


Cool Cube™

Best Practices

Call for
Technical Support
(608) 526-6901



- Prep the PCM (phase change material) panels before use according to one of the described methods provided by VeriCor.
- Ensure all components are clean and free of damage.
- Lay panels flat when turning them solid (to disperse liquid throughout the panel).
- Enable ample air flow around all panel sides.
 - Use spacers (pencils) or racks. →
- Freezing/melting times vary depending on number of panels being prepped and equipment being used.
- Assemble using all six panels for maximum hold time.
 - Using less panels does not change the holding temperature, but does decrease the hold time.
- Panels are reusable (10,000+ cycles).
 - End-of-life disposal: Panels are a plastic #2, typically recycled by businesses/communities. PCM is nontoxic and readily biodegradable.
- Use a calibrated data logger or other temperature monitoring device to observe internal temperature.
- Avoid unnecessary opening of the Cool Cube™ after loading payload. Opening of the Cool Cube™ will decrease hold time.
- An infrared temperature thermometer can assist in ensuring the panels reach a safe pack-out temperature (good for finding out the approximate temperature of each panel).
- The farther the ambient temperatures are from the melting point, the quicker PCM will change states (solidify/liquefy).



Various methods based on type of panel, equipment available & purpose.

Cool Cube™ Room Temp PCM Panels

for FFPE, platelets, biospecimens & more.

All Sizes



Tan Tab/Label



PCM Panel for **ROOM** temps
(20-23°C / 68-73°F)

Remove → Prep (Lab Incubator (see User Guide)) → Assemble → Pack-Out → Close

Part: CC-PCMP-R102 S/N


Do not expose to extreme heat (27°C/81°F) or use otherwise directed.

tel: +1 608 526 6901
VeriCorMed.com/PCM

Video



Prep Method B: Fridge/Room Prep to keep product cool



DO NOT assemble panels directly from a fridge, as they may be initially below 20°C.


Panel Prep

- 2.1 Lay panels flat in a refrigerator until all the PCM (phase change material inside the panel) turns solid. At 4°C/39°F the PCM will solidify in a couple hours.
- 2.2 Spread panels out (enable airflow to all sides) in a room just before use to allow the PCM inside to rise to the appropriate operating temperature. Approximate times:

"03" size = 35 minutes	<i>Times based on a</i>	"28" size = 45 minutes
"08" size = 40 minutes	<i>22°C/72°F room.</i>	"96" size = 50 minutes
- 2.3 Shake panels to verify the PCM is solid with just a little liquid. If a little liquid is heard, it is at 21.5°C. If there is a lot of liquid, restart at step 2.1 to ensure the longest hold time. Using liquid PCM or panels with a solid/liquid combination decreases the hold time. Wipe off condensate and proceed with assembly.



PCM Panel Shake Test



ISTA 7D Thermal Performance Study

Lab-Qualified Hold Times When Starting with Solid PCM

		Qualified Temps: 15-25°C	20-24°C
Cool Cube™ 03	Utilizing Six (6)	91 hrs	47 hrs
Cool Cube™ 08	Lab Freezer Temp	83 hrs	66 hrs
Cool Cube™ 28	PCM Panels	141 hrs	85 hrs
Cool Cube™ 96	(Tan Tab/Label)	143 hrs	91 hrs

Times listed are based on lab-validated, 24-hour cycles of summer & winter profiles (hot & cold ambient temperatures) without the additional thermal mass of a payload, which if conditioned properly, will improve hold times. Actual performance times may vary.