

VeriCor Response Beds

VeriCor Response Beds support your emergency preparation and mitigation plans. Meeting the FEMA guidelines and ADA compliant, VeriCor's line of Response Beds offer EP planners the versatility to choose the level of bed appropriate to the desired care function. From the basic General-Purpose bed to the advanced Emergency Room bed, our Response-Friendly™ beds are made in the USA with quality you and your patients can depend on.



	General-Purpose (RB-GP400)	Medical Response (RB-FN400I)	Active Patient Care (RB-APC400ISR)	Bariatric (RB-B700)	Emergency Room (RB-ER)
Description	Basic bed for rest (used by the general population/staff).	Bed for general patient care (including those with functional needs).	Bed for continuous patient care (standard working height).	Bed for larger individuals or multiple people.	Hard deck for advanced ER care (CPR, surgery, etc.).
Working Capacity	450 lbs.	450 lbs.	450 lbs.	700 lbs.	500 lbs.
Height	22"	18"	28"	22"	Adjustable
Weight	18 lbs.	33 lbs.	38 lbs.	36 lbs.	64 lbs.
Head/Foot Elevator	No	Both	Both	No	Head Only
Included Accessories	-	IV Pole	IV Pole, Side Rails	-	IV Pole, Rails, Casters
Optional Accessories	Hanging Drawer	Side Rails, Hanging Drawer	Hanging Drawer	Hanging Drawer	Hanging Drawer

**Linens sold separately*

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com



General-Purpose Bed

Model: RB-GP400

The General-Purpose Response Bed is engineered for the rigors of long or short term deployment. The lightweight, durable design is comfortable and portable, making it an excellent choice for the general population and staff needs.



Features

- Lightweight and portable
- Spun tube ends
- Waterproof, anti-microbial mattress
- Comfortable spring and deck design
- Heavy duty locking legs

Specifications

- 400+ lbs. working capacity
- Sleeping Surface: 32" x 80"
- Height: 22"
- Folded: 42" x 34" x 4½"
- Weight: 18 lbs.
- Meets FED STD 1633 for flammability

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com



Active Patient Care Bed w/ IV Pole

Model: RB-APC400I

The Active Patient Care (APC) Bed is designed with patient surge treatment care needs in mind. This taller, procedure height bed allows healthcare professionals to perform necessary functions at a working height. In the event of an overflow, the APC doubles effectively as a temporary cot.



Folded (Dual View)



* side rails option



Features

- Lightweight and portable
- 5 position head elevator
- 2 position foot elevator
- Anti-microbial mattress
- Spring and deck design
- Patient restraint belt
- Stow on board IV Pole

Specifications

- 450+ lbs. working capacity
- Sleeping Surface: 32" x 80"
- Height: 28"
- Folded: 34" x 42" x 10"
- Weight: 38 lbs.
- Meets FED STD 1633 for flammability

Options	Model #
w/ IV Pole & Rails	RB-APC400ISR

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com



Functional Needs Bed

Model: RB-FN400

The Functional Needs (FN) Bed is designed with patient surge capacity care needs in mind. Meeting FEMA guidelines and ADA compliant, this bed offers the patient a sense of safety and comfort. It has a wheelchair transferable height and a decontaminable mattress made to reuse over extended periods of time.



* IV pole and side rails options

Folded (Dual View)

Features

- Lightweight and portable
- 5 position head elevator
- 2 position foot elevator
- Anti-microbial mattress
- Spring and deck design
- Patient restraint belt
- Heavy duty locking legs

Specifications

- 450+ lbs. working capacity
- Sleeping Surface: 32" x 80"
- Height: 18"
- Folded: 34" x 42" x 8"
- Weight: 32 lbs.
- Meets FED STD 1633 for flammability



Options	Model #
w/ IV Pole	RB-FN400I
w/ IV Pole & Rails	RB-FN400ISR

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com



Bariatric Bed, Wt Cap 700lb

Model: RB-B700

The Bariatric Response Bed is a lightweight, portable, reusable cot that is highly versatile. The extra sturdy design offers the bariatric population a safe and comfortable bed in times of need. In addition, the taller, longer design has been found to accommodate multiple users, maximizing the limited space of a patient surge situation.



Folded (Dual View)



* linen set sold separately

Features

- Lightweight and portable
- Spun tube ends
- Waterproof, anti-microbial mattress
- Comfortable spring and deck design
- Heavy duty locking legs

Specifications

- 700+ lbs. working capacity
- Sleeping Surface: 38" x 84"
- Height: 22"
- Folded: 42" x 42" x 7"
- Weight: 36 lbs.
- Meets FED STD 1633 for flammability

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com



Emergency Room Bed

Model: RB-ER

The Emergency Room Response Bed is the ideal bed for mobile emergency rooms. The heavy-duty construction withstands the rigors of field medicine, while the features enable safe, efficient treatment. Push-button assembly allows for a quick setup without tools. After use, all parts can be easily disinfected and decontaminated for future use.



Features

- Hard deck for rigorous medical procedures
- Adjustable backrest
- Single-hook IV pole
- Removable side rails
- Anti-microbial mattress
- Set of 4 locking casters

Specifications

- 500 lb. working capacity
- Sleeping Surface: 31" x 81"
- Height: 19" - 32" (adjustable)
- Stored: 82" x 31" x 3"
- Weight: 64 lbs.
- ADA compliant



[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®