

3-Drawer Mobile Medical Case

Model: MC-3D-__

4 Color Choices - Please provide color code when ordering

Red (RD)	OD Green (OD)
Tan (TN)	Black (BK)

The 3-Drawer Mobile Medical Case is a storage container, cabinet and mobile case all-in-one. This durable case provides deep widths and easy access to vital medical gear in the field. It is designed specifically to keep life saving equipment/supplies safe and secure from moisture, dust, sand, temperature and impact.



Features

- Top sliding shelf with mesh & Velcro® closure
- Extra large zippered pouch
- Built-in wheels
- Interchangeable drawers
- Chemical agent decontaminable
- Stackable


Specifications

Dimensions
34" x 21" x 19"

Weight
39 lbs.

Interior Cubic Volume
5.38 cubic feet

Drawer Size (each)
16" x 14³/₄" x 6¹/₂"



Butterfly Latches



[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®