

VeriCor Response Rolling Racks

- Welded Steel Construction
- Laser Cut Components
- Powder Coated Textured Finish
- 2 Fixed & 2 Swivel Casters
- Dust Cover Available

VeriCor's Response Rolling Racks support your emergency preparation and mitigation plans. These racks are constructed of heavy duty steel with a baked on powder coat finish to protect from dings and scratches of exercises and/or deployment. Swivel casters ensure dynamic maneuverability and ease of movement on almost any surface.



	Rolling Rack for 6 Bariatric Beds (RR-6B)	Rolling Rack for 8 APC Beds (RR-8APC)	Rolling Rack for 8 Response Beds (RR-8)	Rolling Rack for 10 Response Beds (RR10)	Rolling Racks for Beds – Custom (RR-CBR)	Rolling Rack for Mattresses – Custom (RR-CMR)
Width x Depth x Height	39" x 43" x 45"	38" x 43" x 80"	38" x 43" x 68"	38" x 43" x 80"	varies	varies
Number of Shelves	1	4	4	5	varies	varies
Bed (Model #)	-----Rolling Rack Capacity-----					
General-Purpose (RB-GP400)	~	~	12	15	limit 30	limit 30
Functional Needs (RB-FN400)	~	~	8	10	limit 30	limit 30
Active Patient Care (RB-APC400I)	~	8	~	~	limit 30	limit 30
Medical Response Bed (RB-MR400)	~	~	8	8	limit 30	limit 30
Bariatric (RB-B700)	6	~	~	~	limit 30	limit 30

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

