

# 42" Workstation w/ 3 Drawers, Table & Refrigerator

**Model:** MW-42R

VeriCor's 42" Workstation with Refrigerator, 3 Drawers & Table comes standard with ample storage, a built-in refrigerator and a detachable table. This single unit station becomes mobile in an instant. Pack in all your medical supplies and the workstation is ready for your next vaccination campaign.

*Integrated Refrigerator*



*Locking Lid Stays*



*Swivel Wheels*



## Features

- Refrigerator without freezer
- Lockable refrigerator & case
- Tilt bins for easy organization
- Easy lock wheels with a wide wheel base for easy and safe maneuverability
- Custom labeling included

## Specifications

*Dimensions (H x W x D)*  
48" x 43" x 25"

*Drawers*  
3 (100 lbs. capacity each)

*Fold-Out Table*  
35" x 43" (30" tall)

*Weight*  
283 lbs.

*Please provide color code when ordering:*

<b>COLOR (CODE)</b>	<b>Red (RD)</b>	<b>Yellow (YL)</b>	<b>Forest Green (FG)</b>	<b>Navy Blue (NB)</b>	<b>Black (BK)</b>
<b>Silver (SL)</b>	<b>Orange (OR)</b>	<b>Green (GR)</b>	<b>Violet Blue (BL)</b>	<b>Wine (WN)</b>	<b>Tan (TN)</b>

**Go to Web Page >>**



info@VeriCorMed.com  
(866) 469-6019

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

**VERICOR** LLC  
Response-Friendly®