

# Cool Cube™ PCM Coolers for Blood Products

A Cool Cube™ at Fridge Temperatures is not an ordinary transport cooler. It is a portable, validated container and pack-out system that keeps temperature sensitive product (whole blood, RBCs, liquid plasma, etc.) chilled between 1 and 10 °C for over 2½ days (or between 1 and 6 °C for over 1½ days) – *without using ice or electricity!*

**Temperature Range** | 1 to 6 °C (34 to 42 °F) / 1 to 10 °C (34 to 50 °F)



	Cool Cube™ 03 at Fridge Temps (BT-03)	Cool Cube™ 08 at Fridge Temps (BT-08)	Cool Cube™ 28 at Fridge Temps (BT-28)	Cool Cube™ 96 at Fridge Temps (BT-96)
<b>Validated Duration</b>	39+ hrs. / 70+ hrs.	53+ hrs. / 83+ hrs.	68+ hrs. / 108+ hrs.	112+ hrs. / 128+ hrs.
<b>Payload Area</b>	5¼" x 5¼" x 5¼"	11¼" x 6½" x 6½"	12" x 12" x 12"	18" x 18" x 18"
<b>Exterior Dimensions</b>	11" x 11" x 11"	16" x 11½" x 11½"	20½" x 20½" x 23"	26½" x 26½" x 28½"
<b>Weight</b>	11 lbs.	16 lbs.	61 lbs.	111 lbs.
<b>Features</b>	<ul style="list-style-type: none"> <li>• Six Refrigerator PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Carrying Case</li> </ul>	<ul style="list-style-type: none"> <li>• Six Refrigerator PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Carrying Case</li> </ul>	<ul style="list-style-type: none"> <li>• Six Refrigerator PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Wheeled Case</li> </ul>	<ul style="list-style-type: none"> <li>• Six Refrigerator PCM Panels</li> <li>• Temp-Shield™ VIP System</li> <li>• User-Friendly Wheeled Case</li> </ul>
	<a href="#">Web/Pricing</a>	<a href="#">Web/Pricing</a>	<a href="#">Web/Pricing</a>	<a href="#">Web/Pricing</a>