


FFPE, Platelets & Biospecimens - Cool Cube™ PCM Transport Coolers



A Cool Cube™ at Room Temperatures is not an ordinary transport cooler. It is a portable, validated container and pack-out system that keeps temperature sensitive product (FFPE, platelets, biospecimens, etc.) between 15 and 25 °C for over 3 days (or between 20 and 24 °C for over 2 days) – *without using ice or electricity!* This simple system streamlines the transport of valuable, room temperature product and protects it from warming or cooling during off-site use or power outages. The validated, patented technology of a Cool Cube™ saves time, saves money and is improving cold chain management practices around the world. The Cool Cube™ is a container and pack-out system so advanced...it's easy!



	Cool Cube™ 03 at Room Temps (RT-03)	Cool Cube™ 08 at Room Temps (RT-08)	Cool Cube™ 28 at Room Temps (RT-28)	Cool Cube™ 96 at Room Temps (RT-96)
Temperature Range	20 to 24 °C / 15 to 25 °C	20 to 24 °C / 15 to 25 °C	20 to 24 °C / 15 to 25 °C	20 to 24 °C / 15 to 25 °C
Validated Duration	47+ hrs. / 91+ hrs. (summer)	66+ hrs. / 83+ hrs. (summer)	85+ hrs. / 141+ hrs. (summer)	91+ hrs. / 143+ hrs. (summer)
Payload Area	5¾" x 5¾" x 5¾"	11¼" x 6½" x 6½"	12" x 12" x 12"	18" x 18" x 18"
Exterior Dimensions	11" x 11" x 11"	16" x 11½" x 11½"	20½" x 20½" x 23"	26½" x 26½" x 28¼"
Weight	11 lbs.	17 lbs.	62 lbs.	114 lbs.
Features	<ul style="list-style-type: none"> • Six Room Temp PCM Panels • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • Six Room Temp PCM Panels • Temp-Shield™ VIP System • User-Friendly Carrying Case 	<ul style="list-style-type: none"> • Six Room Temp PCM Panels • Temp-Shield™ VIP System • User-Friendly Wheeled Case 	<ul style="list-style-type: none"> • Six Room Temp PCM Panels • Temp-Shield™ VIP System • User-Friendly Wheeled Case

[Go to Web Page >>](#)



info@VeriCorMed.com
(608) 526-6901

For more information & pricing: www.VeriCorMed.com

