

# Equipment Chest

Model: MC-EC-\_\_

With an enormous capacity, the Equipment Chest has the space you need to get medical supplies organized. This durable case protects assets during storage and enables the mass movement/deployment of vital medical gear in the field. The removable trays are just the right size to organize smaller items, while the large payload area is ideal for bulky medical equipment. Pack it up and be ready to respond.



Locking Lid Stays



Payload Area



(supplies not included)

## Features

- Large capacity
- Wide wheel base
- Four swivel casters (2 Locking)
- Two removable trays
- Locking lid stays
- Four recessed handles

## Specifications

*Dimensions (L x W x H)*  
50" x 32" x 48"

*Weight*  
230 lbs.

*Casters*  
6" diameter

*Used in Response-Friendly™ Triage & Treatment modules.*

Please provide color code when ordering:

<b>COLOR (CODE)</b>	<b>Red (RD)</b>	<b>Yellow (YL)</b>	<b>Forest Green (FG)</b>	<b>Navy Blue (NB)</b>	<b>Black (BK)</b>
<b>Silver (SL)</b>	<b>Orange (OR)</b>	<b>Green (GR)</b>	<b>Violet Blue (BL)</b>	<b>Wine (WN)</b>	<b>Tan (TN)</b>



[Go to Web Page >>](#)



info@VeriCorMed.com  
(608) 526-6901

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

