

Examination MedKit

Model: MK-E

The Examination MedKit is an UMO (Upright Medical Organizer) loaded with essential components needed for monitoring vital signs and carrying out physical exams. Products include sphygmomanometer kits, stethoscopes, ophthalmoscope set, probes, thermometers, etc.



Features

- The UMO (Upright Medical Organizer)
- Easy to assemble stand
- Items specific for patient examinations
- Packaging choice of rolling hard case or soft carry case
- SmartBook™ with critical data

Specifications

Dimensions On Stand
73" x 39" x 25" @ 100 lbs.

Storage Pockets
21 clear

Rolling Hard Case Dimensions
38" x 23" x 23"

Soft Carry Case Dimensions
30" x 15" x 15"



You choose
the packaging!

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®