

External Pocket Kit of Cool Cube™ 28

Model: CC-EPRK

VeriCor's External Pocket Kit of the Cool Cube™ 28 is a retrofit kit for any tan, hard-cased VeriCor Cool Cube™ that was purchased prior to 2016. This kit enables the addition of an external pocket and flexible grommet to the outside of the case. These two features (now standard on the large Cool Cubes™) enable the user to keep a monitor on the outside of the cooler during transport. Some assembly is required.



Features

- a) 1 - External Pocket
- b) 1 - 1/4"-3/4" Step Drill Bit
- c) 1 - 3/4" Flexible Grommet
- d) 1 - Center Punch
- e) 1 - 9/64" Drill Bit
- f) 3 - Machine Screw
- g) 3 - Rubber Grommet
- h) 1 - Threadlocker
- i) 4 - Acorn Nut (1 extra)
- j) 1 - Pocket Card (not pictured)

Specifications

TOOLS NEEDED (not included)

- Channellock Pliers (or similar)
- Drill
- Hammer
- Phillips-Head Screwdriver
- Small Wrench (Crescent or 5/16")



External Pocket

Flexible Grommet
(former location of
the
pressure release
valve)

Go to Web Page >>



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®