

General-Purpose Bed

Model: RB-GP400

The General-Purpose Response Bed is engineered for the rigors of long or short term deployment. The lightweight, durable design is comfortable and portable, making it an excellent choice for the general population and staff needs.



Features

- Lightweight and portable
- Spun tube ends
- Waterproof, anti-microbial mattress
- Comfortable spring and deck design
- Heavy duty locking legs

Specifications

- 400+ lbs. working capacity
- Sleeping Surface: 32" x 80"
- Height: 22"
- Folded: 42" x 34" x 4½"
- Weight: 18 lbs.
- Meets FED STD 1633 for flammability

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

