

Go-Bag (Bag Only)

Model: MC-GB

The VeriCor Go-Bag gives medical supplies the organization needed for a quick response. This is the same bag used to hold VeriCor's 500 person vaccination system. You get in-line skate wheels, a retractable locking pull handle, padded carry handles, and an adjustable buckle strap to compliment the spacious zippered main compartment, front and side zipper pockets, and the hard-sided separate bottom section.

Features

- Four heavy-duty handles
- Tear-resistant nylon
- Self-repairing nylon zippers
- In-line skate wheels
- Spacious main compartment
- Retractable locking pull handle

Specifications

Dimensions (L x H x D)
31½" x 15" x 14½"

Storage Compartments
5: 1 spacious main, 1 hard-sided bottom, 2 side, 1 front (with smaller pockets)

Weight
8 lbs.



Wheels w/ Curb Guards



Telescopic Handle



[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®