

Litter Equipment Organizer

Model: MC-LEO

The Litter Equipment Organizer protects litters, stands, pads and straps during storage and enables the mass movement/deployment of bulk supplies. Large casters enable this organizer to move easily on carpet, vinyl, wood and concrete flooring. The welded steel construction has no sharp edges and is ideal for the protection of medical assets.



Features

- Capacity: 17 litters, 34 stands, 17 pads, 34+ straps
- Ten hinged doors
- Two steering handles
- Four locking swivel casters
- See-through design
- High quality epoxy paint

Specifications

Dimensions
104" x 48" x 70"

Weight
540 lbs.

Frame Material
1¼" square steel tubing

Casters
6" diameter



Door Latch



Door Handle



Steering Handle



Locking Wheel



[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

