

Medical Organizer w/ 4 Drawers & Table

Model: MC-MO-__

The Medical Organizer w/ 4-Drawers & Table is a storage container, cabinet, table and medical case all-in-one. This durable case provides partitioned drawers for easy access to vital medical gear in the field. The detachable table provides a Response-Friendly™ solution to needing an instant work surface.



Features

- Detachable table
- Deep drawers
- Customizable drawer dividers
- Easy-slide drawer guides
- Drawer label holders
- Butterfly latches
- Wide wheel base

Specifications

Dimensions (H x W x D)
36" x 26" x 23" (closed)

Weight
87 lbs.


Storage
4 partitioned drawers

Detachable Table
35" x 24" (30" tall)

Built-In Wheels



Drawer Dividers



Please provide color code when ordering:

COLOR (CODE)	Red (RD)	Yellow (YL)	Forest Green (FG)	Navy Blue (NB)	Black (BK)
Silver (SL)	Orange (OR)	Green (GR)	Violet Blue (BL)	Wine (WN)	Tan (TN)

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

