

Oxygen Distribution (for Small DUV)

Model: SS-DUV-O2SM

The Oxygen Distribution Specialized Response Module for the Small DUV includes the components necessary to deliver oxygen to patients within a small Dual Use Vehicle. All components are contained within a wheeled, labeled, protective case for safe storage and easy movement. The Cascade System and manifold attach to the L-track within DUV for rapid set-up and safe oxygen distribution. Oxygen tanks not included.

Module Features

- 1 Manifold (8 patients)
- 3 Oxygen Cascade Systems
 - 3 Patient Drops Each
- O₂ Wrench
- O₂ Warning Sign
- 1 Medical Equipment Response Case (MERC)
- 1 SmartBook™ Response Manual
- Dimensions
 - 32" x 21" x 13" @ 58 lbs.

Ask a sales representative about other included items.

Designed as a supplementary set for the VA Dual Use Vehicles.
**SmartBook™ operator's manual included!*

[Go to Web Page >>](#)



Protective MERC
(open view)



SS-SS-DUV-O2SM-012919

info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly™