

Portable Stand for MC-UMO

Model: MC-SUMO

VeriCor's Portable Stand for MC-UMO (the "SUMO") typifies Response-Friendly™ medical equipment. This portable stand, with its snap-together design, sets up in under a minute. Machined aluminum, detention balls, and interchangeable parts enable an easy fit, while the stability straps give it the strength to handle any large, hanging, medical bag.*

**Upright Medical Organizer sold separately.*



Disassembled Stand
w/ Carry Bag



Features

- Lightweight aluminum
- Intuitive design assembles quickly and easily
- High quality powder-coated parts
- Two stability straps
- Carry bag included

Specifications


Dimensions (H x W x D)
73" x 39" x 25"

Weight
9 lbs.

Carry Bag Weight
1.5 lbs.

Carry Bag Dims (L x W)
30" x 24"

Heavy-Duty Welds



[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

