
URGENT MEDICAL DEVICE RECALL

June 19, 2023

Dear Valued McKesson Customer:

Teleflex LLC has notified McKesson Medical-Surgical Inc. (MMS) of an Urgent Medical Device Recall regarding specific lots of their Endotracheal Tubes (ET tube). This notice has been issued due to reports of disconnection of the 15mm connector from the endotracheal tube (ET tube) for the affected products. There is the possibility of oxygen desaturation and in the event any immediate and long-term health consequences are dependent on the degree and duration of desaturation, which may include serious injury or death. Affected product first shipped June 14, 2018.

McKesson Medical-Surgical Inc. has taken appropriate action per this notice.

For questions regarding this notification, please contact Teleflex Customer Service at **(866) 396-2111**.

A review of our records indicates that your company may have purchased items included in this notification. Carefully review the information in this letter and follow the instructions provided below.

Refer to Attachment 1 for a list of affected item(s) distributed by McKesson Medical-Surgical

McKesson Customer Instructions:

- 1.) Immediately discontinue use of any product matching the affected item(s) and lot number(s) listed in Attachment 1. If you have no products matching the affected item(s) and lot number(s), no further action is needed.
- 2.) A copy of the Urgent Medical Device Recall from Teleflex LLC has been included for reference.
- 3.) If you have product affected by this notice, fill out the McKesson Reply Form and return it to our Corporate Customer Service Center via email at MMSRecalls@McKesson.com or fax at **(866) 871-0270**. To ensure timely credit to your account and support the completion of this notice, please respond within 30 days.
 - **Please note:** Any product returned in addition to or in lieu of affected product will be destroyed, without issuance of a credit. The affected product lot numbers are listed in the item table. Once the product is returned, credit will be issued to you.
 - **Please place a new order for replacement product if there is an immediate need.**
- 4.) If you have further distributed any of the item(s) referenced in this notification, provide your accounts with a copy of this Urgent Medical Device Recall and request that they return the affected product directly to you.

We sincerely apologize for any inconvenience this notice may have caused you and your staff. If you have any questions about information provided in this communication, please contact our **McKesson Medical-Surgical Recall Message Center** at MMSRecalls@McKesson.com or call **(800) 688-8840**.

Thank you for your prompt attention,

McKesson Medical-Surgical Inc.

ATTACHMENT 1

MMS #	MFG Catalog #	Description	Affected Lot(s)			
130596	112082070	TUBE, MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX	18FG09	19JT61	KME21C0737	KME22A0699
			18FG12	19JT67	KME21C0822	KME22A0710
			18FG14	19KT14	KME21C0854	KME22A0991
			18FG27	19KT18	KME21C0971	KME22A0992
			18FG33	19KT19	KME21C0973	KME22A1045
			18FT08	19KT20	KME21C1279	KME22A1133
			18GG01	19KT26	KME21C1328	KME22A1776
			18GG02	19KT30	KME21C1655	KME22A2963
			18GG03	19KT31	KME21C1702	KME22A3194
			18GG05	19KT42	KME21C2150	KME22B0630
			18GG07	19LT13	KME21C2403	KME22B1947
			18GG08	19LT31	KME21C2568	KME22B2203
			18GG10	19LT37	KME21C2875	KME22B2316
			18GG12	20AT05	KME21C2954	KME22B2410
			18GG19	20AT11	KME21D0350	KME22B2585
			18GG22	20AT25	KME21D0435	KME22C1342
			18GT07	20AT26	KME21D0491	KME22C2121
			18GT08	20AT27	KME21D0714	KME22C2755
			18GT12	20AT32	KME21D0715	KME22C2854
			18HG28	20AT38	KME21D1218	KME22C3272
			18HG34	20AT49	KME21D1409	KME22D0480
			18HG36	20AT58	KME21D1580	KME22D1175
			18HG40	20BG14	KME21D1652	KME22D1826
			18HT07	20BT02	KME21D1706	KME22D1859
			18HT08	20BT09	KME21D1878	KME22D2302
			18HT31	20BT13	KME21D2020	KME22D2405
			18IG01	20BT22	KME21D2929	KME22D3021
			18IG04	20BT27	KME21D2935	KME22D3274
			18IG07	20BT31	KME21E0001	KME22D3673
			18IG09	20BT37	KME21E0002	KME22D3674
			18IG10	20BT38	KME21E0363	KME22E0177
			18IG11	20BT39	KME21E1390	KME22E0317
			18IG15	20BT46	KME21F0102	KME22E0640
			18IG16	20BT51	KME21F0140	KME22E0655
			18IG18	20BT52	KME21F0239	KME22E0980
			18IG19	20BT53	KME21F0434	KME22E1974
			18IG22	20CT02	KME21F1379	KME22E2108
			18IG23	20CT07	KME21F1887	KME22E2110
			18IG26	20DT05	KME21F2076	KME22E2272
			18IT04	20DT13	KME21F2153	KME22E2559
18IT12	20DT17	KME21G0013	KME22E2560			
18IT16	20DT22	KME21G0662	KME22E2809			
18IT17	20DT23	KME21G1616	KME22E2919			
18IT18	20DT40	KME21G1663	KME22F0123			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
130596	112082070	TUBE, MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX	18IT23	20ET12	KME21G1741	KME22F0279
			18JG04	20ET15	KME21G1923	KME22F0378
			18KG07	20ET17	KME21G1957	KME22F0523
			18KG10	20ET43	KME21G2511	KME22F0792
			18KG25	20ET44	KME21H0975	KME22F0870
			18KG33	20FT27	KME21H1183	KME22F0907
			18KG37	20FT34	KME21H1265	KME22F2117
			18KG40	20FT39	KME21J0157	KME22F2155
			18KG41	20FT43	KME21J0238	KME22F2905
			18KT03	20FT45	KME21J0497	KME22G0012
			18KT30	20FT53	KME21J0639	KME22G0787
			18LG11	20FT54	KME21J0756	KME22G0800
			18LG13	20FT55	KME21J0774	KME22G0809
			18LG16	20FT66	KME21J0855	KME22G0879
			18LG19	20FT71	KME21J1012	KME22G1036
			18LG35	20GG39	KME21J1207	KME22G1949
			18LG40	20GT09	KME21J1256	KME22H1166
			18LT13	20GT12	KME21J1540	KME22H1517
			19AG02	20GT13	KME21J1605	KME22H1995
			19AG03	20GT18	KME21J1885	KME22H2075
			19AT16	20GT31	KME21J2157	KME22H2377
			19BG02	20GT32	KME21K0003	KME22H2379
			19BG14	20GT41	KME21K0201	KME22H2458
			19BG16	KME20H0916	KME21K0241	KME22H2683
			19BG19	KME20H1131	KME21K0388	KME22H2813
			19BG21	KME20H1244	KME21K0553	KME22J0293
			19BG22	KME20H1245	KME21K0635	KME22J0703
			19BG23	KME20H1297	KME21K0931	KME22J0826
			19BG26	KME20H1393	KME21K1129	KME22J1685
			19BG33	KME20H1538	KME21K1202	KME22J3142
			19BG34	KME20H1634	KME21K1334	KME22J3155
			19BG37	KME20H1635	KME21K1415	KME22J3170
			19BG42	KME20H2405	KME21K1557	KME22J3259
			19BT12	KME20H2727	KME21K1618	KME22J3276
			19BT16	KME20J0262	KME21K1667	KME22J3326
			19BT18	KME20J1281	KME21K1816	KME22J3403
			19CG01	KME20J1387	KME21K2169	KME22K0036
			19CG03	KME20J1398	KME21K2270	KME22K0137
			19CG06	KME20J1638	KME21K2533	KME22K0198
			19CT03	KME20J1804	KME21K2571	KME22K0254
			19CT22	KME20J1898	KME21K2638	KME22K0365
			19CT39	KME20J2158	KME21K2719	KME22K0407
			19CT40	KME20J2711	KME21K2746	KME22K0415
			19CT42	KME20J2783	KME21K2798	KME22K0509
			19CT45	KME20J2802	KME21K2864	KME22K2064
19CT72	KME20J2828	KME21K3009	KME22K2087			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
130596	112082070	TUBE, MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX	19CT82	KME20J2874	KME21K3010	KME22K2238
			19DT15	KME20J2896	KME21K3011	KME22K2267
			19DT17	KME20J2901	KME21K3020	KME22K2434
			19DT20	KME20J2942	KME21K3021	KME22L1283
			19DT29	KME20J3023	KME21L0812	KME22L1802
			19DT36	KME20J3160	KME21L1055	KME22L2194
			19DT46	KME20J3266	KME21L1187	KME22L2489
			19DT47	KME20K0061	KME21L1215	KME22L2497
			19ET02	KME20K0437	KME21L1279	KME22L2524
			19ET11	KME20K0462	KME21L1320	KME22L2584
			19ET12	KME20K0585	KME21L1367	KME22L2709
			19ET18	KME20K1007	KME21L1930	KME22M2259
			19ET29	KME20K3013	KME21L2128	KME22M2283
			19ET31	KME20K3014	KME21M0850	KME22M2695
			19ET34	KME20K3145	KME21M1194	KME22M2748
			19ET36	KME20K3365	KME21M1390	KME22M2839
			19ET40	KME20K3448	KME21M1444	KME22M2872
			19ET47	KME20K3518	KME21M1611	KME22M2942
			19ET49	KME20K3560	KME21M1627	KME22M2980
			19ET65	KME20L0276	KME21M1634	KME22M3029
			19ET78	KME20L0371	KME21M1674	KME22M3066
			19ET82	KME20L0785	KME21M1698	KME23A0809
			19FT04	KME20L0928	KME21M1726	KME23A1197
			19FT29	KME20L1691	KME21M1746	KME23A2930
			19FT38	KME20L1752	KME21M1864	KME23A2931
			19FT50	KME20L1858	KME21M1909	KME23A3023
			19GT02	KME20L1874	KME21M2007	KME23B0129
			19GT16	KME20L1923	KME21M2032	KME23B0176
			19GT38	KME20L2111	KME21M2072	KME23B0205
			19GT47	KME20L2210	KME21M2116	KME23B0228
			19GT51	KME20L2273	KME21M2136	KME23B0267
			19GT54	KME20L2350	KME21M2147	KME23B0357
			19GT58	KME20L2373	KME21M2321	KME23B0378
			19GT59	KME20L2463	KME21M2363	KME23B0434
			19GT65	KME20L2478	KME21M2644	KME23B0475
			19HT02	KME20L2499	KME21M2700	KME23B0852
			19HT08	KME20L2699	KME21M2709	KME23B0940
			19HT14	KME20L2826	KME21M2727	KME23B0942
			19HT17	KME20M0490	KME21M2732	KME23B0964
			19HT22	KME20M0664	KME21M2744	KME23B1068
19HT23	KME20M0712	KME21M2783	KME23B1079			
19HT25	KME20M0762	KME21M2968	KME23B1103			
19HT47	KME20M0774	KME21M3075	KME23B1111			
19HT50	KME20M0839	KME21M3138	KME23B1180			
19HT52	KME21A1121	KME21M3181	KME23B1523			
19HT66	KME21A1215	KME21M3224	KME23B1524			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
130596	112082070	TUBE, MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX	19HT69	KME21A1328	KME21M3225	KME23B1569
			19HT78	KME21A2220	KME21M3274	KME23B1662
			19HT82	KME21A2345	KME22A0046	KME23B1663
			19IT02	KME21A2574	KME22A0099	KME23B1701
			19IT04	KME21A2984	KME22A0224	KME23B1800
			19IT15	KME21A3038	KME22A0247	KME23B2428
			19IT31	KME21B0136	KME22A0356	KME23B2554
			19IT42	KME21B0711	KME22A0429	KME23B2980
			19JT05	KME21B2898	KME22A0513	KME23C1764
			19JT47	KME21C0505	KME22A0552	KME23C2325
			19JT55	KME21C0567	KME22A0621	KME23C2356
			19JT60	KME21C0701		
			132026	112080080	TUBE, ENDOTRACH NASORAL 8.0MM (10/BX)112480 TELFLX	18FG02
18GG15	19DT19	19LT47				KME21B2727
18IG13	19GT24	20AT06				KME21C0500
18IG22	19HT26	20BT09				KME21G0866
18JG17	19JT20	20CT14				KME22A2268
18JG38	19KT19	20DT23				KME22B0062
18KG22	19KT38	20GG39				KME22F2858
19BG21	19KT50					
134105	112082080	TUBE, MURPHY ENDOTRACH CUFF 8.0MM (10/BX) TELFLX	18FG16	20AT49	KME21B2489	KME22A2492
			18FT08	20BG15	KME21C0705	KME22A2693
			18FT24	20BG16	KME21C0706	KME22A2932
			18FT26	20BG19	KME21C1063	KME22A2947
			18GG06	20BT05	KME21C1325	KME22A2956
			18GG08	20BT09	KME21C1441	KME22A2966
			18GG11	20BT10	KME21C1745	KME22A2967
			18GG12	20BT13	KME21C1785	KME22A3223
			18GG13	20BT24	KME21C2068	KME22A3241
			18GG16	20BT27	KME21C3283	KME22B0118
			18GT06	20BT31	KME21C3320	KME22B0500
			18GT08	20BT37	KME21D1653	KME22B0525
			18GT11	20BT51	KME21D1946	KME22B0737
			18GT22	20BT52	KME21D2011	KME22B0898
			18HG28	20CT42	KME21D2012	KME22B0915
			18HG37	20CT48	KME21E0584	KME22B0991
			18HG40	20CT54	KME21E0687	KME22B1148
			18HT27	20DT02	KME21E0770	KME22B1149
			18HT31	20DT04	KME21E1452	KME22C1286
			18IG04	20DT10	KME21F0981	KME22C1362
			18IG07	20DT23	KME21F1006	KME22C1525
			18IG11	20DT28	KME21F1215	KME22C1987
			18IG12	20ET12	KME21F1530	KME22C2072
			18IG13	20ET16	KME21F1604	KME22C2219
			18IG23	20ET17	KME21F2016	KME22C2400
			18IT06	20ET39	KME21G0736	KME22C3376

MMS #	MFG Catalog #	Description	Affected Lot(s)			
13410 5	112082080	TUBE, MURPHY ENDOTRACH CUFF 8.0MM (10/BX) TELFLX	18IT10	20FT17	KME21G1100	KME22D0004
			18IT16	20FT18	KME21H0887	KME22D1589
			18IT18	20FT19	KME21H0919	KME22D1689
			18JG07	20FT21	KME21H0962	KME22D1897
			18JG09	20FT43	KME21H0965	KME22D2014
			18JG12	20FT46	KME21H1260	KME22D2505
			18JG13	20FT49	KME21J0109	KME22D2636
			18JG26	20FT58	KME21J0484	KME22D3062
			18JG37	20FT61	KME21J0563	KME22D3275
			18JG38	20FT64	KME21J1295	KME22D3421
			18JT04	20GG48	KME21J1364	KME22D3423
			18JT06	20GT10	KME21J1415	KME22D3427
			18JT10	20GT11	KME21J1468	KME22D3712
			18KG23	20GT12	KME21J1532	KME22E2014
			18LG08	20GT23	KME21J2055	KME22E2271
			18LG11	KME20G0204	KME21K0022	KME22E2901
			18LT05	KME20G0436	KME21K0187	KME22F0278
			18LT07	KME20H0953	KME21K0408	KME22F0621
			19AG35	KME20H1298	KME21K0449	KME22F1270
			19AG38	KME20H1394	KME21K0636	KME22F2238
			19AT20	KME20H1530	KME21K0637	KME22G1511
			19BG06	KME20H2356	KME21K2457	KME22G1515
			19BG07	KME20H2404	KME21K2534	KME22G1526
			19BG14	KME20H2871	KME21K2657	KME22J1751
			19BG18	KME20J2507	KME21K2757	KME22J1773
			19BG22	KME20J2562	KME21K2884	KME22J2000
			19BG23	KME20J2805	KME21K3013	KME22J3346
			19BG25	KME20J2830	KME21L0744	KME22K0540
			19BG30	KME20J2895	KME21L0850	KME22K0541
			19BG37	KME20J2945	KME21L1438	KME22K0654
			19BT08	KME20J2950	KME21L1533	KME22K0691
			19BT12	KME20J2991	KME21L1661	KME22K0694
			19DT36	KME20J2998	KME21L1787	KME22K0715
			19DT43	KME20J3341	KME21L2012	KME22K0777
			19ET22	KME20K3520	KME21L2177	KME22K1034
			19ET27	KME20L0194	KME21L2881	KME22L2710
			19ET58	KME20L1157	KME21M2059	KME22L2784
			19ET66	KME20L1237	KME21M2134	KME22M0921
			19ET70	KME20L1706	KME21M2327	KME22M0933
			19ET82	KME20L2207	KME21M2412	KME22M0968
19ET83	KME20L2333	KME21M2636	KME22M1121			
19FT04	KME20L2375	KME21M2915	KME23A0132			
19FT09	KME20L2607	KME21M3032	KME23A0174			
19FT13	KME20L2696	KME21M3072	KME23A0310			
19FT25	KME20M0577	KME21M3137	KME23A0459			
19FT26	KME20M1043	KME21M3180	KME23A0498			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
134105	112082080	TUBE, MURPHY ENDOTRACH CUFF 8.0MM (10/BX) TELFLX	19FT28	KME20M1069	KME21M3236	KME23A0510
			19GT02	KME20M1260	KME21M3273	KME23A0556
			19GT25	KME20M1342	KME22A0044	KME23A0562
			19GT36	KME20M1780	KME22A0109	KME23A0582
			19GT55	KME20M2363	KME22A0394	KME23A0604
			19IT13	KME21A0082	KME22A0437	KME23A0612
			19IT24	KME21A1435	KME22A0547	KME23A0667
			19IT33	KME21A1522	KME22A0626	KME23A0714
			19IT37	KME21A1524	KME22A0793	KME23A0761
			19IT39	KME21A2014	KME22A0855	KME23B1442
			19IT41	KME21A2058	KME22A1014	KME23B1443
			19IT43	KME21A2576	KME22A1043	KME23B1674
			19IT45	KME21B0049	KME22A1132	KME23B2427
			19KT14	KME21B0137	KME22A1255	KME23B3074
			19LT31	KME21B0248	KME22A1701	KME23B3262
			19LT35	KME21B0428	KME22A1752	KME23C0459
			19LT37	KME21B0616	KME22A1787	KME23C0460
			19LT46	KME21B0709	KME22A1991	KME23C0916
			19LT47	KME21B0863	KME22A1992	KME23C1244
			19LT51	KME21B1214	KME22A2130	KME23C1291
			20AT01	KME21B1284	KME22A2203	KME23C1339
			20AT22	KME21B2488	KME22A2429	KME23C1406
144243	112082065	TUBE, MURPHY ENDOTRACH CUFF 6.5MM (10/BX) TELFLX	18FG23	19JT18	KME21A0404	KME22A1234
			18FG27	19JT19	KME21A0695	KME22A1388
			18FG30	19JT20	KME21A1588	KME22A1505
			18FG31	19JT26	KME21A1779	KME22A2005
			18FT06	19JT32	KME21A2102	KME22A2059
			18FT08	19JT36	KME21A2230	KME22A2198
			18GG19	19JT37	KME21A2572	KME22A2316
			18GT32	19JT41	KME21B1352	KME22A2350
			18GT34	19JT47	KME21B2607	KME22A2401
			18HG02	19JT49	KME21B2789	KME22A2416
			18HG04	19JT63	KME21B2856	KME22A2481
			18HG15	19JT67	KME21C0504	KME22A2694
			18HG17	19JT68	KME21C0981	KME22A2726
			18HG30	19KT20	KME21C0983	KME22A2950
			18IG02	19KT21	KME21C1035	KME22A2959
			18IG07	19KT26	KME21C1036	KME22A3224
			18IG10	19LT07	KME21C1097	KME22B0070
			18IG13	19LT31	KME21C1925	KME22B0466
			18IG16	19LT46	KME21C2997	KME22B0566
			18IT04	20AT05	KME21C3166	KME22B0946
			18IT25	20AT06	KME21C3456	KME22B1049
			18JG01	20AT20	KME21E0607	KME22B1945
18JG02	20AT47	KME21E0685	KME22B2102			
18JG04	20AT49	KME21F0384	KME22B2322			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
14424 3	112082065	TUBE, MURPHY ENDOTRACH CUFF 6.5MM (10/BX) TELFLX	18JG07	20BT13	KME21F0516	KME22B2644
			18JG12	20CT02	KME21F0694	KME22C2076
			18JG20	20CT13	KME21F0834	KME22D0407
			18KG35	20CT23	KME21F0930	KME22D0573
			18KG37	20CT37	KME21F1554	KME22D0810
			18KG38	20CT47	KME21F1885	KME22D1482
			18KG44	20CT54	KME21F2156	KME22D1780
			18KT28	20DT01	KME21F2203	KME22D2011
			18LG14	20DT02	KME21G0532	KME22D2630
			19AG09	20DT13	KME21H0202	KME22D3020
			19AG15	20ET15	KME21H0916	KME22E1854
			19AG32	20ET16	KME21H1233	KME22E2488
			19AG33	20ET18	KME21H1264	KME22F1264
			19AT20	20GG27	KME21J0302	KME22F2116
			19BG14	20GT05	KME21J0389	KME22F2172
			19BG18	20GT06	KME21J0528	KME22F2296
			19BG40	20GT13	KME21J0637	KME22F2720
			19CG02	20GT22	KME21J0725	KME22G0552
			19CG06	20GT24	KME21J0909	KME22G1173
			19CG07	20GT32	KME21J1020	KME22G2426
			19CT26	KME20H0882	KME21J1056	KME22G2579
			19CT39	KME20H2726	KME21J1173	KME22G2943
			19CT40	KME20H2934	KME21J1205	KME22G3248
			19CT76	KME20J1386	KME21J1253	KME22H0766
			19DT10	KME20J1399	KME21J1753	KME22H1917
			19DT16	KME20J1656	KME21J1833	KME22H2054
			19DT27	KME20K0036	KME21J1871	KME22H2504
			19DT47	KME20K0259	KME21J1945	KME22J1150
			19ET02	KME20K0436	KME21J1998	KME22J1860
			19ET05	KME20K0504	KME21K0018	KME22J2225
			19ET40	KME20K0677	KME21K0570	KME22J3273
			19ET49	KME20K0881	KME21K0587	KME22J3274
			19ET65	KME20L1421	KME21K0829	KME22K0003
			19GT19	KME20L1436	KME21L0344	KME22K0995
			19GT29	KME20L2349	KME21L0408	KME22K2411
			19GT30	KME20L2447	KME21L0447	KME22K2444
			19GT31	KME20L2448	KME21L0566	KME22K2445
			19GT38	KME20M0710	KME21L0672	KME22K2561
			19HT26	KME20M0892	KME21L1306	KME22L0264
			19HT61	KME20M0928	KME21L2132	KME22L0386
19HT66	KME20M1806	KME21L2653	KME23B1438			
19HT68	KME20M1923	KME21M0853	KME23B1830			
19IT20	KME20M2249	KME21M1380	KME23B1977			
19IT33	KME20M2450	KME21M1469	KME23B2917			
19IT59	KME20M2708	KME22A0676	KME23B3102			
19IT61	KME21A0338					

MMS #	MFG Catalog #	Description	Affected Lot(s)			
144244	112082075	TUBE, MURPHY ENDOTRACH CUFF 7.5MM (10/BX) TELFLX	18FG16	19KT41	KME21C2026	KME22A2602
			18FG29	19KT48	KME21C2051	KME22A2603
			18GG13	19KT49	KME21C2066	KME22A2898
			18GG15	19LT07	KME21C2195	KME22A2937
			18GG16	19LT23	KME21C2499	KME22A3082
			18GG22	19LT31	KME21C2615	KME22A3196
			18GG24	19LT35	KME21C2680	KME22A3219
			18GG26	19LT54	KME21C2766	KME22A3222
			18GG27	20AT11	KME21C2908	KME22A3238
			18GT16	20AT14	KME21C2964	KME22B0266
			18GT19	20AT20	KME21C3394	KME22C2126
			18GT32	20AT22	KME21D0029	KME22C2350
			18GT34	20AT26	KME21D0058	KME22C2437
			18HG14	20AT32	KME21D0186	KME22C2439
			18HG21	20AT38	KME21D0586	KME22C2604
			18HG23	20AT47	KME21D0648	KME22C2756
			18HG24	20AT54	KME21D0775	KME22C2855
			18HG27	20AT58	KME21D1103	KME22C3145
			18HG33	20BT01	KME21D1400	KME22D0324
			18HT15	20BT09	KME21D1423	KME22D0464
			18HT18	20BT24	KME21D1705	KME22D1176
			18HT22	20BT25	KME21D1780	KME22D1481
			18HT24	20BT31	KME21D1923	KME22D1660
			18HT31	20BT37	KME21D1956	KME22D1688
			18IG08	20BT39	KME21D2014	KME22D2635
			18IG22	20BT50	KME21D2018	KME22D2853
			18IG23	20BT51	KME21D2262	KME22D3116
			18IG25	20CT42	KME21E0021	KME22D3236
			18IG29	20CT43	KME21E0288	KME22D3419
			18IT04	20CT51	KME21E0469	KME22D3420
			18IT18	20CT53	KME21E0498	KME22D3711
			18IT23	20CT54	KME21E1150	KME22E0179
			18IT25	20DT01	KME21E1178	KME22E0319
			18JG01	20DT02	KME21E1186	KME22E0321
			18JG02	20DT05	KME21E1264	KME22E0322
			18JG04	20DT08	KME21F0184	KME22E0641
			18JG07	20DT20	KME21F0351	KME22E0642
			18JG10	20DT32	KME21F0385	KME22E1861
			18JG24	20DT37	KME21F0557	KME22E2010
			18JG32	20ET12	KME21F0594	KME22E2011
18JG38	20ET34	KME21F0725	KME22E2647			
18JT04	20ET37	KME21F1292	KME22E3002			
18JT06	20ET39	KME21F1380	KME22F0185			
18JT22	20FT27	KME21F1529	KME22F0317			
18KG05	20FT35	KME21F1583	KME22F0775			
18KG35	20FT36	KME21F1607	KME22F0776			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
144244	112082075	TUBE, MURPHY ENDOTRACH CUFF 7.5MM (10/BX) TELFLX	18KT04	20FT39	KME21F1696	KME22F1265
			18KT11	20FT43	KME21F1751	KME22F2028
			18KT21	20FT49	KME21F2182	KME22F2237
			18LG19	20FT53	KME21G0114	KME22F2721
			18LG38	20FT55	KME21G0169	KME22H1514
			18LG40	20FT58	KME21G0898	KME22H1855
			18LT07	20FT61	KME21G1942	KME22H1898
			19AG21	20FT71	KME21H0829	KME22H1953
			19AG36	20GG41	KME21J0200	KME22H1987
			19AT20	20GG42	KME21J0359	KME22H1993
			19BG01	20GT04	KME21J0418	KME22H2002
			19BG04	20GT10	KME21J0562	KME22H2172
			19BG06	20GT11	KME21J0647	KME22H2247
			19BG37	20GT14	KME21J0671	KME22H2300
			19BT04	20GT18	KME21J0726	KME22H2359
			19BT12	20GT25	KME21J0763	KME22H2417
			19BT18	20GT31	KME21J0854	KME22H2448
			19CG04	20GT41	KME21J1212	KME22H2607
			19CG06	20GT51	KME21J1709	KME22H2865
			19CT03	20GT57	KME21J1842	KME22J0150
			19CT04	KME20H2118	KME21J2080	KME22J0511
			19CT20	KME20H2346	KME21K0305	KME22J1107
			19CT26	KME20H2728	KME21K0986	KME22J1258
			19CT34	KME20H3115	KME21K1712	KME22J1490
			19CT42	KME20H3120	KME21K1855	KME22J3017
			19CT50	KME20J0048	KME21K2171	KME22J3143
			19CT57	KME20J1051	KME21K2475	KME22J3344
			19CT64	KME20J1455	KME21K2531	KME22K0071
			19CT69	KME20J1469	KME21K2745	KME22K0116
			19DT10	KME20J1657	KME21L0628	KME22K0138
			19DT34	KME20J1971	KME21L0880	KME22K0277
			19DT36	KME20J1997	KME21L1025	KME22K0364
			19DT42	KME20J2159	KME21L1118	KME22K0430
			19DT45	KME20J2538	KME21L1443	KME22K0580
			19DT49	KME20J2789	KME21L1581	KME22K0581
			19DT50	KME20J2836	KME21L1663	KME22L1286
			19ET02	KME20J2918	KME21L1698	KME22L1559
			19ET06	KME20J2944	KME21L1992	KME22L1716
			19ET08	KME20J2988	KME21L2222	KME22L2285
			19ET12	KME20J2992	KME21L2455	KME22M1576
19ET22	KME20J3022	KME21L2775	KME22M1645			
19ET29	KME20J3165	KME21M0193	KME22M1646			
19ET34	KME20J3284	KME21M0310	KME22M1682			
19ET36	KME20K0052	KME21M0314	KME22M1688			
19ET49	KME20K0181	KME21M0382	KME22M1717			
19ET67	KME20K1113	KME21M0500	KME22M1730			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
144244	112082075	TUBE, MURPHY ENDOTRACH CUFF 7.5MM (10/BX) TELFLX	19ET70	KME20K3044	KME21M0690	KME22M1764
			19ET72	KME20K3577	KME21M0802	KME22M1963
			19ET77	KME20L0484	KME21M0852	KME22M2037
			19ET82	KME20L0755	KME21M0930	KME22M2165
			19FT04	KME20L1092	KME21M1006	KME22M2260
			19FT25	KME20L1156	KME21M1065	KME22M2454
			19FT33	KME20L1298	KME21M1070	KME22M2521
			19GT02	KME20L1299	KME21M1085	KME22M2556
			19GT14	KME20L1704	KME21M1245	KME22M2565
			19GT16	KME20L1924	KME21M1268	KME22M2630
			19GT17	KME20L1983	KME21M1385	KME22M2637
			19GT25	KME20L2605	KME21M1503	KME22M2649
			19GT29	KME20L2798	KME21M1583	KME22M2722
			19GT30	KME20M0034	KME21M1613	KME22M2726
			19GT31	KME20M0321	KME21M1638	KME22M2773
			19GT32	KME20M1067	KME21M1749	KME22M2794
			19GT35	KME20M1507	KME21M2135	KME22M2875
			19GT40	KME20M1890	KME21M2478	KME22M2897
			19GT64	KME20M2212	KME21M2499	KME22M2908
			19HT03	KME21A3004	KME21M2646	KME22M2928
			19HT09	KME21A3259	KME21M2699	KME22M2943
			19HT15	KME21B0179	KME21M2708	KME22M3002
			19HT22	KME21B0823	KME21M2745	KME22M3017
			19HT71	KME21B2480	KME21M2790	KME22M3046
			19IT12	KME21B2482	KME21M2871	KME22M3088
			19IT13	KME21C0239	KME21M2967	KME22M3100
			19IT41	KME21C0264	KME22A0545	KME22M3142
			19IT42	KME21C0614	KME22A0854	KME22M3188
			19IT43	KME21C0615	KME22A0940	KME23A2932
			19IT45	KME21C0661	KME22A1136	KME23A3025
			19JT11	KME21C0702	KME22A1224	KME23B0026
			19JT14	KME21C0730	KME22A1367	KME23B0240
			19JT59	KME21C0826	KME22A1392	KME23B0703
			19JT61	KME21C0855	KME22A1754	KME23B1328
			19KT06	KME21C1093	KME22A1990	KME23B1570
			19KT09	KME21C1223	KME22A1993	KME23B1672
			19KT13	KME21C1323	KME22A2115	KME23B1673
			19KT14	KME21C1399	KME22A2129	KME23B1718
			19KT18	KME21C1601	KME22A2315	KME23B2081
			19KT19	KME21C1704	KME22A2351	KME23B2098
19KT20	KME21C1746	KME22A2364	KME23C1322			
19KT26	KME21C1944	KME22A2414	KME23C1441			
19KT35	KME21C2008	KME22A2537	KME23C1689			
19KT38						

MMS #	MFG Catalog #	Description	Affected Lot(s)			
144476	112082060	TUBE, MURPHY ENDOTRACH CUFF 6.0MM (10/BX) TELFLX	18FG02	19IT12	KME20L1143	KME22A1518
			18FG20	19IT33	KME20L1420	KME22A2004
			18FG24	19IT43	KME20L2089	KME22A2431
			18FT33	19IT55	KME20L2783	KME22A2484
			18GG03	19IT59	KME20M0497	KME22A2892
			18GG05	19JT30	KME20M0663	KME22A2958
			18GG06	19KT01	KME20M1064	KME22B0288
			18GG13	19KT02	KME20M1143	KME22B0499
			18GG15	19KT07	KME20M1218	KME22B0740
			18GG23	19KT08	KME21A1521	KME22B1509
			18GG25	19KT17	KME21A1803	KME22B1847
			18GT11	19KT31	KME21A2066	KME22B1895
			18HG07	19KT38	KME21A2344	KME22B2593
			18HG30	19KT50	KME21B0252	KME22C0832
			18HG33	19LT04	KME21B1151	KME22D2180
			18HG37	19LT08	KME21B1592	KME22D3181
			18HT05	20AT06	KME21B2606	KME22E0328
			18IG06	20AT41	KME21B2729	KME22E0977
			18IG14	20AT49	KME21C0483	KME22F0005
			18IG27	20BT39	KME21C1094	KME22F0189
			18IT04	20BT46	KME21C2681	KME22F0762
			18IT06	20BT51	KME21C2996	KME22F1398
			18IT10	20BT52	KME21D0703	KME22F2154
			18IT18	20CT13	KME21D1159	KME22F3471
			18IT25	20CT54	KME21D1326	KME22G0140
			18JG16	20ET12	KME21E0660	KME22G0188
			18JG24	20ET15	KME21E1391	KME22G2351
			18JT06	20ET16	KME21G0572	KME22G2656
			18KG23	20ET17	KME21G1924	KME22G3065
			18KG25	20FT23	KME21G1976	KME22G3236
			18KG33	20FT45	KME21G2250	KME22H0773
			18LG06	20FT58	KME21G2260	KME22H0828
			18LT25	20GT09	KME21H0729	KME22H1512
			19BG22	20GT34	KME21H0838	KME22H1589
			19BT04	20GT38	KME21H0911	KME22J0226
			19CG03	KME20H0123	KME21H1123	KME22J1998
			19CT25	KME20H2930	KME21J0646	KME22K0733
			19CT29	KME20J0548	KME21J1838	KME22K0898
			19CT45	KME20J1446	KME21K0411	KME22K0991
			19CT50	KME20J1457	KME21K0550	KME23B0758
			19CT61	KME20J1458	KME21L0214	KME23B0877
			19CT76	KME20J2801	KME21L0952	KME23B0878
			19DT05	KME20K0435	KME21L2133	KME23B1439
			19DT14	KME20K1112	KME21L2271	KME23B2302
			19GT07	KME20K2168	KME21M2791	KME23C0917
			19HT04	KME20K3530	KME21M2792	KME23C1809
			19HT61	KME20L0196	KME22A1454	KME23C2256
			19IT05			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
14447 7	112082090	TUBE, MURPHY ENDOTRACH CUFF 9.0MM(10/BX) TELFLX	18FG23	19AT19	KME20J2421	KME22C0133
			18FT17	19CG02	KME20K0180	KME22C0202
			18FT24	19ET11	KME20K0212	KME22C0454
			18FT25	19ET12	KME20K0600	KME22C0573
			18FT49	19JT61	KME20K0978	KME22C0675
			18GG01	19KT09	KME20K1434	KME22D3053
			18GG12	19KT19	KME20K3021	KME22D3492
			18GG26	19KT35	KME20K3525	KME22D3565
			18GT08	19LT35	KME20L0805	KME22E1281
			18HG24	19LT46	KME20L1142	KME22E1282
			18HG27	20AT57	KME21A2016	KME22E2749
			18HT07	20BT02	KME21A2676	KME22F2444
			18IG28	20BT38	KME21A3263	KME22G1323
			18IG29	20CT29	KME21B0632	KME22G1463
			18IT25	20DT04	KME21B1283	KME22H0649
			18JG38	20DT20	KME21K0162	KME22H2380
			18KG04	20DT22	KME21K0478	KME22K0996
			18KG05	20DT23	KME21L1274	KME23A2189
			18KG23	20DT28	KME21L1371	KME23A2224
			18KG26	20ET17	KME22A2199	KME23A2288
			18KG28	KME20H0415	KME22A2286	KME23A2435
			18KG38	KME20H2140	KME22A2317	KME23A2734
			18KT04	KME20H2314	KME22B1140	KME23B2546
			18KT27	KME20H2913	KME22B1997	KME23C0173
18LG13	KME20J1999	KME22B2180				
14447 9	100382030	TUBE, MURPHY ENDOTRACH UNCUFF 3.0MM (10/BX) TELFLX	18FT13	19JT21	KME20M3220	KME22E0560
			18GG16	19KT13	KME20M3221	KME22E2506
			18GT16	19LT27	KME21A0366	KME22E2684
			18GT32	20AT05	KME21A0764	KME22E2925
			18HG08	20AT37	KME21A2350	KME22F1606
			18HG25	20AT50	KME21B1132	KME22G0726
			18IG14	20BT43	KME21B1875	KME22G1285
			18IG28	20CT07	KME21B2182	KME22G3009
			18IG29	20CT12	KME21B2185	KME22H0691
			18JG35	20CT53	KME21C0890	KME22H1054
			18KG08	20DT02	KME21C1020	KME22H1223
			18KG33	20DT22	KME21C1033	KME22H1909
			18LG03	20ET30	KME21C1572	KME22H2746
			18LG23	20ET33	KME21C2686	KME22J1478
			18LG26	20ET56	KME21C2852	KME22J2117
			19AG08	20FT43	KME21D1293	KME22K0057
			19AG23	20FT55	KME21D1448	KME22K0878
			19AG38	20FT68	KME21D1743	KME22K2488
			19BG12	20GT18	KME21D1913	KME22L2920
			19BG35	20GT22	KME21E1429	KME22M0084

MMS #	MFG Catalog #	Description	Affected Lot(s)			
14447 9	100382030	TUBE, MURPHY ENDOTRACH UNCUFF 3.0MM (10/BX) TELFLX	19CT37	KME20H0952	KME21F0962	KME22M0212
			19CT42	KME20H1089	KME21F0984	KME22M0213
			19CT64	KME20H1644	KME21F1546	KME22M0265
			19DT36	KME20H2418	KME21J1148	KME22M0415
			19ET51	KME20J0795	KME21J2065	KME22M2425
			19FT08	KME20J1896	KME22A2454	KME23A0108
			19FT28	KME20J2817	KME22C1461	KME23A0308
			19GT05	KME20K0276	KME22C1763	KME23A0413
			19GT29	KME20K0325	KME22C1862	KME23A0663
			19GT35	KME20K0846	KME22D0815	KME23B0502
			19HT26	KME20K3209	KME22D1065	KME23B0537
			19HT61	KME20K3303	KME22D2425	KME23B0538
			19HT80	KME20K3401	KME22D2823	KME23B0725
			19IT11	KME20L1517	KME22D3265	KME23B1910
			19IT12	KME20L1893	KME22D3266	KME23C0161
			19IT38	KME20M1110	KME22E0277	KME23C1149
			19IT56	KME20M1899	KME22E0291	KME23C1602
19JT01						
14448 0	100382035	TUBE, MURPHY ENDOTRACH UNCUFF 3.5MM (10/BX) TELFLX	18FG14	20AT05	KME21A0354	KME22E2685
			18FT13	20AT06	KME21A0766	KME22E2686
			18GG15	20AT50	KME21A1487	KME22F1707
			18GT32	20BG19	KME21B0039	KME22F2556
			18HG08	20BT43	KME21B1553	KME22G2888
			18HG40	20CT12	KME21C0891	KME22G3193
			18IG13	20CT27	KME21C0930	KME22G3229
			18IG14	20DT01	KME21C1580	KME22H0114
			18IG28	20DT02	KME21C2242	KME22H1215
			18JG08	20DT04	KME21C2687	KME22H2298
			18KG09	20DT24	KME21C3125	KME22H2747
			18KG24	20ET30	KME21G1006	KME22J0247
			18KG25	20ET56	KME21J2129	KME22J1043
			18KG33	20FT21	KME21K0796	KME22J2124
			18KG34	20FT43	KME21K0947	KME22K0712
			18LG01	20FT55	KME21L2524	KME22K0713
			18LG03	20FT58	KME22A2459	KME22K1577
			18LG15	20FT68	KME22B0422	KME22K2171
			19AG21	20GT18	KME22B0727	KME22K2382
			19AG38	20GT22	KME22B0933	KME22K2545
			19BG12	KME20G0092	KME22B0990	KME22L1846
			19BG35	KME20H0244	KME22C0851	KME22L2147
			19CG05	KME20H0292	KME22C0852	KME22L2341
			19CT36	KME20H0780	KME22C2632	KME22L2372
			19CT40	KME20H1090	KME22C2804	KME22L2411
			19CT44	KME20H1757	KME22D0304	KME22M1269
			19CT64	KME20H1841	KME22D0880	KME22M2219

MMS #	MFG Catalog #	Description	Affected Lot(s)			
14448 0	100382035	TUBE, MURPHY ENDOTRACH UNCUFF 3.5MM (10/BX) TELFLX	19CT72	KME20H2506	KME22D1106	KME22M2426
			19DT13	KME20H2729	KME22D1895	KME22M2624
			19FT08	KME20J0673	KME22D2022	KME22M2625
			19FT12	KME20J1311	KME22D2157	KME22M2710
			19FT13	KME20J1772	KME22D2241	KME22M2900
			19GT04	KME20J2002	KME22D2242	KME22M3152
			19GT08	KME20J2096	KME22D2320	KME23A0664
			19GT29	KME20K0870	KME22D2462	KME23A0732
			19GT46	KME20K3235	KME22D2650	KME23A2851
			19HT05	KME20L0693	KME22D2810	KME23A2852
			19HT32	KME20L1581	KME22D3220	KME23A3132
			19HT61	KME20L1894	KME22D3334	KME23B0485
			19JT26	KME20M3043	KME22E0292	KME23C0162
			19KT20	KME20M3045	KME22E0669	KME23C1052
			19LT27	KME20M3223	KME22E0906	KME23C1234
14448 1	100382040	TUBE, MURPHY ENDOTRACH UNCUFF 4.0MM (10/BX) TELFLX	18FG31	19KT20	KME20L0503	KME22D0811
			18FG33	19KT37	KME20M0209	KME22D0943
			18FT13	19KT38	KME20M1518	KME22D1009
			18GG15	19LT23	KME21A1348	KME22E0606
			18HG08	20AT37	KME21B0084	KME22E1274
			18HG10	20BG19	KME21B0573	KME22E2288
			18HG25	20CT12	KME21B1133	KME22E2289
			18IG13	20CT53	KME21B1830	KME22F0571
			18IG14	20DT04	KME21C0220	KME22F1710
			18IG28	20DT24	KME21C1573	KME22G2306
			18JG08	20ET56	KME21C1675	KME22G2356
			18JG35	20FT21	KME21C2245	KME22G3111
			18KG09	20FT43	KME21C3124	KME22H0140
			18KG33	20GG42	KME21C3173	KME22H1124
			18LG01	20GT18	KME21C3313	KME22H1225
			18LG03	KME20G0126	KME21D0525	KME22H1500
			18LG05	KME20H0286	KME21D0945	KME22H1963
			18LG15	KME20H1091	KME21E1431	KME22H2405
			19BG12	KME20H2041	KME21F0921	KME22J0383
			19BG35	KME20H2419	KME21J0629	KME22J1148
			19ET09	KME20H2420	KME21K0797	KME22J1529
			19FT12	KME20J0796	KME21L2516	KME22J2122
			19FT13	KME20J1073	KME21L2526	KME22L1570
			19GT05	KME20J1844	KME21M0868	KME23A0854
			19GT35	KME20J2613	KME21M3069	KME23B1454
			19HT07	KME20J2688	KME22A2460	KME23B1756
			19HT34	KME20K0326	KME22A2706	KME23B2161
			19HT79	KME20K0601	KME22B0531	KME23C1126
19IT11	KME20K3305	KME22C2637	KME23C1289			
19JT21	KME20L0502	KME22D0032				

MMS #	MFG Catalog #	Description	Affected Lot(s)			
144482	100382050	TUBE, MURPHY ENDOTRACH UNCUFF 5.0MM (10/BX) TELFLX	18FG02	19DT39	20CT38	KME21D0736
			18FG27	19ET15	20DT13	KME21J0154
			18FT26	19ET77	20DT17	KME21K1131
			18FT35	19FT13	20DT22	KME21M2476
			18GG26	19FT19	20DT23	KME22A1524
			18GG27	19GT17	20DT32	KME22B2166
			18HG11	19GT58	20GG45	KME22D0053
			18HT18	19HT26	20GT14	KME22D0857
			18IG04	19IT05	20GT24	KME22D2267
			18IG10	19IT20	20GT32	KME22D2479
			18JG36	19IT41	20GT41	KME22D2540
			18JT06	19IT59	20GT44	KME22D3022
			18JT10	19JT18	KME20H0293	KME22E0992
			18JT21	19JT22	KME20H1490	KME22F2156
			18KG23	19JT26	KME20H2049	KME22F2426
			18KT30	19JT46	KME20H2478	KME22F2427
			18LG06	19JT49	KME20J1793	KME22G0200
			18LG19	19JT61	KME20K3053	KME22G1191
			18LG22	19KT29	KME20L1755	KME22G1749
			18LT16	19KT41	KME20L2383	KME22G1843
			18LT18	19LT46	KME20L2596	KME22H2870
			19BG03	20AT26	KME20M2056	KME22H3306
			19BG16	20AT40	KME21A0768	KME22J3090
			19BG21	20AT41	KME21D0189	KME22K0205
			19CG06	20BT08	KME21D0453	KME22K0782
			19DT16	20CT11	KME21D0647	
144483	100382060	TUBE, MURPHY ENDOTRACH UNCUFF 6.0MM (10/BX) TELFLX	18FG02	19ET09	KME20H2043	KME22A2339
			18FG28	19ET72	KME20H2479	KME22B0949
			18FT08	19GT21	KME20H2730	KME22B1778
			18FT33	19IT05	KME20H2738	KME22C2310
			18GT11	19JT14	KME20J0379	KME22C2707
			18HG20	19JT18	KME20J1224	KME22D1804
			18HG23	19JT46	KME20J3207	KME22F0417
			18HG33	20BT17	KME20K3366	KME22F0556
			18IG26	20CT07	KME21A1430	KME22F1739
			18JG15	20DT02	KME21D0667	KME22F1875
			18KG23	20DT13	KME21D0828	KME22G3211
			18KG39	20DT19	KME21F0561	KME22H0518
			18LG01	KME20H1758	KME21J0120	KME23A1071
			18LG14	KME20H2042	KME21J1583	
146251	112082085	TUBE, MURPHY ENDOTRACH CUFF 8.5MM (10/BX) TELFLX	18FG09	19LT51	KME20M1849	KME21L1369
			18FT26	20AT01	KME20M3122	KME21L1586
			18FT30	20AT06	KME20M3248	KME21M2500
			18GG02	20AT49	KME21A0697	KME21M2728
			18GG15	20BT08	KME21A1093	KME22A0787

MMS #	MFG Catalog #	Description	Affected Lot(s)			
146251	112082085	TUBE, MURPHY ENDOTRACH CUFF 8.5MM (10/BX) TELFLX	18HG15	20BT13	KME21A2265	KME22A1549
			18HG37	20BT37	KME21A2577	KME22A2287
			18HT11	20BT51	KME21B0932	KME22A2494
			18HT15	20CT02	KME21B0958	KME22A2586
			18HT24	20CT22	KME21B1215	KME22A2607
			18HT35	20CT28	KME21C3024	KME22A2957
			18IG04	20CT32	KME21C3324	KME22A2968
			18IG10	20CT42	KME21D0061	KME22B0518
			18IG25	20DT02	KME21E0902	KME22B0519
			18IG27	20DT17	KME21E1291	KME22B0750
			18IG29	20DT18	KME21F0741	KME22B0844
			18JG12	20ET12	KME21F1214	KME22B0899
			18JG15	20ET16	KME21F1531	KME22B0942
			18JG18	20ET17	KME21F2006	KME22B0993
			18JG26	20GG42	KME21F2074	KME22B1129
			18JG31	20GG43	KME21H0920	KME22B1154
			18JG33	20GT34	KME21H1219	KME22B1343
			18JG34	20GT40	KME21J0417	KME22B1696
			18JT12	20GT50	KME21J0485	KME22B1775
			19BG18	KME20H0750	KME21J0640	KME22B2179
			19BT08	KME20J2539	KME21J1304	KME22B2320
			19CT69	KME20J2680	KME21J1481	KME22B2379
			19CT73	KME20J2999	KME21J1609	KME22B2477
			19CT74	KME20K0923	KME21J1634	KME22B2797
			19CT76	KME20K1114	KME21J1712	KME22C0572
			19CT77	KME20K1273	KME21J1755	KME22C2074
			19DT05	KME20K1513	KME21J1864	KME22D1891
			19DT10	KME20K3020	KME21J1872	KME22D1892
			19ET08	KME20K3034	KME21J2082	KME22D3186
			19KT09	KME20K3146	KME21K0095	KME22D3488
			19KT14	KME20L0670	KME21K0279	KME22D3564
			19KT30	KME20L0756	KME21K0412	KME22E0663
			19KT31	KME20L0757	KME21K1193	KME22E0988
			19KT38	KME20L1437	KME21K2275	KME22E1279
			19KT46	KME20L1549	KME21K2575	KME22E1583
			19KT49	KME20L1678	KME21K2839	KME22E1811
			19KT50	KME20L1707	KME21K3015	KME22E2253
			19KT55	KME20L1867	KME21L0394	KME22F1373
			19LT04	KME20L2274	KME21L0571	KME22H0102
			19LT35	KME20L2831	KME21L0799	KME22H1602
			19LT38	KME20M1511	KME21L1115	KME22J2933
			19LT47	KME20M1548	KME21L1276	KME23B1630
146252	112082055	TUBE, MURPHY ENDOTRACH CUFF 5.5MM (10/BX) TELFLX	18FG20	19AG07	KME20K0464	KME22D0603
			18FG33	19AG19	KME20K0683	KME22D1799
			18FT24	19AG32	KME20L1703	KME22D3018

MMS #	MFG Catalog #	Description	Affected Lot(s)			
146251	112082085	TUBE, MURPHY ENDOTRACH CUFF 8.5MM (10/BX) TELFLX	18FT26	19AG33	KME20M0996	KME22E0652
			18GG02	19BG02	KME21A0492	KME22E1164
			18GG13	19BG10	KME21A1433	KME22E1423
			18GG15	19BG37	KME21C0502	KME22E2359
			18GG25	19ET31	KME21C1347	KME22F0277
			18GT11	19ET73	KME21C1496	KME22F0418
			18GT34	19GT25	KME21C2745	KME22F0898
			18HG24	19GT58	KME21D0835	KME22F1662
			18HG26	19JT13	KME21E0864	KME22G0318
			18HG28	19JT36	KME21E1230	KME22G2350
			18HG33	19JT46	KME21F0946	KME22H0635
			18HG34	19KT02	KME21F1967	KME22H1516
			18IG05	19KT14	KME21G0850	KME22H1612
			18IG14	19KT18	KME21J0319	KME22J1044
			18IG19	19KT20	KME21K0161	KME22J1464
			18IG22	19KT21	KME21K1627	KME22J2331
			18IG23	20AT32	KME21M0312	KME22J2332
			18IG26	20AT40	KME21M2559	KME22K0403
			18JT04	20BT01	KME22A0149	KME22K0534
			18KG10	20BT10	KME22A2333	KME22K0635
			18KG25	20BT13	KME22A2942	KME22K0780
			18KG26	20BT38	KME22B0352	KME22K1440
			18KG28	20BT51	KME22B0563	KME22L1218
			18KG33	20CT54	KME22B1357	KME23A2300
			18KG37	20ET17	KME22C0323	KME23A2432
			18KT22	20ET18	KME22C1515	KME23A2845
			18KT31	KME20G0395	KME22C2218	KME23A3010
			18LG19	KME20H0640	KME22C2923	KME23B0428
			18LG40	KME20J2395	KME22D0055	KME23C0594
			19AG02	KME20J2829	KME22D0282	KME23C1669
			147586	100382025	TUBE, MURPHY ENDOTRACH UNCUFF 2.5MM (10/BX) TELFLX	18FG22
18FG31	19GT04	KME20L1514				KME22F2559
18FT13	19GT29	KME20L2483				KME22F2561
18GG16	19GT31	KME20M3046				KME22G0176
18GG26	19HT32	KME21A2200				KME22G0436
18GG34	19IT11	KME21A3288				KME22G1235
18GT12	19IT56	KME21B0567				KME22H0629
18HG08	19IT58	KME21C1019				KME22H0630
18IG13	19JT51	KME21C1802				KME22H1372
18IG14	19KT08	KME21C2684				KME22H1542
18IG28	19KT65	KME21D0352				KME22H1766
18IG29	19LT03	KME21D0601				KME22K0062
18IT06	19LT27	KME21D1281				KME22K0094
18JG09	20AT06	KME21D1818				KME22K0425
18JG14	20AT37	KME21E1427				KME22K2437

MMS #	MFG Catalog #	Description	Affected Lot(s)			
147586	100382025	TUBE, MURPHY ENDOTRACH UNCUFF 2.5MM (10/BX) TELFLX	18JG35	20BT15	KME21K0765	KME22K2588
			18KG08	20CT12	KME21M1976	KME22L2985
			18KG24	20CT54	KME21M3021	KME22M0174
			18KG25	20DT02	KME21M3076	KME22M0558
			18LG05	20DT04	KME21M3077	KME22M0560
			18LG12	20DT22	KME22A2458	KME22M1206
			18LG15	20ET03	KME22B0930	KME22M1625
			19AG04	20ET33	KME22B0932	KME23A0661
			19AG07	20ET56	KME22B1025	KME23A1137
			19AG21	20FT21	KME22B1089	KME23A1600
			19AG38	20FT40	KME22C0099	KME23A2103
			19BG12	20FT43	KME22C3367	KME23A2228
			19CT37	20GG41	KME22D0179	KME23A2314
			19CT41	20GT37	KME22D2727	KME23B1565
			19CT44	KME20G0261	KME22D2824	KME23B1728
			19CT64	KME20H2887	KME22D3221	KME23B2474
			19ET09	KME20J1072	KME22E0297	KME23C0307
			19ET53	KME20J1895	KME22E2283	KME23C1051
			19ET57	KME20K2027	KME22E2285	KME23C1276
			19FT08	KME20K3208	KME22F0538	KME23C1709
19FT25	KME20K3400					
149619	112082050	TUBE, MURPHY ENDOTRACH CUFF 5.0MM (10/BX) TELFLX	18FG24	19BG26	20GG47	KME22B2692
			18FT39	19BG42	20GT32	KME22C0066
			18GG01	19BT04	KME20H0638	KME22C0571
			18GG02	19CT03	KME20H2725	KME22C3132
			18GG19	19CT25	KME20J0260	KME22D0281
			18GG23	19CT32	KME20J0421	KME22D2010
			18GT06	19CT82	KME20J1385	KME22D2503
			18HG30	19DT05	KME20J2065	KME22D2504
			18IG18	19FT38	KME20K0211	KME22D2648
			18IG23	19GT07	KME20K0463	KME22E0316
			18IG25	19HT26	KME21A2017	KME22E1277
			18IG30	19HT79	KME21A2343	KME22E1422
			18IT12	19IT26	KME21A2686	KME22E2743
			18IT18	19JT18	KME21B2491	KME22E2879
			18JG04	19KT17	KME21F0321	KME22G0063
			18JG10	19KT20	KME21F0727	KME22J0616
			18JG15	19LT07	KME21F2114	KME22J1331
			18JG16	19LT13	KME21F2115	KME22J1463
			18JG18	19LT17	KME21G1704	KME22K0660
			18JG31	19LT50	KME21G1949	KME22K1433
18JT06	20AT13	KME21J0561	KME22L1657			
18JT16	20AT40	KME21K0151	KME23A2133			
18KG33	20BT37	KME21K0292	KME23A2145			
18KG35	20BT38	KME21M2075	KME23A2148			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
149619	112082050	TUBE, MURPHY ENDOTRACH CUFF 5.0MM (10/BX) TELFLX	18KG37	20CT02	KME21M2306	KME23A2214
			18KT22	20CT12	KME21M2649	KME23A2221
			18LG06	20CT18	KME22A1589	KME23A2299
			19AG33	20CT37	KME22A2197	KME23A2358
			19AT20	20CT43	KME22A2604	KME23A3027
			19BG02	20DT13	KME22B1848	KME23A3166
			19BG03	20ET12	KME22B2108	KME23C0172
			19BG10	20ET17	KME22B2478	KME23C0323
			19BG15	20ET56	KME22B2538	KME23C0713
			19BG22	20FT23	KME22B2578	KME23C1115
149620	112082100	TUBE, ENDOTRACH NASORAL 10.0MM(10/BX) TELFLX	18GG13	20AT01	KME21A2668	KME22G2264
			18JG24	20AT22	KME21C0017	KME22H0209
			18JG36	20BG18	KME21C0068	KME22H2251
			18JT06	20DT18	KME21D0827	KME22J2931
			18KG10	20DT19	KME21D1219	KME22J3329
			18LG06	20ET17	KME21F2252	KME22L0218
			19AG36	KME20H2141	KME21F2345	KME22L0219
			19CG02	KME20K0977	KME22B1777	KME22L2730
			19GT16	KME20K1809	KME22B1881	KME22L2731
			19GT59	KME20K2553	KME22C0019	KME23A2196
			19IT33	KME20L0335	KME22C0063	KME23A2539
			19JT56	KME20M1341	KME22C0327	KME23A3169
			19KT14	KME20M2113	KME22E2679	KME23C2078
19KT49	KME21A2342	KME22F0401				
152534	112082095	TUBE, ENDOTRACH NASORAL 9.5MM (10/BX) TELFLX	18HG08	19LT08	KME20L2462	KME22C2606
			18HT04	20AT01	KME20L2557	KME22C2757
			18IG29	20AT13	KME20M2115	KME22D1593
			18IT04	20CT02	KME20M2198	KME22D1778
			18JG36	20CT07	KME21A2346	KME22E1285
			18JT16	20CT29	KME21A3264	KME22E2362
			18KG37	20DT18	KME21B1895	KME22F1401
			18KT27	20DT32	KME21C3326	KME22F2452
			18LT16	20GT23	KME21E1393	KME22F2784
			19BG02	20GT44	KME21J0116	KME22J1579
			19CG02	KME20H1486	KME21K1818	KME23A0557
			19CT04	KME20J1899	KME22B1259	KME23A0763
			19CT72	KME20K1583	KME22B1420	KME23A1256
			19GT25	KME20K1831	KME22B1854	KME23B0010
			19HT19	KME20K2166	KME22C0677	KME23B0168
			19HT47	KME20K2574	KME22C0713	KME23B0169
			19JT05	KME20L1041	KME22C0844	KME23B1752
19KT55	KME20L1592	KME22C1612	KME23B1964			
159586	100382055	TUBE, MURPHY ENDOTRACH UNCUFF 5.5MM (10/BX) TELFLX	18FG21	18KG23	19JT49	KME21J0727
			18FG27	18LE50	20CT47	KME21J0885
			18FG33	18LG35	20DT01	KME21J1363

MMS #	MFG Catalog #	Description	Affected Lot(s)			
159586	100382055	TUBE, MURPHY ENDOTRACH UNCUFF 5.5MM (10/BX) TELFLX	18FT17	18LT04	20GG45	KME22B2278
			18FT31	18LT16	20GT40	KME22C2285
			18FT49	19AG05	KME20H0294	KME22C2351
			18GT06	19AT23	KME20H0885	KME22E0326
			18HG06	19BG10	KME20H3036	KME22F2548
			18HG15	19BG21	KME20J0600	KME22H0306
			18HG23	19BG26	KME20J2166	KME22H0308
			18HG30	19BT08	KME20K1116	KME22H1209
			18HT04	19ET72	KME20K1206	KME22H1373
			18HT07	19FT29	KME20K1272	KME22H1783
			18IT06	19HT26	KME20K3277	KME22H2171
			18JG14	19IT42	KME20L1040	KME22H3008
			18JT12	19JT16	KME20L1224	KME22J0513
			18JT21	19JT18	KME21A0136	KME22J1254
			18KG10	19JT19	KME21H0922	KME22J2057
18KG11	19JT25	KME21J0638				
159587	100382045	TUBE, MURPHY ENDOTRACH UNCUFF 4.5MM (10/BX) TELFLX	18GG15	19HT32	KME20J3084	KME22D2158
			18GG26	19HT61	KME20K0469	KME22D2478
			18GG34	19HT80	KME20K2028	KME22E0295
			18HG10	19IT11	KME20K3205	KME22E0608
			18HG25	19IT58	KME20K3307	KME22E0909
			18IG04	19IT59	KME20L1508	KME22E0911
			18IG14	19JT01	KME20M0210	KME22E1313
			18IG28	19KT13	KME20M1075	KME22E2291
			18IG29	19LT50	KME20M3225	KME22F3304
			18JG08	20AT50	KME21A1488	KME22G2358
			18JG35	20BG15	KME21A2309	KME22H2502
			18KG08	20BT43	KME21B0085	KME22H3167
			18KG33	20BT51	KME21B1134	KME22J3174
			18LG01	20CT11	KME21B1828	KME22K2663
			18LG12	20DT04	KME21C1577	KME22K2686
			18LG30	20DT22	KME21C1677	KME22M0003
			19BG12	20FT43	KME21C2244	KME22M1623
			19CG05	20GT22	KME21C2444	KME23A2310
			19CT16	20GT37	KME21D2950	KME23B0782
			19DT13	KME20H0783	KME21D2951	KME23B1757
19ET09	KME20J0674	KME21G1008	KME23B2163			
19FT08	KME20J2818	KME21G1167	KME23C1119			
19FT28	KME20J2819	KME22C2577	KME23C1128			
19GT02	KME20J2846	KME22D2082	KME23C1336			
19GT35						
173546	100380055	TUBE, MAGILL ENDOTRACH UNCUFF 5.5MM (10/BX) TELFLX	19BG37	KME20H2315	KME21G2361	KME22F2818
			19ET11	KME21A0571	KME22A0155	KME22H2362
			19GT51	KME21A2778	KME22A2335	KME22K0054
			19JT20	KME21C2155	KME22E0288	KME22K1481

MMS #	MFG Catalog #	Description	Affected Lot(s)			
173546	100380055	TUBE, MAGILL ENDOTRACH UNCUFF 5.5MM (10/BX) TELFLX	19LT51	KME21C2572	KME22E1272	KME22M2009
			20AT13	KME21G1496	KME22F1551	KME23C3069
			20ET07			
177517	100380035	TUBE, MAGILL ENDOTRACH UNCUFF 3.5MM (10/BX) TELFLX	18HG08	19FT13	KME20L1365	KME22G3264
			18IG13	19HT07	KME20M3212	KME22H0689
			18IG14	19HT79	KME21C1679	KME22H1774
			18JG35	19LT28	KME21K0768	KME22J0819
			18KG24	20ET03	KME21M2863	KME22J2899
			18LG31	20GT37	KME22D0945	KME22J2900
			19BG15	KME20H2590	KME22E1145	KME22K2123
			19CT16	KME20L0690	KME22G2935	KME23A3131
189776	100380025	TUBE, MAGILL ENDOTRACH UNCUFF 2.5MM (10/BX) TELFLX	19ET31			
			19ET51	19LT28	KME21D0327	KME22H0686
			19FT29	20AT05	KME21F0137	KME22H1127
			19GT05	20AT37	KME21K0799	KME22H1360
			19GT35	20AT54	KME21K0800	KME22J2603
			19HT05	20BG15	KME22A2704	KME22K0674
			19HT32	20DT24	KME22C0981	KME23A1837
			19HT79	20ET03	KME22C1568	KME23B2585
			19JT51	KME21B0529	KME22D0029	KME23C0306
189777	100380030	TUBE, MAGILL ENDOTRACH UNCUFF 3.0MM (10/BX) TELFLX	19KT06	KME21C0887		
			19ET31	19HT32	20FT68	KME21M2737
			19FT28	19HT60	20GT37	KME22G1609
			19GT08	19IT38	KME20J3083	KME22H1076
			19GT31	20DT09	KME21A0758	KME22J1610
			19GT35	20DT24	KME21C2233	KME22J2700
			19GT48	20ET03	KME21J1139	KME23B1755
189778	100380040	TUBE, MAGILL ENDOTRACH UNCUFF 4.0MM (10/BX) TELFLX	19HT05	20FT55	KME21L0959	
			19FT28	20BG15	KME21A0067	KME22B2857
			19GT04	20ET03	KME21A0858	KME22C3373
			19GT51	KME20H2609	KME21C1765	KME22D0876
			19HT32	KME20J3088	KME21K0893	KME22G2187
			19HT61	KME20K2749	KME21K1264	KME22L2330
			19JT26	KME20L0640	KME21L0962	KME23A1889
			19KT13	KME20M1726	KME21M2918	KME23B2199
189780	100380050	TUBE, MAGILL ENDOTRACH UNCUFF 5.0MM (10/BX) TELFLX	19LT28	KME20M2538		
			18FG05	19HT39	20GG27	KME22A0422
			18FT26	19IT20	KME20H2359	KME22D2223
			18GT27	19KT21	KME20K1713	KME22D3339
			18IG06	19LT51	KME20L1811	KME22E0284
			18JG20	20AT26	KME20M2219	KME22F3118
			18LG03	20AT41	KME21A0129	KME22G2639
			18LT24	20CT43	KME21A2253	KME22H2816
			19BG30	20DT13	KME21B0984	KME22J3320
			19CT50	20DT32	KME21K1153	KME23A1612
	19DT17	20ET07	KME21L0988			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
189781	100380060	TUBE, MAGILL ENDOTRACH UNCUFF 6.0MM (10/BX) TELFLX	18FG31	19ET72	KME21A0131	KME22D1803
			18HG33	19JT30	KME21D0600	KME22G3107
			19AG09	20AT13	KME21H0921	KME22H0501
			19BT11	20BT13	KME21J0790	KME22K0061
			19CT26	20DT13	KME22A0448	KME22K1283
189782	100380065	TUBE, MAGILL ENDOTRACH UNCUFF 6.5MM (10/BX)	19CT32	20AT20	KME22A2337	KME22H2352
			19GT16	KME21D0517	KME22G0433	KME22K0063
			19KT41	KME21K2919	KME22H0503	KME22K1651
189783	100382065	TUBE, MURPHY ENDOTRACH UNCUFF 6.5MM (10/BX) TELFLX	18GG26	19BT19	19LT13	KME22C0085
			18HG16	19CT22	KME20J0549	KME22C2640
			18HG28	19CT39	KME20L1582	KME22C2849
			18KG25	19ET06	KME20L1583	KME22E0482
			18LG08	19ET77	KME21A1615	KME22E2541
			18LT04	19GT16	KME21D2952	KME22F1554
			19BG10	19IT59	KME21L2054	KME22F2582
			19BG16	19KT02	KME22A2341	
191305	112080070	TUBE, MAGILL ENDOTRACH NASAL 7.0MM (10/BX) TELFLX	18FG19	19HT26	20ET12	KME21B2490
			18GG19	19HT69	20ET17	KME21C1184
			18HG28	19IT11	20FT11	KME21F1585
			18IG06	19JT21	20FT19	KME21J0644
			18IG19	19KT14	20GG40	KME21J1078
			18JG13	19LT31	20GT38	KME22A1109
			18JG17	19LT51	KME20K2577	KME22A1781
			18KG13	20AT22	KME20L1985	KME22B0401
			18LT07	20AT58	KME20M1771	KME22B0402
			19BG19	20BT13	KME20M1772	KME22F1174
			19CT26	20BT42	KME20M2856	KME22F2171
			19DT20	20CT13	KME21A2355	KME22G1295
			19ET12	20DT10	KME21B1624	KME22G2999
			19ET36			
191306	112080075	TUBE, MAGILL ENDOTRACH NASAL 7.5MM (10/BX) TELFLX	18FT26	19AT23	20BT42	KME21F1534
			18GG17	19CT64	20CT02	KME22A2269
			18GG19	19ET12	20CT54	KME22B0741
			18HG15	19ET24	20DT09	KME22B0896
			18HG37	19GT31	20FT21	KME22B1031
			18IG04	19HT69	20FT60	KME22C0096
			18IG29	19IT11	20GT18	KME22E2348
			18JG17	20AT22	KME20H0334	KME22F2118
			18KG10	20BT13	KME20K0110	KME22J0567
191307	112080090	TUBE, MAGILL ENDOTRACH NASAL 9.0MM (10/BX) TELFLX	18FT24	19GT29	19KT58	KME20M2248
			18FT26	19HT79	KME20K1808	KME22E2355
			19ET83	19KT38	KME20L2352	KME22H0310
191450	112080060	TUBE, MAGILL ENDOTRACH NASAL 6.0MM (10/BX) TELFLX	18GG26	19ET31	20CT13	KME21K0521
			18GT08	19GT25	20FT66	KME22A2893
			18GT11	19IT20	KME20J1463	KME22B2215

MMS #	MFG Catalog #	Description	Affected Lot(s)			
191450	112080060	TUBE, MAGILL ENDOTRACH NASAL 6.0MM (10/BX) TELFLX	18IG22	19KT14	KME20M0616	KME22C2044
			18IT18	19KT47	KME20M2112	KME22G2998
			18JG33	20AT58	KME21A2685	KME22J0859
			19AG39			
204182	112080065	TUBE, MAGILL ENDOTRACH NASAL 8.5MM (10/BX) TELFLX	18GT07	19CT64	20BT13	KME22A2546
			18HG20	19ET12	20DT40	KME22B2739
			18IG16	19ET31	20GT05	KME22C2271
			18IG22	19HT69	KME20K2581	KME22E0189
			18JG13	19LT47	KME21C1701	KME22F1550
			19AT23	20AT41	KME22A2193	
204183	112080085	TUBE, MAGILL ENDOTRACH NASAL 8.5MM (10/BX) TELFLX	18IG19	19KT38	20ET17	KME21C0067
			18KG26	19LT17	20GT57	KME22A1989
			19CT61	20CT02	KME20M0909	KME22E2141
			19GT31	20DT23	KME21B0519	KME22G2263
			19JT19			
229777	111781035	TUBE, MURPHY ENDOTRACH CUFF 3.5MM (10/BX) TELFLX	18FG20	20CT11	KME22A0560	KME22G0600
			18FG31	20CT32	KME22A0755	KME22G0898
			18GG23	20DT19	KME22A1376	KME22G1642
			18GG34	20DT32	KME22C0493	KME22G1861
			18GT32	20ET34	KME22D1697	KME22G2270
			18HG17	20ET49	KME22D2626	KME22G3150
			18IG15	20FT18	KME22D2999	KME22G3268
			18JG01	20FT64	KME22D3173	KME22H2027
			18JG16	20GT11	KME22D3253	KME22H3065
			18JG35	KME20H0466	KME22E0727	KME22H3272
			18KG17	KME20H2769	KME22E0728	KME22H3273
			18KG18	KME20H2998	KME22E1011	KME22H3400
			18KG36	KME20H3129	KME22E1155	KME22J0633
			18LG15	KME20M0650	KME22E1392	KME22J0901
			18LG30	KME20M1283	KME22E2151	KME22J1975
			18LG39	KME20M2799	KME22E2239	KME22J2073
			19AG15	KME20M2822	KME22E2744	KME22K1778
			19BG33	KME21A2028	KME22E2753	KME22K1842
			19CG01	KME21A3020	KME22F0075	KME22K1924
			19CG02	KME21B1044	KME22F0229	KME22K2089
			19CT26	KME21B1197	KME22F0540	KME22L0809
			19DT03	KME21B1960	KME22F0540	KME22L1004
			19DT30	KME21B1961	KME22F0600	KME22L1046
			19DT38	KME21C0324	KME22F0880	KME22L1130
			19ET58	KME21C0426	KME22F1315	KME22L1463
			19ET69	KME21C1119	KME22F1736	KME23A1767
			19FT37	KME21L0024	KME22F2068	KME23A2684
			19FT46	KME21L0026	KME22F2353	KME23B0204
			19IT53	KME21L1782	KME22F2627	KME23B0380
			19JT10	KME21L2301	KME22F2949	KME23B1912

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229777	111781035	TUBE, MURPHY ENDOTRACH CUFF 3.5MM (10/BX) TELFLX	19KT14	KME21L2785	KME22F3119	KME23B3058
			19KT46	KME21M1710	KME22F3242	KME23B3059
			19LT16	KME21M2600	KME22F3425	KME23C0339
			20AT22	KME21M2601	KME22G0214	KME23C0586
			20AT26			
229778	111781040	TUBE, NASAL MURPHY ENDOTRACH CUFF 40MM (10/BX) TELFLX	18FT35	19JT63	KME21K1714	KME22E0945
			18GG16	19KT02	KME21K1841	KME22E1105
			18GG19	19KT41	KME21K2300	KME22E1394
			18GG23	19KT46	KME21K3006	KME22E1499
			18GG34	19LT21	KME21L2290	KME22E1933
			18GT32	19LT36	KME21L2744	KME22E2376
			18HG17	20AT22	KME21L2745	KME22E2659
			18HG21	20AT53	KME21L2760	KME22E2810
			18HG22	20BT01	KME21L2761	KME22E2889
			18HG29	20BT02	KME21L2787	KME22F0001
			18HG30	20CT05	KME21L2869	KME22F0541
			18IG15	20DT23	KME21M0203	KME22F1109
			18JG01	20DT29	KME21M0801	KME22F2070
			18JG03	20DT32	KME21M0809	KME22F2950
			18JG16	20DT35	KME21M0950	KME22F2997
			18JG19	20ET04	KME21M1706	KME22F3120
			18JG35	20ET33	KME21M1904	KME22F3243
			18JG37	20ET34	KME21M1905	KME22F3423
			18KG17	20ET38	KME21M2908	KME22G0815
			18KG19	20ET39	KME22A0561	KME22G1098
			18KG38	20FT17	KME22A2283	KME22G1167
			18LG12	20FT19	KME22A2404	KME22G1228
			18LG14	20FT60	KME22A2849	KME22G1229
			18LG18	20FT64	KME22A2915	KME22G1330
			18LG30	20FT68	KME22A3199	KME22G2240
			19AG01	20GT11	KME22B0002	KME22G2307
			19AG19	20GT34	KME22B0003	KME22G2576
			19AG29	KME20H1786	KME22B0004	KME22G2629
			19BG18	KME20H2999	KME22B0043	KME22G2924
			19BG32	KME20J0367	KME22B0123	KME22G3134
			19BG33	KME20J0494	KME22B0124	KME22H0036
			19CG01	KME20J1209	KME22B0165	KME22H0098
			19CG02	KME20J1388	KME22B0192	KME22H0545
			19CT26	KME20M2796	KME22B0545	KME22H0600
19CT33	KME20M2906	KME22B0585	KME22H0761			
19CT53	KME20M2907	KME22B1999	KME22H1361			
19CT73	KME20M2994	KME22C1877	KME22H1725			
19CT80	KME20M3092	KME22C1989	KME22H1932			
19CT82	KME21A1005	KME22C1990	KME22H2056			
19DT05	KME21A1006	KME22C2189	KME22H2520			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229778	111781040	TUBE, NASAL MURPHY ENDOTRACH CUFF 40MM (10/BX) TELFLX	19DT30	KME21A2030	KME22C2311	KME22H3068
			19DT35	KME21A2037	KME22C2479	KME22J0197
			19DT38	KME21A3023	KME22C2480	KME22J0848
			19DT42	KME21A3024	KME22C2481	KME22J2169
			19ET18	KME21A3028	KME22C2659	KME22J2411
			19ET34	KME21A3029	KME22C2660	KME22J2625
			19ET39	KME21B0259	KME22C2760	KME22J2840
			19ET40	KME21B0788	KME22C2839	KME22J2964
			19ET58	KME21B0994	KME22C2840	KME22K1626
			19ET60	KME21B1198	KME22C2866	KME22K1715
			19ET61	KME21B1286	KME22C2966	KME22K1779
			19ET65	KME21B1288	KME22C3124	KME22K1850
			19ET69	KME21B1620	KME22C3126	KME22K1851
			19FT01	KME21B2627	KME22C3454	KME22K2013
			19FT41	KME21B2628	KME22D0284	KME22K2507
			19GT13	KME21B2705	KME22D1566	KME22K2510
			19GT16	KME21C0428	KME22D1674	KME23A1733
			19GT63	KME21C1121	KME22D1698	KME23A2846
			19HT31	KME21C1199	KME22D1813	KME23A2847
			19HT82	KME21C1879	KME22D1890	KME23A2941
			19IT01	KME21C3336	KME22D2708	KME23B0383
			19IT14	KME21F0980	KME22D3011	KME23B0384
			19JT35	KME21K1094	KME22E0730	KME23B0882
			19JT56	KME21K1655	KME22E0731	
229779	111781045	TUBE, MURPHY ENDOTRACH CUFF 4.5MM (10/BX) TELFLX	18FG33	19KT14	KME21K2459	KME22C1991
			18GG16	19KT38	KME21K2495	KME22C1992
			18HG17	19KT41	KME21K2579	KME22C2190
			18HG21	19LT16	KME21K2998	KME22C2968
			18HG29	19LT21	KME21K3000	KME22D0091
			18HG30	19LT36	KME21K3001	KME22D0153
			18HG31	20AT22	KME21L0074	KME22D0248
			18IG10	20BT02	KME21L0890	KME22D0285
			18IG14	20CT02	KME21L1033	KME22D0321
			18IG17	20CT05	KME21L1359	KME22D0442
			18IG20	20CT11	KME21L1364	KME22D0444
			18JG02	20CT14	KME21L1365	KME22D0536
			18JG03	20CT32	KME21L1435	KME22D0717
			18JG16	20DT19	KME21L1525	KME22D1659
			18JG19	20DT20	KME21L2143	KME22D2057
			18JG20	20DT23	KME21L2788	KME22D2469
			18JG35	20DT29	KME21M0949	KME22D3211
			18JG37	20DT32	KME21M1654	KME22E1501
			18KG01	20DT35	KME21M1655	KME22E1932
			18KG17	20ET33	KME21M1656	KME22E2026
			18KG19	20ET34	KME21M1906	KME22E2122

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229779	111781045	TUBE, MURPHY ENDOTRACH CUFF 4.5MM (10/BX) TELFLX	18KG36	20ET38	KME21M1907	KME22E2690
			18KG39	20ET39	KME21M2271	KME22E2805
			18KG41	20ET49	KME21M2424	KME22E2888
			18LG12	20FT17	KME21M2479	KME22E2928
			18LG14	20FT18	KME21M2674	KME22E3030
			18LG30	20FT60	KME21M2706	KME22F0230
			18LG39	20GT34	KME21M2909	KME22F2072
			19AG01	KME20H0612	KME22A0562	KME22F2355
			19AG17	KME20H1787	KME22A1471	KME22F2625
			19AG19	KME20H2570	KME22A2179	KME22F3383
			19AG29	KME20H2572	KME22A2917	KME22F3502
			19BG02	KME20H2573	KME22A3198	KME22G0215
			19BG18	KME20H3001	KME22A3264	KME22G0234
			19BG19	KME20H3127	KME22B0005	KME22G0262
			19BG32	KME20H3128	KME22B0167	KME22G0476
			19BG33	KME20J1211	KME22B0191	KME22G0537
			19CG01	KME20J1389	KME22B0239	KME22G0688
			19CG02	KME21A1007	KME22B0277	KME22G0735
			19CG04	KME21A1008	KME22B0358	KME22G0788
			19CT26	KME21A1182	KME22B0399	KME22G1047
			19CT33	KME21A2025	KME22B0502	KME22G1109
			19CT53	KME21A2026	KME22B0580	KME22G1332
			19CT56	KME21A2189	KME22B1038	KME22G1466
			19CT80	KME21A3059	KME22B1144	KME22G2049
			19DT05	KME21A3060	KME22B1145	KME22G2244
			19DT30	KME21A3083	KME22B1224	KME22G2774
			19DT46	KME21A3084	KME22B1265	KME22H1040
			19ET18	KME21B0258	KME22B1309	KME22H1202
			19ET34	KME21B0790	KME22B1341	KME22H1424
			19ET39	KME21B1083	KME22B1418	KME22H1619
			19ET40	KME21B1259	KME22C0060	KME22H1648
			19ET61	KME21B1289	KME22C0494	KME22H2057
			19ET65	KME21B1545	KME22C0791	KME22H2522
			19ET69	KME21B1605	KME22C0792	KME22H2724
			19FT01	KME21B1754	KME22C1080	KME22H3066
			19FT41	KME21B2217	KME22C1081	KME22J0547
			19GT13	KME21B2530	KME22C1082	KME22J1265
			19HT68	KME21B2672	KME22C1236	KME22J2026
			19IT01	KME21C0424	KME22C1237	KME22K0152
			19IT15	KME21C0425	KME22C1472	KME22K0249
			19IT53	KME21C2416	KME22C1710	KME22K0340
			19IT54	KME21J1495	KME22C1711	KME22K0559
			19JT06	KME21K1633	KME22C1724	KME22K0703
			19JT10	KME21K1715	KME22C1850	KME22K0721
			19JT14	KME21K1842	KME22C1852	KME22K0752
			19JT56	KME21K1866	KME22C1878	KME22M2486
			19JT63	KME21K2302	KME22C1879	KME23A0494
			19KT02	KME21K2458	KME22C1880	KME23B1758

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229780	111781050	TUBE, MURPHY ENDOTRACH CUFF 5.0MM (10/BX) TELFLX	18FG31	19LT21	KME21K1079	KME22D0805
			18FT35	20AT22	KME21K1716	KME22D0854
			18FT39	20BT53	KME21K1826	KME22D0936
			18GG17	20CT07	KME21K1868	KME22D1008
			18GG18	20CT12	KME21K2303	KME22D1453
			18GT32	20CT14	KME21K2335	KME22D1957
			18HG17	20CT32	KME21K3007	KME22D2058
			18HG19	20DT18	KME21L0889	KME22D2293
			18HG21	20DT29	KME21L0985	KME22D2643
			18HG22	20ET34	KME21L1526	KME22D2765
			18HG29	20ET38	KME21L1597	KME22D3013
			18IG15	20ET49	KME21L1776	KME22E0947
			18IG17	20FT17	KME21L2153	KME22E0956
			18IG30	20FT19	KME21L2154	KME22E1500
			18JG16	20FT21	KME21L2299	KME22E1885
			18JG19	20FT64	KME21L2885	KME22E2219
			18JG35	20GT11	KME21M0162	KME22E2220
			18JG37	KME20J0370	KME21M0337	KME22E2440
			18KG17	KME20J0663	KME21M0342	KME22E2512
			18KG39	KME20J0665	KME21M0502	KME22E2872
			18KG40	KME20J1212	KME21M0733	KME22F0074
			18LG14	KME20K3628	KME21M0951	KME22F0542
			18LG15	KME20L2213	KME21M1978	KME22F0987
			18LG18	KME20L2286	KME21M1979	KME22F1112
			18LG39	KME20L2666	KME21M2427	KME22F2173
			19AG25	KME20L2668	KME21M2675	KME22F2206
			19BG02	KME20M0644	KME21M2676	KME22F2951
			19BG04	KME20M1284	KME21M2711	KME22G0207
			19BG18	KME20M2643	KME21M2807	KME22G0380
			19BG19	KME20M2645	KME22A0564	KME22G0831
			19BG32	KME20M2956	KME22A0565	KME22G1050
			19CG01	KME20M3015	KME22A1888	KME22G1227
			19CG02	KME21A1024	KME22A2257	KME22G1267
			19CT28	KME21A1065	KME22A2369	KME22G1467
			19CT33	KME21A1273	KME22A2370	KME22G1643
			19CT36	KME21A2079	KME22A2850	KME22G2688
			19CT56	KME21A2235	KME22A2919	KME22G3007
			19CT57	KME21A3066	KME22A2970	KME22G3031
			19CT80	KME21B0255	KME22B0006	KME22G3234
			19DT05	KME21B0256	KME22B0174	KME22H0492
19DT30	KME21B0378	KME22B0360	KME22H0763			
19ET34	KME21B2209	KME22B0400	KME22H0785			
19ET36	KME21B2274	KME22B0665	KME22H0829			
19ET39	KME21B2275	KME22B0754	KME22H1041			
19ET40	KME21B2535	KME22B2315	KME22H1074			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229780	111781050	TUBE, MURPHY ENDOTRACH CUFF 5.0MM (10/BX) TELFLX	19ET61	KME21C0433	KME22B2554	KME22H1240
			19ET65	KME21C0643	KME22B2557	KME22J1082
			19ET69	KME21C1123	KME22B2631	KME22J1118
			19FT01	KME21C1125	KME22B2633	KME22J3075
			19FT45	KME21C1164	KME22B2635	KME22K0453
			19GT13	KME21E0006	KME22B2678	KME22K0605
			19GT63	KME21F0902	KME22B2793	KME22K0709
			19HT24	KME21F1216	KME22B2794	KME22K0751
			19HT31	KME21F2139	KME22B2920	KME22M0926
			19HT53	KME21F2198	KME22C0553	KME22M2488
			19IT54	KME21J0892	KME22C1122	KME22M2542
			19JT06	KME21J1496	KME22C3230	KME23A2849
			19JT10	KME21J1690	KME22C3394	KME23A2963
			19JT11	KME21J1813	KME22C3396	KME23B0513
			19JT14	KME21J1815	KME22D0097	KME23B0535
			19JT35	KME21J2074	KME22D0108	KME23B0881
			19JT63	KME21K0004	KME22D0244	KME23B1193
			19KT14	KME21K0115	KME22D0286	KME23B3061
			19KT38	KME21K0296	KME22D0322	KME23B3062
			19KT41	KME21K0875	KME22D0538	KME23B3318
19KT46	KME21K0877	KME22D0539	KME23C0879			
19LT17	KME21K1029	KME22D0602				
229781	111781055	TUBE, MURPHY ENDOTRACH CUFF 5.5MM (10/BX) TELFLX	18FT35	19LT17	KME21K1046	KME22E1502
			18FT39	20BT01	KME21K1078	KME22E1580
			18GG17	20BT02	KME21K2304	KME22E1581
			18GG18	20CT07	KME21K3005	KME22E1629
			18GT32	20DT18	KME21L0631	KME22E1778
			18IG10	20DT29	KME21L0756	KME22E1851
			18IG16	20DT32	KME21L1428	KME22E1931
			18IG17	20ET34	KME21L1429	KME22E2306
			18JG02	20ET49	KME21L1596	KME22F0073
			18JG13	20FT18	KME21L2141	KME22F0543
			18JG20	20FT64	KME21L2291	KME22F1077
			18JG35	20GT11	KME21L2736	KME22G1333
			18JG37	20GT34	KME22A1889	KME22G2578
			18KG18	20GT49	KME22A1890	KME22H1365
			18KG39	KME20H0487	KME22A2018	KME22H1508
			18KG41	KME20H1788	KME22A2281	KME22H1647
			18LG12	KME20H2574	KME22A2921	KME22H2499
			18LG14	KME20H2768	KME22A2969	KME22H2500
			18LG30	KME20H3002	KME22B0359	KME22H2812
			18LG39	KME20J1213	KME22C0207	KME22H3174
19AG29	KME20L0905	KME22C2843	KME22H3384			
19BG03	KME20L2170	KME22C2927	KME22H3402			
19BG04	KME20L2191	KME22C2969	KME22J0431			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229781	111781055	TUBE, MURPHY ENDOTRACH CUFF 5.5MM (10/BX) TELFLX	19BG18	KME20L2670	KME22C3213	KME22J0529
			19BG19	KME20M0304	KME22C3304	KME22J0977
			19BG33	KME20M2647	KME22C3375	KME22J1104
			19CG01	KME20M3094	KME22C3379	KME22J1305
			19CG02	KME21A1082	KME22D0052	KME22J1499
			19CG04	KME21A2080	KME22D0587	KME22J1601
			19CT28	KME21A2926	KME22D0719	KME22J1723
			19CT33	KME21B0543	KME22D0808	KME22J1977
			19CT56	KME21B1199	KME22D0844	KME22J2150
			19DT03	KME21B1630	KME22D0846	KME22J2412
			19DT05	KME21B1716	KME22D0855	KME22J2510
			19DT30	KME21B2541	KME22D0899	KME22J2755
			19DT35	KME21B2670	KME22D1001	KME22K0002
			19DT42	KME21C0592	KME22D1093	KME22K0136
			19DT46	KME21C0593	KME22D1115	KME22K0250
			19ET16	KME21C3340	KME22D1337	KME22K0251
			19ET18	KME21D2002	KME22D1348	KME22K1436
			19ET34	KME21E0008	KME22D1454	KME22K2261
			19FT01	KME21E0009	KME22D1568	KME22K2568
			19GT63	KME21F0901	KME22D1675	KME22L0033
			19HT31	KME21F1424	KME22D1820	KME22L0237
			19HT55	KME21F1459	KME22D2489	KME22L0610
			19IT01	KME21F2371	KME22D2490	KME23A0718
			19IT12	KME21J1239	KME22D3014	KME23A3005
			19IT14	KME21J1487	KME22D3071	KME23A3006
			19IT15	KME21J1691	KME22E1106	KME23B0514
			19IT19	KME21J1692	KME22E1398	KME23B1775
			19KT02	KME21J1854	KME22E1402	KME23B3051
			19KT41	KME21K0178		
			229782	111781060	TUBE, MURPHY ENDOTRACH CUFF 6.0MM (10/BX) TELFLX	18FT39
18GT32	19JT35	KME21C3344				KME22D2491
18HG19	19KT41	KME21D2006				KME22E2435
18IG10	19LT21	KME21D2007				KME22E2436
18IG16	19LT36	KME21F1423				KME22E2874
18IG30	20AT22	KME21J1549				KME22G0208
18JG35	20BT02	KME21K3004				KME22G0690
18LG15	20CT32	KME21L1432				KME22G0862
18LG30	20DT35	KME21L2739				KME22G1268
18LG39	20ET34	KME21M0028				KME22G2247
19AG15	20ET39	KME21M0163				KME22G2660
19AG25	20FT18	KME21M0799				KME22G3252
19CT33	20FT64	KME21M2429				KME22H1630
19CT36	KME20H0431	KME21M2911				KME22H2409
19CT80	KME20H2972	KME21M3080				KME22J0850
19DT25	KME20L2192	KME22A0809				KME22J1083

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229782	111781060	TUBE, MURPHY ENDOTRACH CUFF 6.0MM (10/BX) TELFLX	19DT35	KME20M0306	KME22A1219	KME22J1978
			19DT46	KME20M1285	KME22B1419	KME22J2703
			19ET16	KME20M2828	KME22C1123	KME22J2876
			19ET34	KME20M2905	KME22C1555	KME22K0908
			19ET69	KME21A1085	KME22C1556	KME22K1023
			19FT01	KME21A2932	KME22C2192	KME22K1050
			19GT13	KME21A3056	KME22C2845	KME22L0970
			19GT17	KME21B0722	KME22C2970	KME23A2886
			19GT63	KME21B1047	KME22D0540	KME23A3101
			19HT68	KME21B1628	KME22D1570	KME23B0515
			19IT01	KME21B2706	KME22D1620	KME23B0516
			19IT53	KME21C1054	KME22D1841	KME23B3050
			229783	111781065	TUBE, MURPHY ENDOTRACH CUFF 6.5MM (10/BX) TELFLX	18FG20
18FG31	19JT60	KME21C0646				KME22D1843
18GG17	19KT14	KME21C1045				KME22D3133
18HG21	19KT46	KME21C3350				KME22E0948
18HG30	19LT17	KME21D1229				KME22E0949
18IG16	19LT21	KME21E0011				KME22F0366
18IG30	20AT22	KME21F0146				KME22F0988
18JG16	20BT01	KME21J0899				KME22F3121
18JG20	20CT11	KME21K1657				KME22G0692
18KG01	20CT32	KME21K1658				KME22G1335
18LG18	20DT19	KME21K1713				KME22G3008
18LG30	20ET39	KME21K1748				KME22G3063
19AG01	20FT18	KME21K2308				KME22H0521
19AG25	20FT64	KME21L2302				KME22H1633
19AG29	KME20H0432	KME21L2792				KME22H3274
19BG18	KME20H2681	KME21L2886				KME22J0144
19CG02	KME20K1682	KME22A1379				KME22J1123
19DT25	KME20K3648	KME22A1884				KME22J2413
19DT34	KME20L2673	KME22A1885				KME22J2479
19DT35	KME20M1288	KME22A2922				KME22K0090
19DT46	KME20M2831	KME22C0332				KME22K1413
19ET36	KME20M2908	KME22C0413				KME22K1843
19FT01	KME20M2957	KME22C0504				KME23A0055
19FT02	KME21A2081	KME22C1125				KME23A0357
19FT45	KME21A2935	KME22C1127				KME23A2168
19GT63	KME21B0719	KME22C1581				KME23B0519
19HT55	KME21B1201	KME22C2971	KME23C0313			
19HT68	KME21B1825	KME22D0178				
229784	111781070	TUBE, MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX	18GG17	19HT31	KME21A2083	KME22C0505
			18GG25	19HT68	KME21A2187	KME22D2264
			18GT32	19HT69	KME21A2929	KME22D3154
			18HG19	19IT19	KME21B0718	KME22E1996
			18IG10	19JT11	KME21B1715	KME22E2224

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229784	111781070	TUBE, MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX	18IG16	19JT35	KME21B1946	KME22E2513
			18JG13	19JT63	KME21C0552	KME22F0202
			18JG21	19KT02	KME21C1043	KME22G1652
			18KG01	19LT21	KME21C1044	KME22G1698
			18LG12	20AT22	KME21C1910	KME22H0792
			18LG15	20BG11	KME21C1912	KME22H1588
			18LG30	20CT11	KME21F1682	KME22H2923
			19AG15	20DT22	KME21F1864	KME22J1124
			19AG25	20DT35	KME21K1635	KME22J2482
			19AG29	20FT64	KME21M0798	KME22J3020
			19BG03	20FT68	KME21M2432	KME22K1296
			19CT56	20GT34	KME22A0810	KME22M2817
			19DT25	20GT49	KME22A2177	KME23A0424
			19ET15	KME20H0436	KME22A2883	KME23A0625
			19ET18	KME20H2683	KME22A2924	KME23B1602
			19ET58	KME20K1677	KME22A3263	KME23B1972
			19GT13	KME20M0445	KME22B2796	KME23B3008
			19GT17	KME20M2650	KME22C0042	KME23C0514
			19HT02	KME21A2042	KME22C0206	
			229785	111781075	TUBE, MURPHY ENDOTRACH CUFF 7.5MM (10/BX) TELFLX	18FG21
18GG18	19JT11	KME21C1041				KME22E1853
18GG19	19JT63	KME21C1165				KME22E2748
18HG29	19LT21	KME21C2485				KME22F0118
18IG17	19LT36	KME21C3404				KME22F0603
18IG29	20AT26	KME21F0143				KME22G0210
18JG03	20BG11	KME21F0144				KME22G1168
18JG20	20CT11	KME21F1428				KME22H1123
18JG35	20DT19	KME21J1568				KME22H3070
18KG01	20FT64	KME21K1054				KME22J1084
18KG43	20GT11	KME21K2543				KME22K2080
18LG18	20GT34	KME21L2839				KME22K2470
19AG25	KME20J1216	KME22A0811				KME22K2512
19BG03	KME20J1852	KME22A1220				KME22L0201
19BG18	KME20J1853	KME22A2928				KME22L0202
19BG20	KME20J2131	KME22B0589				KME23A1196
19CT80	KME20K3678	KME22B2683				KME23B1974
19DT38	KME20L0953	KME22C0030				KME23B2026
19ET36	KME20M1291	KME22C2313				KME23B3049
19FT02	KME21A1263	KME22D0105				KME23C0156
19GT13	KME21A2084	KME22D0296	KME23C0157			
19HT55	KME21B1948	KME22D0589	KME23C0315			
19HT68	KME21B1950	KME22D2493				
229786	111781080	TUBE, MURPHY NASAL CUFF 8.0MM (10/BX) TELFLX	18GG18	19IT01	KME20M0409	KME22D2401
			18HG21	19IT15	KME21A1111	KME22E2898
			18IG16	19JT35	KME21A2095	KME22H2600

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229786	111781080	TUBE, MURPHY NASAL CUFF 8.0MM (10/BX) TELFLX	18KG01	19KT14	KME21B2212	KME22H3295
			18KG19	19LT21	KME21C2483	KME22J2877
			18LG18	20AT26	KME21C3405	KME22K2861
			19BG18	20CT11	KME22B0921	KME22K3117
			19BG33	20DT18	KME22B0923	KME22M1510
			19FT02	20FT68	KME22D0289	KME22M2546
			19GT13	20GT49	KME22D1811	KME23B0521
			19HT31	KME20J2105	KME22D1855	KME23B3047
			19HT55	KME20K3653		
229787	100181030	TUBE, ENDOTRACH NASAL UNCUFF 3.0MM (10/BX) TELFLX	18GG23	18LG18	19GT16	KME22F2174
			18HG17	19AG25	19IT53	KME22G1753
			18HG21	19BG19	19LT21	KME22G3017
			18IG17	19BG32		
229789	100181040	TUBE, ENDOTRACH NASAL UNCUFF 4.0MM (10/BX) TELFLX	18GG17	19HT31	KME21B2620	KME22F0607
			18GG34	19IT01	KME21B2621	KME22F2175
			18IG15	19IT53	KME21C1148	KME22G1862
			18JG03	19KT02	KME21C3436	KME22H0489
			18JG21	19LT21	KME21D0720	KME22K0562
			18KG18	20AT26	KME21D2148	KME22K1052
			18LG14	20CT11	KME21M2369	KME22L0811
			18LG15	20DT19	KME22B0041	KME23A1885
			19AG01	20DT32	KME22B0278	KME23A2176
			19AG25	KME20J0427	KME22B0588	KME23A2309
			19BG02	KME20M0303	KME22C0507	KME23B0501
			19ET39	KME21A2096	KME22D3254	KME23B2935
			19FT46	KME21A3064		
229790	100181045	TUBE, ENDOTRACH NASAL UNCUFF 4.5MM (10/BX) TELFLX	18FG31	19CT26	KME21D2150	KME22F1313
			18FT35	19CT53	KME21L2312	KME22F2624
			18GG17	19CT56	KME21M0208	KME22G0321
			18GG25	19CT57	KME21M0957	KME22G1706
			18GG34	19DT03	KME21M0958	KME22G2273
			18HG17	19DT34	KME21M0959	KME22H0491
			18IG14	19DT42	KME21M2377	KME22H1075
			18IG30	19ET61	KME21M2379	KME22H1324
			18JG16	19GT13	KME22B0281	KME22J0385
			18JG35	19GT16	KME22B1307	KME22J0386
			18JG37	19GT60	KME22B1308	KME22J1310
			18KG17	19HT59	KME22B1464	KME22K0563
			18KG36	19LT36	KME22C0411	KME22K0907
			18KG41	20AT26	KME22C2193	KME22K1053
			18LG14	20CT05	KME22D1004	KME22K1241
			18LG39	20CT32	KME22D1005	KME22K3166
			19AG25	KME21C1880	KME22E1405	KME23A1147
			19BG02	KME21D0722	KME22E1503	KME23A1148
			19BG32	KME21D1235	KME22E2767	KME23A2177
			19CG01	KME21D2149	KME22F0777	KME23B2938
19CG02						

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229791	100181050	TUBE, ENDOTRACH NASAL UNCUFF 5.0MM (10/BX) TELFLX	18FG31	19DT30	KME21L2790	KME22F3241
			18FT35	19DT35	KME21M2313	KME22G0696
			18GG23	19ET36	KME21M2314	KME22H0787
			18GT32	19ET39	KME21M2317	KME22H1645
			18HG17	19ET40	KME21M2318	KME22H2935
			18JG19	19ET61	KME21M2380	KME22H3404
			18JG21	19ET69	KME21M2481	KME22J2045
			18JG37	19FT01	KME22A0760	KME22K0981
			18KG18	19FT46	KME22A2297	KME22L0396
			18KG19	20AT26	KME22B0755	KME22L0444
			18KG21	20DT18	KME22B2679	KME22L0735
			18KG38	KME21A2039	KME22C0412	KME22L1059
			18KG39	KME21B2623	KME22C2194	KME23A0603
			18LG14	KME21C0543	KME22C2195	KME23A0617
			19AG01	KME21C1149	KME22D1145	KME23A1167
			19BG02	KME21C3398	KME22D2262	KME23A1169
			19BG18	KME21C3399	KME22D3177	KME23A1240
			19CG01	KME21D1236	KME22F1334	KME23A1249
			19CG02	KME21F0888	KME22F2203	KME23A3143
			19CT28	KME21K2310	KME22F3000	KME23B2940
19DT05	KME21L2303	KME22F3239	KME23C2375			
229792	100181055	TUBE, ENDOTRACH NASAL UNCUFF 5.5MM (10/BX) TELFLX	18FT35	19DT35	KME22A1382	KME22H0984
			18GG18	19ET36	KME22A2298	KME22H1510
			18GG23	19ET61	KME22A2299	KME22H1511
			18GT32	19ET65	KME22B1465	KME22H3258
			18IG16	19FT46	KME22B1666	KME22J1450
			18JG13	19KT02	KME22B2680	KME22K1154
			18JG35	20AT26	KME22C2315	KME22K1984
			18KG01	20CT11	KME22C2316	KME22K3163
			18KG39	20DT19	KME22D1671	KME22K3164
			18LG15	KME21A2087	KME22D2352	KME22L0576
			18LG18	KME21B1774	KME22E1006	KME22L0611
			18LG39	KME21B2624	KME22E1007	KME22L0734
			19AG01	KME21B2710	KME22E1406	KME23A1086
			19AG25	KME21B2711	KME22E1930	KME23B1606
			19CG02	KME21C0549	KME22F2204	KME23B3006
			19CT33	KME22A0758	KME22F2482	KME23C2632
19DT05						
229793	100181060	TUBE, ENDOTRACH NASAL UNCUFF 6.0MM (10/BX) TELFLX	18HG21	19ET40	KME22A2300	KME22F3384
			18KG41	19ET61	KME22C2820	KME22F3472
			18KG43	19ET65	KME22D1551	KME22G3018
			18LG15	19FT46	KME22D2353	KME22G3213
			19AG01	20CT32	KME22E1009	KME22K2121
			19BG33	KME21C1936	KME22E1144	KME22L0551
			19CT33	KME21K2311	KME22E1407	KME23B3004
19ET39	KME21M2381	KME22F2205				

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229794	100181065	TUBE, ENDOTRACH NASAL UNCUFF 6.5MM (10/BX) TELFLX	18FT35	19AG25	20CT11	KME22H0137
			18GG25	19CT33	KME22B1670	KME22J0202
			18JG35	19ET69	KME22E1010	KME22K2262
			18KG21	20BT02	KME22E1408	KME23B3131
			18LG18			
229795	111780035	TUBE, MURPHY ENDOTRACH CUFF 3.5MM (10/BX) TELFLX	18GG03	19FT22	KME20M1519	KME22G0575
			18GG24	19FT49	KME21A0965	KME22H0031
			18GT26	19GT43	KME21A2395	KME22H0797
			18HG36	19HT03	KME21B1964	KME22H1789
			18IG11	19JT01	KME21B2528	KME22H2875
			18IG24	19JT20	KME21C0337	KME22H3073
			18JG07	19JT55	KME22A1674	KME22J1158
			18JG14	19LT32	KME22B1172	KME22J1777
			18JG36	20AT09	KME22B1719	KME22K2438
			18KG22	20CT24	KME22C0647	KME22L1579
			18KG34	20DT15	KME22D0114	KME22L1617
			18LG31	20ET04	KME22D2229	KME22L1659
			19AG21	20ET56	KME22E0680	KME22L1725
			19AG36	20GG43	KME22E1094	KME23A2277
			19BG12	20GT18	KME22E2895	KME23A3144
			19CG07	KME20H2627	KME22E2896	KME23B0437
			19CT36	KME20K2726	KME22F1454	KME23B1557
19CT69	KME20L2489	KME22F1458	KME23C0152			
19DT42	KME20M0355	KME22F2737				
229796	111780040	TUBE, MURPHY ENDOTRACH CUFF 4.0MM (10/BX) TELFLX	18FG24	19BG12	20FT36	KME22F2739
			18GG03	19CG07	20FT60	KME22G0790
			18GG04	19CT69	20GT40	KME22H1304
			18GG24	19DT42	KME21C0405	KME22H2878
			18GT26	19ET58	KME21C1057	KME22J1447
			18HG24	19FT18	KME22C0550	KME22J3205
			18IG11	19GT43	KME22C0648	KME22K0482
			18IG14	19GT44	KME22C1816	KME23A1241
			18IG24	19JT55	KME22C2806	KME23A2001
			18JG07	19KT38	KME22D1073	KME23A2181
			18KG22	19LT32	KME22D2150	KME23A2683
			18KG34	20AT09	KME22E0683	KME23A2864
			18KG44	20AT49	KME22E0684	KME23A3032
			18LG31	20BT53	KME22E2814	KME23B1558
18LG33	20DT24	KME22E2897	KME23C0871			
19AG36	20ET50	KME22F1760				
229797	111780045	TUBE, AGT MURPHY ENDOTRACH CUFF 4.5MM (10/BX) TELFLX	18FG12	19GT43	KME20M1522	KME22J0267
			18GG26	19JT02	KME21A0198	KME22J1449
			18GT26	19KT37	KME21A0967	KME22J2193
			18HG24	19LT32	KME21A2397	KME22J2455
			18HG36	20AT13	KME21B2435	KME22K0574

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229797	111780045	TUBE, AGT MURPHY ENDOTRACH CUFF 4.5MM (10/BX) TELFLX	18IT05	20BT53	KME21C0325	KME22M1047
			18JG07	20CT21	KME22C1814	KME22M1305
			18JG13	20DT15	KME22C1815	KME22M1323
			18JG36	20ET50	KME22C2815	KME22M2819
			18KG44	20FT36	KME22D1075	KME22M2882
			18LG31	20FT60	KME22D3468	KME22M3040
			18LG33	20GG43	KME22E0686	KME23A2233
			19AG11	20GT40	KME22E2786	KME23A2332
			19AG21	KME20K2722	KME22F1459	KME23A2894
			19CG07	KME20K2733	KME22F2736	KME23A2895
			19DT06	KME20L0659	KME22G0186	KME23A3034
			19DT42	KME20L1620	KME22G0338	KME23B1307
			19FT18	KME20M0897	KME22H1790	KME23B1560
229798	111780050	TUBE, AGT MURPHY ENDOTRACH CUFF 5.0MM (10/BX) TELFLX	18FG12	19CT36	KME20M0952	KME22H1792
			18FG24	19DT06	KME21A0202	KME22H1890
			18GG03	19FT18	KME21A2486	KME22J2277
			18GT26	19GT44	KME21B0747	KME22J2836
			18HG24	19JT01	KME21C0326	KME22K0470
			18IG24	19KT37	KME22A1676	KME22K0876
			18IT05	19LT04	KME22B0553	KME22M1319
			18JG07	19LT32	KME22B1015	KME22M1954
			18KG34	20AT09	KME22B1725	KME23A0157
			18KG44	20BT53	KME22C2788	KME23A1282
			18LG06	20CT24	KME22C2818	KME23A1422
			18LG11	20DT24	KME22D0291	KME23A1483
			18LG31	20ET50	KME22D0293	KME23A2002
			19AG09	20ET56	KME22D2681	KME23A2278
			19AG11	20FT36	KME22F0733	KME23B0505
			19AG21	20FT58	KME22G0349	KME23B0506
			19AG36	KME20J1117	KME22G3043	KME23B1320
19CG07						
229799	111780055	TUBE, AGT MURPHY ENDOTRACH CUFF 5.5MM (10/BX) TELFLX	18FG23	19CT36	KME20K0629	KME22F1274
			18FG24	19DT06	KME20K2746	KME22G0791
			18GG03	19FT18	KME20L0661	KME22H0940
			18GT26	19HT80	KME20L1622	KME22H2230
			18IG24	19LT04	KME20M0351	KME22H2331
			18IT05	20AT09	KME21A0978	KME22H3266
			18IT06	20AT49	KME21A2488	KME22J2275
			18JG13	20BT53	KME21B1917	KME22J2875
			18JG36	20DT14	KME21C1735	KME22K2770
			18KG08	20FT59	KME21E1424	KME23A0296
			18KG34	20FT60	KME22C0548	KME23A0482
			18KG44	KME20H2632	KME22C0549	KME23A0601
			18LG11	KME20J1036	KME22D0160	KME23A2865
			18LG31	KME20J1207	KME22D0263	KME23A3176
			19AG21	KME20J1863	KME22E1097	KME23B0508
19AG36						

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229800	111780060	TUBE, AGT MURPHY ENDOTRACH CUFF 6.0MM (10/BX) TELFLX	18FG12	19AG21	20AT09	KME20K3201
			18FG23	19DT06	20CT23	KME20L0662
			18GG03	19DT42	20DT14	KME20L1712
			18GT26	19ET34	20DT24	KME20M0352
			18HG24	19FT18	20DT29	KME20M1524
			18IG11	19FT29	20ET04	KME21A2399
			18IG24	19GT44	20FT36	KME21B0095
			18JG07	19HT66	20FT58	KME21C0926
			18JG14	19JT01	KME20H1915	KME21J1233
			18KG08	19JT20	KME20H2130	KME21K2858
			18KG22	19KT17	KME20H2633	KME22C0547
			18LG06	19LT04	KME20J1035	KME22D1074
			18LG31	19LT08	KME20K1794	KME22D3373
			19AG09	19LT31	KME20K2067	KME22G3153
			19AG11	19LT32	KME20K2745	
229801	111780065	TUBE, AGT MURPHY ENDOTRACH CUFF 6.5MM (10/BX) TELFLX	18FG12	19ET12	KME20M0906	KME22C1836
			18FG23	19FT18	KME20M2652	KME22C2817
			18FG24	19FT49	KME20M3194	KME22D1076
			18GT26	19HT03	KME21B0072	KME22D2087
			18HG24	19HT66	KME21C0332	KME22F1763
			18IG11	19IT43	KME21C0923	KME22G1315
			18IG14	19JT20	KME21C3299	KME22G3080
			18IG24	19KT21	KME21D0466	KME22H0004
			18IT05	19LT08	KME21D0468	KME22H2233
			18JG07	19LT31	KME21D1903	KME22J1231
			18JG36	20AT49	KME21E1460	KME22J2727
			18KG10	20BT53	KME21L1006	KME22K0999
			18LG11	20DT14	KME21L1007	KME22K1158
			18LG31	20DT24	KME21M1776	KME22K1240
			19AG11	20FT59	KME21M3025	KME22M3068
			19AG36	KME20H0325	KME21M3026	KME23A2003
			19CG07	KME20J1141	KME22A0377	KME23A2192
19CT36	KME20L0727	KME22B1011	KME23A2962			
19DT10						
229802	111780070	TUBE, AGT MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX	18FG12	19IT12	KME20L0731	KME22D2351
			18FG24	19IT43	KME20L1681	KME22D2363
			18GG01	19JT01	KME20M0276	KME22D3366
			18GG26	19JT55	KME20M0955	KME22D3379
			18GT26	19KT21	KME20M2655	KME22E1102
			18HG36	19LT08	KME20M3190	KME22E2443
			18IG14	19LT31	KME21A0979	KME22F1333
			18IT05	20AT13	KME21B0096	KME22G0171
			18JG07	20AT49	KME21B0154	KME22G1728
			18JG36	20BT53	KME21B1518	KME22G3156
			18KG10	20CT24	KME21B1963	KME22H0815

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229802	111780070	TUBE, AGT MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX	18KG34	20DT09	KME21C0334	KME22H1321
			18LG11	20DT24	KME21C2295	KME22J0749
			18LG31	20ET50	KME21J1237	KME22K1286
			18LG35	20FT39	KME21K0732	KME22K2440
			19AG36	20FT59	KME21M1046	KME22L1550
			19CG07	20GG27	KME22A0450	KME22L1660
			19CT36	KME20G0273	KME22A1594	KME22L2544
			19CT69	KME20H2634	KME22B0425	KME22M1758
			19ET12	KME20J2813	KME22B1002	KME22M1953
			19FT14	KME20J3155	KME22C0157	KME22M2759
			19GT39	KME20K1111	KME22C0473	KME22M2837
			19GT44	KME20K2068	KME22D1110	KME22M2886
			19HT04	KME20K3244	KME22D2152	KME23A1734
			19HT80	KME20L0730		
229803	111780075	TUBE, AGT MURPHY ENDOTRACH CUFF 7.5MM (10/BX) TELFLX	18FG21	19ET39	KME20K3625	KME22D3251
			18FG31	19FT37	KME20L2680	KME22E1151
			18GG19	19GT65	KME21A1113	KME22E1506
			18HG23	19HT55	KME21A2088	KME22F1667
			18HG30	19IT01	KME21A3085	KME22G0264
			18IG11	19IT19	KME21B0714	KME22G0386
			18IG16	19JT38	KME21B2708	KME22G3151
			18IG30	19JT63	KME21C0596	KME22G3270
			18JG35	19LT21	KME21C1167	KME22H2700
			18JG36	20AT26	KME21C2517	KME22J0257
			18KG19	20BG11	KME21L0761	KME22J1105
			18KG41	20CT12	KME21L0859	KME22J1335
			18LG18	20DT22	KME22A1383	KME22K3101
			19AG03	20ET04	KME22B2543	KME22M1594
			19BG04	20ET49	KME22B2906	KME22M2312
			19BG19	20FT65	KME22C0033	KME23A0829
			19CG03	KME20H0627	KME22C0034	KME23A1588
			19CT36	KME20H2734	KME22D0652	KME23C0868
			19DT25	KME20J0490	KME22D1158	KME23C1106
			19DT34	KME20J1392	KME22D1372	KME23C1182
19ET20	KME20K0991	KME22D2254				
229804	111780080	TUBE, AGT MURPHY ENDOTRACH CUFF 8.0MM (10/BX) TELFLX	18FG33	19ET20	KME21A1299	KME22D2122
			18GG23	19ET39	KME21A2090	KME22D2509
			18GG34	19ET40	KME21B0716	KME22D2752
			18HG30	19GT65	KME21B1783	KME22E1153
			18IG11	19HT68	KME21B1962	KME22G0816
			18IG18	19IT15	KME21B2204	KME22H0192
			18IG30	19LT21	KME21B2676	KME22H2094
			18JG21	20AT26	KME21C0599	KME22H2235
			18JG35	20CT12	KME21C1953	KME22H3271
			18KG19	20CT13	KME21D2027	KME22J0148
			18KG41	20CT32	KME21J1694	KME22J0532

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229804	111780080	TUBE, AGT MURPHY ENDOTRACH CUFF 8.0MM (10/BX) TELFLX	19AG15	20DT35	KME21J1930	KME22J3210
			19AG19	20FT19	KME21K1106	KME22K2844
			19BG04	KME20J0491	KME22B0591	KME22M0776
			19BG19	KME20K1107	KME22C2317	KME22M1210
			19DT05	KME20M3089	KME22D0156	KME23A3014
			19DT25	KME20M3090	KME22D1159	KME23B1497
			19DT38			
229805	111780085	TUBE, MURPHY ENDOTRACH CUFF 8.5MM (10/BX) TELFLX	18FG21	18LG18	20CT32	KME22B0760
			18GG25	19AG19	20GT34	KME22F2487
			18IG29	19DT36	KME21B2675	KME22H1797
			18IG30	19FT02	KME21F0887	KME23A1524
			18LG12			
229807	100180030	TUBE, ENDOTRACH ORAL UNCUFF 3.0MM (10/BX) TELFLX	18GG26	19FT19	20ET56	KME22D1514
			18GT26	19FT29	KME20K1810	KME22F2666
			18JG14	19HT04	KME20M0946	KME22G0172
			19AG21	19IT43	KME20M2783	KME22G1751
			19CG07	19KT17	KME21A2389	KME22L1527
			19CT69	19LT08	KME21B0075	KME22M1326
			19DT42	20DT15	KME21B1923	KME23A1328
			19ET34	20DT24	KME21L1114	
229808	100180035	TUBE, ENDOTRACH ORAL UNCUFF 3.5MM (10/BX) TELFLX	18IG24	19FT19	20DT24	KME22D1515
			18JG07	19FT29	20ET56	KME22G0175
			18LG06	19HT34	KME20K1811	KME22K0607
			18LG33	19IT43	KME20M0947	KME22K3167
			19AG36	19JT02	KME21A0060	KME23A1605
			19CT69	19LT08	KME21A2390	KME23A1884
			19DT42	19LT32	KME21B0599	KME23C0905
			19ET34			
229809	100180040	TUBE, ENDOTRACH ORAL UNCUFF 4.0MM (10/BX) TELFLX	18GT26	19CT36	KME21A0061	KME22H0033
			18IG24	19DT42	KME21B1925	KME22H3078
			18JG07	19FT29	KME21M1736	KME22J0843
			18KG10	19FT49	KME22C0700	KME22K0856
			18LG33	20DT24	KME22F0734	KME23A0056
			19AG21	KME20L2493	KME22G0184	KME23B2624
			19BG12	KME20M0278		
229810	100180045	TUBE, ENDOTRACH ORAL UNCUFF 4.5MM (10/BX) TELFLX	18GT26	19GT43	20ET04	KME22J3224
			18KG10	19HT66	20GG27	KME22K0489
			18LG06	19JT55	KME20H2636	KME22L0279
			18LG33	19KT17	KME20K1812	KME22L1109
			19AG11	19LT08	KME20K2573	KME23A1186
			19AG21	19LT32	KME21L2619	KME23A3133
			19CT69	20AT49	KME21M1773	KME23B1541
			19DT06	20BT53	KME22E2900	KME23B2654
			19FT29	20DT15	KME22G0347	KME23C0302
			19FT49			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
229811	100180050	TUBE, ENDOTRACH ORAL UNCUFF 5.0MM (10/BX) TEFLX	18GT26	19CG07	19LT32	KME22G0655
			18IG24	19DT06	20DT14	KME22H3083
			18KG10	19ET12	20GG27	KME22K0927
			18KG22	19ET58	KME20J1208	KME22L0442
			18LG06	19GT44	KME21B0077	KME23A1609
			18LG33	19HT04	KME21C0441	KME23B1319
			19AG21	19HT66	KME22D3452	KME23C1958
			19AG36	19IT12	KME22F0736	
229812	100180055	TUBE, ENDOTRACH ORAL UNCUFF 5.5MM (10/BX) TEFLX	18IG24	19FT22	20AT49	KME21A0310
			18JG36	19FT49	20BT53	KME21B0602
			18LG33	19HT04	20CT24	KME22F0738
			19AG11	19HT34	20DT24	KME22L0860
			19AG21	19HT80	20GG27	KME23A1607
			19DT06	19IT43	KME20K2673	KME23B0500
			19ET58	19JT02	KME20M0950	KME23C1479
229813	100180060	TUBE, ENDOTRACH ORAL UNCUFF 6.0MM (10/BX) TEFLX	18KG08	19HT34	20GG27	KME22D3370
			18LG33	19IT12	KME20K1813	KME22F2673
			19AG21	19IT43	KME20L2063	KME22J0208
			19BG12	19LT32	KME20M1586	KME22L0861
			19ET12	20AT49	KME21A2393	KME22M3021
			19FT19	20DT15	KME21C3189	
229815	100180070	TUBE, ENDOTRACH ORAL UNCUFF 7.0MM (10/BX) TEFLX	18HG40	19FT19	KME21A2555	KME22K1045
			18LG33	19GT46	KME21L0242	KME22L1142
			19CT36	19LT31	KME22H2313	
242472	112080050	TUBE, MAGILL ENDOTRACH NASAL 5.0MM (10/BX) TELFLX	18FG07	19KT38	KME20K1512	KME22B0593
			18JG33	19LT07	KME20L1145	KME22E2347
			18JG36	20AT41	KME20M2383	KME22H1178
			18LT16	20CT13	KME21B0357	KME22H3021
			19CT26	20ET17		
242473	112080055	TUBE, MAGILL ENDOTRACH NASAL 5.5MM (10/BX) TELFLX	18FG02	19CT64	KME20K2579	KME22B0079
			18IG29	19GT31	KME21A2569	KME22B0080
			18KG28	19IT33	KME21C1978	KME22E2802
			19AG39	20DT18	KME21K0764	KME22H2250
			19BG02	KME20G0435	KME21K1278	
252280	100380045	TUBE, MAGILL ENDO UNCUFF 4.5MMTELFLX	18FG22	19JT51	KME20L0425	KME22G1456
			18GT16	19KT13	KME21A2196	KME22G2369
			18IG04	19LT31	KME21A3268	KME22H2833
			18JG08	20AT06	KME21C0914	KME22J2485
			19CT42	20DT04	KME22D0877	KME22K2783
			19HT05	KME20J0756	KME22E0896	KME23B2478
			19HT61	KME20K2144	KME22F1601	
265412	112480035	TUBE, MURPHY ENDOTRACH CUFF 3.5MM (10/BX) TELFLX	18FG14	19LT27	KME21B1099	KME22G0614
			18FG22	20AT54	KME21C1766	KME22G1180
			18GG15	20BT42	KME21C2690	KME22G1181
			18GG16	20BT51	KME21C3242	KME22G1876
			18GG34	20CT07	KME21D0295	KME22H2913

MMS #	MFG Catalog #	Description	Affected Lot(s)			
265412	112480035	TUBE, MURPHY ENDOTRACH CUFF 3.5MM (10/BX) TELFLX	18HG08	20CT54	KME21D2960	KME22H3442
			18IG01	20DT20	KME21K1922	KME22J2322
			18IG03	20DT22	KME21M0820	KME22J3348
			18IG04	20DT24	KME21M2050	KME22L1190
			18IG13	20ET30	KME21M3073	KME22L1191
			18LG30	20FT43	KME22A0416	KME22L1547
			19BG16	KME20H2894	KME22A1499	KME23A2353
			19BG35	KME20J0754	KME22A2444	KME23A2534
			19DT35	KME20J3097	KME22B2850	KME23B0177
			19ET51	KME20M2591	KME22C1453	KME23B0178
			19GT51	KME20M3049	KME22C1455	KME23B0590
			19IT38	KME20M3050	KME22C2694	KME23B1301
			19JT20	KME21A0349	KME22D1761	KME23B2516
			19JT50	KME21A3147	KME22D2536	KME23B2606
			19KT37	KME21A3239	KME22E0601	KME23C0082
			19KT60	KME21B0749	KME22F2453	KME23C0083
			19LT03			
265413	112480040	TUBE, MURPHY ENDOTRACH CUFF 4.0MM (10/BX) TELFLX	18FG14	20AT06	KME21K1187	KME22F1692
			18FG22	20AT09	KME21K1188	KME22F2455
			18GG16	20AT39	KME21K1977	KME22G1300
			18IG03	20AT54	KME21M1084	KME22G2312
			18IG29	20BG14	KME21M1842	KME22G2367
			18JG08	20BT01	KME21M1845	KME22H1126
			18JG09	20BT53	KME21M2733	KME22H1212
			18KG25	20CT12	KME21M2734	KME22H2302
			18LG04	20CT27	KME21M2735	KME22J1615
			18LG30	20CT54	KME21M2736	KME22K0283
			19AG38	20DT04	KME22B2019	KME22K0632
			19CG05	20ET30	KME22C1693	KME22K1096
			19DT35	20FT58	KME22D0057	KME22M1155
			19DT36	KME20G0182	KME22D0180	KME22M1478
			19ET09	KME20H0372	KME22D0819	KME22M1659
			19ET31	KME20H0373	KME22D1066	KME23A1486
			19ET51	KME20H1889	KME22D1832	KME23B1677
			19FT25	KME20J0924	KME22D3269	KME23B2354
			19HT60	KME20J2043	KME22E0602	KME23B3174
			19IT32	KME21C2368	KME22E2648	KME23C0081
19IT43	KME21D1287	KME22E2681	KME23C0692			
19KT60	KME21D2961	KME22F1440	KME23C1442			
19LT03						
265414	112480045	TUBE, MURPHY ENDOTRACH CUFF 4.5MM (10/BX) TELFLX	18FG02	19JT19	KME20M0861	KME22D1836
			18FG14	19JT21	KME20M1576	KME22D2080
			18GT16	19JT50	KME20M1925	KME22D2309
			18GT32	19KT58	KME20M1926	KME22D2310
			18HG08	19LT09	KME20M3207	KME22D2481
			18HG37	19LT27	KME21A0769	KME22F0391

MMS #	MFG Catalog #	Description	Affected Lot(s)			
265414	112480045	TUBE, MURPHY ENDOTRACH CUFF 4.5MM (10/BX) TELFLX	18IG13	20AT06	KME21B0450	KME22F1537
			18IG29	20BG15	KME21B1090	KME22F2459
			18JG09	20BT01	KME21C0314	KME22F2526
			18JG35	20BT42	KME21C1807	KME22G1183
			18KG25	20BT43	KME21D0652	KME22H0051
			18KG34	20CT07	KME21D1746	KME22H3521
			18LG03	20CT12	KME21K1189	KME22J2886
			19AG04	20CT54	KME21K1920	KME22K1019
			19AG07	20DT01	KME21M1895	KME22M2104
			19AG08	20DT09	KME21M2015	KME23A1068
			19AG23	20DT39	KME21M2776	KME23A1487
			19AG38	20ET30	KME21M3038	KME23A2313
			19BG35	20FT58	KME22B2939	KME23B1676
			19CG05	20GT22	KME22C0808	KME23B2166
			19DT09	KME20H0774	KME22C0928	KME23B3254
			19FT12	KME20J0792	KME22D0939	KME23B3255
			19IT38	KME20J2687	KME22D0940	KME23B3257
			278648	112080100	TUBE, MAGILL ENDOTRACH CUFF 10.0MM (10/BX)	18FT26
18IT18	19IT31	20GT57				KME22C0947
18KG10	20BT13	KME20K1005				
349769	506525	TUBE, ENDOTRACH UNCUFF W/STYLETTE 2.5MM (10/BX) TELFLX	18FG05	19AG18	20AT46	KME22B0302
			18FG13	19AG21	20BT13	KME22B1071
			18GG17	19AG37	20BT27	KME22C3357
			18GG28	19CT32	20CT12	KME22D0784
			18HG16	19CT55	20CT22	KME22D1062
			18HG25	19ET45	20CT36	KME22D3648
			18HG29	19GT05	20CT51	KME22E1130
			18HG35	19GT20	20DT37	KME22F0161
			18HG38	19GT38	20ET16	KME22F1501
			18IG06	19GT59	20ET28	KME22G0561
			18IG12	19HT05	20FT44	KME22G1100
			18IG24	19HT15	KME20L0615	KME22H2626
			18IG27	19HT40	KME20M2047	KME22K0013
			18JG11	19HT80	KME21A1889	KME22K0331
			18JG14	19IT29	KME21A1890	KME22K2365
			18JG26	19KT08	KME21A2804	KME22L0760
			18JG34	19KT25	KME21B1143	KME22M0741
			18KG23	19KT41	KME21B1672	KME23A0852
			18LG05	19KT48	KME21C0063	KME23A1053
			18LG12	19LT04	KME21C0301	KME23B2941
18LG22	20AT25	KME21M1833	KME23C0925			
349770	506530	TUBE, ENDOTRACH UNCUFF W/STYLETTE 3.0MM (10/BX) TELFLX	18FG05	19CT65	20BT47	KME22B2326
			18GG12	19ET72	20CT05	KME22B2327
			18GG17	19FT01	20CT12	KME22C1075
			18GG31	19GT05	20CT36	KME22C1826

MMS #	MFG Catalog #	Description	Affected Lot(s)			
349770	506530	TUBE, ENDOTRACH UNCUFF W/STYLETTE 3.0MM (10/BX) TELFLX	18GT32	19GT20	20DT37	KME22C3192
			18HG12	19GT59	20ET16	KME22D0787
			18HG16	19GT65	20ET28	KME22E0223
			18IG02	19HT40	KME20J0318	KME22E1029
			18IG06	19HT61	KME20K0152	KME22F0563
			18IG15	19HT80	KME20K0153	KME22G0910
			18IG21	19IT29	KME20K0154	KME22G0911
			18IG27	19IT39	KME20K3587	KME22H0780
			18JG02	19JT18	KME20L1467	KME22H2909
			18JG23	19JT22	KME20L1468	KME22J0577
			18JG32	19JT29	KME20L1838	KME22K0015
			18JG34	19JT49	KME20M2048	KME22K0016
			18KG15	19JT57	KME20M2480	KME22K0202
			18KG23	19KT08	KME21A0669	KME22K0335
			18LG02	19LT04	KME21A1897	KME22K0918
			18LG22	19LT48	KME21A2807	KME22K1688
			18LG28	20AT09	KME21B1674	KME22K2774
			19AG21	20AT32	KME21B1675	KME22K3077
			19AG24	20AT44	KME21C0489	KME23A1501
			19AG27	20AT46	KME21M0865	KME23A1648
			19AG37	20AT53	KME21M2740	KME23A2751
			19BG25	20BT23	KME21M2741	KME23B2942
			19BG40	20BT37	KME22B2144	KME23C1505
349771	506535	TUBE, ENDOTRACH UNCUFF 3.5MM (10/BX) TELFLX	18FG05	19GT20	20DT14	KME22B0539
			18FG20	19GT45	20DT37	KME22B0725
			18GG12	19GT59	20ET16	KME22B1028
			18GG17	19GT65	20ET21	KME22B1562
			18GG28	19HT05	20ET28	KME22C0916
			18GT32	19HT40	20FT44	KME22C1821
			18HG12	19HT61	KME20K1624	KME22C1940
			18HG16	19HT80	KME20K3589	KME22D0350
			18JG17	19IT29	KME20L0616	KME22D1078
			18KG07	19JT49	KME20L0617	KME22F0162
			18LG02	19KT08	KME20M2049	KME22F2162
			18LG22	19LT04	KME20M3274	KME22G3101
			18LG28	19LT48	KME20M3275	KME22G3233
			19AG21	20AT32	KME20M3278	KME22H2624
			19BG40	20AT53	KME21A1899	KME22J0803
			19CT58	20BT05	KME21A2809	KME22J3077
			19ET38	20BT23	KME21B0575	KME22K1550
			19ET60	20BT47	KME21C0027	KME22K2775
			19ET72	20CT12	KME21C0491	KME23A3203
			19FT01	20CT22	KME22A0057	KME23C0591
19GT05	20CT36	KME22A0206	KME23C0959			
19GT10	20CT44					

MMS #	MFG Catalog #	Description	Affected Lot(s)			
349772	506540	TUBE, ENDOTRACH UNCUFF W/STYLETTE 4.0MM (10/BX) TELFLX	18FG05	19BG21	20AT32	KME21A2814
			18FG13	19BG40	20AT53	KME21B1689
			18FG22	19ET38	20BT23	KME21B1692
			18GG17	19ET60	20CT12	KME22B1078
			18GG28	19ET72	20CT22	KME22C1972
			18GT32	19FT01	20CT36	KME22C2370
			18HG16	19FT18	20ET49	KME22D0541
			18JG02	19FT35	20ET56	KME22D0788
			18JG17	19GT05	20GG45	KME22E0224
			18JG32	19GT28	KME20J2322	KME22E1238
			18JG34	19HT15	KME20K0940	KME22F0048
			18KG29	19HT40	KME20K2247	KME22G2817
			18KG35	19HT80	KME20L0621	KME22H2454
			18LG02	19IT39	KME20L1470	KME22K0336
			18LG12	19JT57	KME20L2029	KME23C0353
			18LG18	19LT04	KME20M2093	KME23C0482
			18LG28	19LT17	KME20M2481	KME23C0960
			19AG21	20AT09	KME21A0670	
349773	506545	TUBE, ENDOTRACH UNCUFF 4.5MM (10/BX) TELFLX	18FG05	19BG21	20AT09	KME20L2030
			18FG13	19CT49	20BG16	KME20M2482
			18GG20	19FT35	20BT05	KME21A1892
			18HG16	19GT20	20CT05	KME21A1893
			18IG06	19GT28	20CT36	KME21A2811
			18IG12	19GT38	20CT44	KME21B0576
			18JG04	19GT59	20DT37	KME22G3102
			18JG14	19GT65	20GG45	KME22J2956
			18JG21	19HT24	KME20H2927	KME22K2778
			18JG34	19HT80	KME20H2928	KME22L3057
			18KG35	19IT29	KME20K2248	KME23A2083
			18LG02	19KT25	KME20L0620	KME23B1203
			18LG18	19KT48	KME20L1471	KME23C0961
			18LG28	19LT38	KME20L1472	
349774	506550	TUBE, ENDOTRACH UNCUFF W/STYLETTE 5.0MM (10/BX) TELFLX	18FG01	19DT33	20AT25	KME20J2114
			18FG05	19ET54	20AT44	KME20J3032
			18FG13	19ET77	20BT05	KME20K1344
			18GG17	19FT24	20BT27	KME20L1842
			18GG31	19GT59	20BT37	KME20L1843
			18HG12	19GT65	20BT51	KME20L2299
			18HG19	19HT24	20CT05	KME20M2094
			18JG21	19HT61	20CT12	KME20M2483
			18JG32	19IT29	20CT36	KME21A2812
			18KG41	19JT42	20DT09	KME21B1693
			18LG18	19JT49	20ET56	KME21C0020
			18LG28	19JT57	20GT06	KME22B1976
			19AG21	19KT25	20GT38	KME22H3204
19DT25	20AT14	KME20H1790	KME23C2070			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
349775	504550	TUBE, ENDOTRACH CUFF 5.0MM (10/BX) TELFLX	18FG03	19GT20	KME20K0467	KME22B0172
			18FG13	19GT38	KME20K1581	KME22B1395
			18GG12	19GT65	KME20L0618	KME22B1632
			18GG17	19HT05	KME20L1090	KME22C0436
			18GG31	19HT40	KME20L2027	KME22C1654
			18HG12	19HT66	KME20M0019	KME22C3003
			18IG02	19JT70	KME21A2076	KME22D0145
			18IG12	19KT08	KME21A2787	KME22D0825
			18IG21	19LT17	KME21B1140	KME22E1129
			18IG27	19LT32	KME21C0059	KME22H0163
			18JG02	20AT14	KME21C2363	KME22H1169
			18KG15	20AT25	KME21D1964	KME22H1703
			18LG25	20BT27	KME21F1551	KME22H1709
			18LG37	20BT37	KME21G1707	KME22J3162
			19AG37	20CT36	KME21K0545	KME22J3163
			19BG07	20CT54	KME21K1249	KME22K0743
			19BG21	20DT14	KME21M0011	KME22K1246
			19CT71	20ET04	KME21M1831	KME22K1247
			19ET54	20ET21	KME21M2074	KME23B0167
			19FT14	20GT38	KME21M2647	KME23C0967
			19FT24	KME20H0128	KME22B0171	KME23C2505
19GT10	KME20K0001					
349776	504555	TUBE, ENDOTRACH CUFF 5.5MM (10/BX) TELFLX	18FG20	19GT10	20CT12	KME22A2940
			18FG22	19GT20	20CT36	KME22A3267
			18FG25	19GT38	20ET21	KME22A3268
			18FT32	19GT59	KME20G0388	KME22B1633
			18GG14	19GT65	KME20J0319	KME22B1751
			18HG12	19HT61	KME20J2321	KME22C3207
			18HG27	19IT25	KME20K0765	KME22D0590
			18IG06	19JT29	KME20K0766	KME22D1785
			18IG24	19JT42	KME20K0935	KME22E0619
			18IG27	19JT49	KME20K1582	KME22E0620
			18JG02	19KT30	KME20L0619	KME22F1666
			18JG14	19KT41	KME20L1857	KME22G0225
			18KG07	19LT04	KME20L2028	KME22G2294
			18KG23	19LT50	KME20M0020	KME22H3105
			18KG35	20AT09	KME21A1902	KME22J1458
			18LG31	20AT44	KME21A2796	KME22J1459
			19AG27	20BT05	KME21B2681	KME22J2714
			19BG04	20BT13	KME21C2707	KME22L1104
			19CT55	20BT27	KME21H1189	KME22L1105
			19ET72	20BT37	KME22A0148	KME23C0571
			19FT30	20CT05	KME22A2939	KME23C0590
19GT05						

MMS #	MFG Catalog #	Description	Affected Lot(s)			
349777	504560	TUBE, ENDOTRACH CUFF W/STYLETTE 6.0MM (10/BX) TELFLX	18FG03	19FT24	KME20H0129	KME21K1165
			18FG15	19FT35	KME20H0906	KME21L0442
			18FG28	19GT05	KME20H1791	KME21L0443
			18GG24	19GT45	KME20H1847	KME21L1900
			18IG27	19HT05	KME20J0848	KME21L2539
			18JG02	19HT61	KME20K1550	KME21M1391
			18JG17	19JT42	KME20L0336	KME22A1432
			18JG26	19KT17	KME20L2298	KME22B1411
			18JG34	19KT21	KME20M0021	KME22C1068
			18KG09	19KT65	KME20M0492	KME22C2858
			18KG35	19LT04	KME21A0380	KME22D0479
			18LG05	19LT23	KME21A0817	KME22D0482
			18LG12	20AT09	KME21A1160	KME22D0591
			18LG18	20AT37	KME21A1770	KME22D0592
			18LG37	20AT53	KME21B1137	KME22D2846
			19AG27	20BT51	KME21B1138	KME22F1495
			19AG37	20CT05	KME21B1139	KME22G2465
			19BG12	20CT12	KME21C3255	KME22G2466
			19CT32	20CT36	KME21E1412	KME23C0829
			19CT41	20DT37	KME21G2207	KME23C0830
			19CT55	20ET04	KME21H0839	KME23C0831
			19CT58	20FT10	KME21H0972	KME23C1796
			19ET66	20FT35	KME21H1124	KME23C2255
			19ET77	20GT06	KME21K0660	KME23C2503
			19FT14			
349778	504565	TUBE, ENDOTRACH CUFF WSTYLETTE6.5MM (10/BX) TELFLX	18FG05	19HT10	KME20K1551	KME22A3028
			18FG13	19HT66	KME20K1628	KME22A3156
			18FG25	19IT14	KME20K2246	KME22A3313
			18GG20	19IT25	KME20L0488	KME22B0067
			18HG07	19JT18	KME20L1840	KME22B1218
			18HG12	19JT29	KME20M0022	KME22B1409
			18HG19	19JT42	KME20M0854	KME22C0499
			18HG32	19JT57	KME20M1102	KME22C0863
			18HG38	19KT08	KME20M1579	KME22C1073
			18IG06	19KT41	KME21A0085	KME22C1225
			18IG09	19KT55	KME21A1040	KME22C1655
			18IG15	19LT32	KME21A1292	KME22C1941
			18JG04	20AT09	KME21A2794	KME22D0747
			18JG21	20AT14	KME21B1141	KME22D0755
			18JG34	20AT25	KME21B1669	KME22F1616
			18KG35	20AT53	KME21B2155	KME22F1617
			18LG39	20BT27	KME21B2501	KME22G0632
			19AG01	20CT05	KME21B2502	KME22G1224
			19AG24	20CT12	KME21B2855	KME22G2467
			19AG27	20CT36	KME21C0025	KME22H1900
19AG37	20CT44	KME21C1202	KME22H3348			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
349778	504565	TUBE, ENDOTRACH CUFF WSTYLETTE6.5MM (10/BX) TELFLX	19CT32	20DT09	KME21C1331	KME22J1135
			19CT41	20ET04	KME21C1790	KME22J1136
			19CT58	20ET16	KME21C1897	KME22J1928
			19CT65	20GT14	KME21C2188	KME22J3352
			19DT06	KME20G0189	KME21C2262	KME22K2413
			19DT20	KME20G0316	KME21C2835	KME22K2414
			19DT33	KME20H2312	KME21K1128	KME22L0241
			19ET32	KME20H2585	KME21L0328	KME22L0242
			19ET54	KME20H2586	KME21L2340	KME22L1792
			19ET60	KME20H2587	KME21M2283	KME23B2021
			19FT24	KME20J0320	KME21M2446	KME23C0112
			19FT41	KME20J0570	KME21M2568	KME23C0572
			19GT28	KME20K0465	KME22A0675	KME23C1198
			19GT38	KME20K0939	KME22A2415	KME23C2268
349779	504570	TUBE, ENDOTRACH CUFF W/STYLETTE 7.0MM (10/BX)	18FG03	20BT27	KME21B2460	KME22B0522
			18FG13	20BT37	KME21B2719	KME22B0523
			18FG21	20BT47	KME21B2721	KME22B0536
			18FG22	20CT05	KME21C0123	KME22B0537
			18FG30	20CT12	KME21C0376	KME22B0624
			18FT32	20CT22	KME21C0421	KME22B0633
			18GG02	20CT27	KME21C0620	KME22B0684
			18GG08	20CT36	KME21C0634	KME22B0849
			18GG20	20CT44	KME21C1016	KME22B0943
			18GT32	20DT14	KME21C1333	KME22B0985
			18HG07	20DT37	KME21C1416	KME22B1396
			18HG16	20ET04	KME21C1434	KME22B1408
			18HG19	20ET16	KME21C1489	KME22B1472
			18HG32	20FT10	KME21C1593	KME22B1560
			18HG35	20FT39	KME21C1653	KME22B1634
			18HG38	20FT54	KME21C1703	KME22B1636
			18IG02	20FT61	KME21C1786	KME22B1694
			18IG06	20GG45	KME21C1896	KME22B1752
			18IG09	20GG47	KME21C2065	KME22B1753
			18IG12	20GT06	KME21C2146	KME22B1793
			18IG17	20GT25	KME21G0508	KME22B1794
			18IG24	20GT31	KME21G1162	KME22B1873
			18IG27	20GT38	KME21G1179	KME22B2124
			18JG02	KME20H0972	KME21G1354	KME22B2125
			18JG04	KME20H1935	KME21G1355	KME22B2752
			18JG21	KME20H2588	KME21G1379	KME22B2884
			18JG32	KME20H2783	KME21G1385	KME22B2901
			18JG37	KME20J0058	KME21G1705	KME22C0067
			18KG23	KME20J0060	KME21G1782	KME22C0241
			18KG35	KME20J0156	KME21G2060	KME22C0242
			18KG38	KME20J0321	KME21G2061	KME22C0532
18LG05	KME20J0573	KME21G2090	KME22C0864			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
349779	504570	TUBE, ENDOTRACH CUFF W/STYLETTE 7.0MM (10/BX)	18LG08	KME20J0574	KME21H0604	KME22C1020
			18LG12	KME20J0803	KME21H0605	KME22C1317
			18LG18	KME20J0804	KME21H0607	KME22C1860
			18LG31	KME20J0805	KME21H0776	KME22C1942
			18LG39	KME20J0946	KME21H0777	KME22C2352
			19AG05	KME20J1112	KME21J1183	KME22C2353
			19AG09	KME20J1233	KME21J1206	KME22C2649
			19AG31	KME20J1692	KME21J1244	KME22C2654
			19AG34	KME20J1694	KME21J1255	KME22C2723
			19BG04	KME20J2175	KME21J1258	KME22C2859
			19CG06	KME20J2577	KME21J1482	KME22C2987
			19CT17	KME20J2957	KME21J1604	KME22C2988
			19CT23	KME20J3035	KME21J1747	KME22C3185
			19CT32	KME20K0026	KME21J2053	KME22D0594
			19CT37	KME20K0360	KME21K0016	KME22D0711
			19CT49	KME20K0361	KME21K0261	KME22D0716
			19CT55	KME20K0768	KME21K0271	KME22D0831
			19CT58	KME20K0941	KME21K0272	KME22D0836
			19CT65	KME20K1148	KME21K0314	KME22D1498
			19CT71	KME20K1202	KME21K0554	KME22D1499
			19DT25	KME20K1304	KME21K0662	KME22D1964
			19DT33	KME20K1552	KME21K0728	KME22D1965
			19DT45	KME20K1658	KME21K0830	KME22D2545
			19ET05	KME20K2057	KME21K1041	KME22D2547
			19ET15	KME20K2187	KME21K1044	KME22D2548
			19ET22	KME20K3161	KME21K1331	KME22D2640
			19ET28	KME20K3355	KME21K1413	KME22D2738
			19ET45	KME20L0057	KME21K1503	KME22D2847
			19ET60	KME20L0059	KME21K1517	KME22D2848
			19FT14	KME20L0490	KME21K1556	KME22F2400
			19FT24	KME20L0623	KME21K1620	KME22F2402
			19FT30	KME20L0943	KME21K1845	KME22F2868
			19FT41	KME20L0944	KME21K2355	KME22F2869
			19GT05	KME20L1091	KME21K2356	KME22G1064
			19GT10	KME20L1122	KME21L0093	KME22G1066
			19GT20	KME20L1124	KME21L0094	KME22G1068
			19GT59	KME20L1475	KME21L0768	KME22G1069
			19GT65	KME20L1844	KME21L1053	KME22G1089
			19HT40	KME20L2039	KME21L1054	KME22G1629
			19HT53	KME20L2258	KME21L1216	KME22G1632
			19IT25	KME20L2827	KME21L1217	KME22G1633
			19IT32	KME20M0378	KME21L1277	KME22M2258
19IT39	KME20M0379	KME21L1966	KME22M2534			
19IT44	KME20M0380	KME21L2129	KME22M2940			
19IT54	KME20M0498	KME21L2690	KME22M3055			
19JT22	KME20M1105	KME21L2811	KME23A2908			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
349779	504570	TUBE, ENDOTRACH CUFF W/STYLETTE 7.0MM (10/BX)	19JT57	KME20M1153	KME21M1888	KME23A3001
			19JT64	KME20M1156	KME21M2117	KME23A3199
			19JT70	KME21A0228	KME21M2320	KME23B0009
			19KT17	KME21A0273	KME21M2562	KME23B0377
			19KT25	KME21A0304	KME21M2845	KME23B1200
			19KT48	KME21A0305	KME21M3200	KME23B1201
			19KT55	KME21A1769	KME22A0150	KME23B1202
			19LT04	KME21A2227	KME22A0245	KME23B1311
			19LT10	KME21A2228	KME22A0390	KME23B1658
			19LT32	KME21A2229	KME22A0698	KME23B1659
			19LT38	KME21A2263	KME22A0797	KME23B1660
			19LT50	KME21A3119	KME22A0798	KME23B1661
			20AT09	KME21A3120	KME22A0990	KME23B3075
			20AT14	KME21A3121	KME22A1135	KME23B3146
			20AT25	KME21B0261	KME22A1161	KME23C1395
			20AT32	KME21B0338	KME22A1257	KME23C1538
			20AT44	KME21B0559	KME22A3253	KME23C1562
			20AT53	KME21B0750	KME22A3254	KME23C2218
			20AT57	KME21B1073	KME22A3312	KME23C2308
			20BG12	KME21B1751	KME22B0301	KME23C2438
			20BT05	KME21B2107	KME22B0354	KME23C3007
			20BT23	KME21B2153		
			349780	504575	TUBE, ENDOTRACH CUFF W/STYLETTE 7.5MM (10/BX) TELFLX	18FG05
18FG13	19IT06	KME20M2485				KME22A1786
18FG30	19IT14	KME20M2556				KME22A1877
18GG02	19IT25	KME21A0084				KME22A1909
18GG08	19IT29	KME21A0303				KME22A1912
18GG12	19IT39	KME21A0671				KME22A2057
18GG17	19IT44	KME21A0672				KME22A3084
18GG20	19JT22	KME21A0673				KME22A3086
18GG28	19JT57	KME21A0674				KME22A3087
18HG01	19JT70	KME21A0812				KME22A3151
18HG12	19KT08	KME21A1041				KME22A3270
18HG22	19KT25	KME21A1464				KME22B0570
18HG27	19KT48	KME21A1605				KME22B0649
18HG35	19LT04	KME21A2799				KME22B0850
18HG38	19LT10	KME21A2815				KME22B1220
18IG09	19LT23	KME21B0753				KME22B1407
18IG17	19LT50	KME21B0814				KME22B2126
18IG27	20AT09	KME21B0888				KME22B2127
18IG30	20AT14	KME21B0889				KME22B2128
18JG14	20AT25	KME21B0931				KME22B2162
18JG17	20AT32	KME21B1240				KME22B2333
18KG07	20AT37	KME21B2461				KME22B2471
18KG09	20AT46	KME21B2462				KME22B2472
18KG18	20AT53	KME21B2465	KME22B2497			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
349780	504575	TUBE, ENDOTRACH CUFF W/STYLETTE 7.5MM (10/BX) TELFLX	18KG29	20CT05	KME21B2716	KME22B2703
			18KG41	20CT12	KME21B2717	KME22B2903
			18KG44	20CT27	KME21C0700	KME22C0243
			18LG02	20CT36	KME21C0972	KME22C1859
			18LG05	20CT44	KME21C1203	KME22C2294
			18LG12	20DT09	KME21C1336	KME22C2402
			18LG22	20DT37	KME21C1435	KME22C2578
			18LG31	20ET04	KME21C1893	KME22D1474
			19AG05	20ET16	KME21C2565	KME22D2638
			19AG18	20ET37	KME21C2610	KME22D3063
			19AG24	20ET56	KME21G1164	KME22D3064
			19AG27	20FT61	KME21G1165	KME22D3677
			19AG31	20FT65	KME21G1356	KME22G1070
			19AG34	20GT06	KME21G2407	KME22G1071
			19BG01	20GT31	KME21J0947	KME22H2369
			19BG07	20GT38	KME21J1716	KME22H2372
			19CT20	20GT49	KME21J1841	KME22H2630
			19CT32	KME20H0980	KME21K0659	KME22H2905
			19CT37	KME20H0981	KME21K1347	KME22J0801
			19CT49	KME20H1654	KME21K1765	KME22J0802
			19CT55	KME20H1794	KME21K1844	KME22J2720
			19CT58	KME20H1848	KME21K2349	KME22J2721
			19CT65	KME20J0061	KME21K2974	KME22J3353
			19DT38	KME20J0157	KME21L0634	KME22K0311
			19DT45	KME20J0322	KME21L0876	KME22K0313
			19ET05	KME20J0323	KME21L1024	KME22M1776
			19ET22	KME20J0576	KME21L1392	KME22M1819
			19ET28	KME20J0829	KME21L2221	KME22M2535
			19ET45	KME20J0831	KME21L2348	KME22M2755
			19ET50	KME20J0834	KME21L2765	KME22M2756
			19ET54	KME20J1114	KME21M0086	KME22M2874
			19ET60	KME20J2544	KME21M0286	KME22M3030
			19FT01	KME20K0769	KME21M0363	KME22M3114
			19FT14	KME20K0770	KME21M0546	KME23A1028
			19FT24	KME20K1150	KME21M0962	KME23A3200
			19FT35	KME20K1204	KME21M1889	KME23B0861
			19FT41	KME20K1665	KME21M2037	KME23B0900
			19FT49	KME20K2188	KME21M2325	KME23B1023
			19GT01	KME20K2920	KME21M2564	KME23B1898
			19GT05	KME20K3016	KME21M2848	KME23B2846
19GT10	KME20K3018	KME22A0204	KME23B2919			
19GT28	KME20K3159	KME22A0546	KME23C1269			
19GT38	KME20K3575	KME22A1057	KME23C1323			
19GT59	KME20L1476	KME22A1198	KME23C1324			
19HT05	KME20L2374	KME22A1514	KME23C1396			
19HT10	KME20M0383	KME22A1588	KME23C1568			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
349780	504575	TUBE, ENDOTRACH CUFF W/STYLETTE 7.5MM (10/BX) TELFLX	19HT61	KME20M1151	KME22A1725	KME23C1790
			19HT70	KME20M2051	KME22A1727	KME23C2992
			19HT75	KME20M2344		
349781	504580	TUBE, ENDOTRACH CUFF W/STYLETTE 8.0MM (10/BX) TELFLX	18FG03	19JT42	KME20L1277	KME22A3271
			18FG05	19JT49	KME20L1479	KME22A3272
			18FG13	19KT17	KME20L2302	KME22A3283
			18GG14	19LT04	KME20L2828	KME22A3284
			18HG22	19LT17	KME20M1752	KME22A3285
			18HG38	20AT09	KME20M1753	KME22B1406
			18IG09	20AT14	KME21A0381	KME22B1996
			18JG06	20AT25	KME21A0820	KME22B2164
			18JG14	20AT32	KME21A0823	KME22B2329
			18JG26	20AT37	KME21A1291	KME22B2330
			18KG09	20AT46	KME21A1604	KME22B2496
			18KG15	20AT53	KME21B0339	KME22B2504
			18KG23	20BT05	KME21B1241	KME22B2595
			18LG05	20BT23	KME21B1855	KME22C0247
			18LG12	20BT47	KME21B2108	KME22C0865
			18LG18	20CT05	KME21B2154	KME22C1074
			19AG05	20CT12	KME21B2697	KME22C1321
			19AG31	20CT27	KME21C0377	KME22C1322
			19BG01	20CT36	KME21C0623	KME22H3393
			19BG07	20CT44	KME21C2407	KME22J1538
			19CT32	20DT14	KME21C2738	KME22J1713
			19CT41	20DT37	KME21C3448	KME22J3328
			19CT55	20ET16	KME21D1784	KME22K0537
			19DT38	20ET21	KME21E1286	KME22K0602
			19DT45	20GT06	KME21F1330	KME22K0744
			19DT52	20GT25	KME21J0409	KME22K0760
			19ET05	KME20H0590	KME21K2402	KME22L1503
			19ET32	KME20H0973	KME21L0372	KME22L1515
			19ET60	KME20H0974	KME21L1534	KME22M2049
			19ET72	KME20H1289	KME21L1899	KME23A0707
			19ET77	KME20H1290	KME21L2812	KME23B0978
			19FT24	KME20H2924	KME21M2036	KME23B1899
			19GT10	KME20H2926	KME21M2053	KME23B1900
19GT20	KME20K0449	KME21M2870	KME23B2022			
19GT38	KME20K0450	KME22A0395	KME23C0595			
19GT59	KME20K1414	KME22A1237	KME23C0968			
19HT05	KME20K2622	KME22A3025	KME23C3126			
		19IT06				
349782	504585	TUBE, ENDOTRACH CUFF 8.5MM (10/BX) TELFLX	18FG13	19FT14	KME20L0624	KME21M2519
			18FG15	19FT24	KME20L1482	KME22B0987
			18FG25	19GT52	KME20L2303	KME22B1754
			18GG20	19HT61	KME20L2832	KME22B2155
			18GG28	19HT80	KME21A0891	KME22C0679

MMS #	MFG Catalog #	Description	Affected Lot(s)			
349782	504585	TUBE, ENDOTRACH CUFF 8.5MM (10/BX) TELFLX	18HG12	19IT39	KME21A1043	KME22C0866
			18HG19	19KT17	KME21A2802	KME22C2001
			18IG06	19KT48	KME21A2803	KME22C2381
			18IG27	19LT50	KME21B0933	KME22D0823
			18JG02	20AT32	KME21B1250	KME22F1553
			18JG23	20AT53	KME21B1251	KME22F2230
			18JG26	20BT05	KME21B2029	KME22H0145
			18KG18	20CT12	KME21C0506	KME22J0252
			18KG23	20CT44	KME21C2107	KME22J0259
			19AG05	20DT09	KME21C2339	KME22J0470
			19AG31	20DT37	KME21C2341	KME22J2352
			19AG37	20GG27	KME21C2553	KME22L2186
			19CT17	20GG45	KME21D0210	KME22L2295
			19CT58	20GT38	KME21E1295	KME22L2804
			19CT71	KME20H0532	KME21K1874	KME23B2023
			19DT02	KME20J2404	KME21M1728	KME23C1449
			19ET32	KME20J2405	KME21M2447	KME23C1450
19ET45	KME20K2625					
349783	504590	TUBE, ENDOTRACH CUFF 9.0MM (10/BX) TELFLX	18FG18	19AG27	20BG16	KME21C2105
			18FG20	19BG04	20BT05	KME21E1114
			18GG17	19DT38	20BT47	KME22B1404
			18GT32	19DT45	20CT36	KME22B1405
			18HG12	19ET32	20CT44	KME22B1755
			18HG27	19ET72	KME20H0238	KME22C0915
			18HG38	19FT14	KME20H0239	KME22C1943
			18IG15	19GT20	KME20H0591	KME22C2372
			18IG30	19HT40	KME20J2420	KME22D1717
			18KG03	19IT39	KME20L0626	KME22D1718
			18KG07	19IT58	KME20M0856	KME22G1443
			18KG35	19KT17	KME20M2091	KME23A2760
			18LG18	19KT21	KME21C0299	KME23C0137
19AG24	19LT50	KME21C2104				
401083	150020	TUBE, ENDOTRACH UNCUFF W/STYL 2.0MM (10/BX) TELFLX	18FG01	18HG22	18IG02	18JG17
			18FG05	18HG27	18JG04	18JG21
			18HG07	18HG38	18JG06	18JG23
401084	150025	TUBE, ENDOTRACH UNCUFF 2.5MM (10/BX) TELFLX	18GT32	18HG29	18IG17	18IG24
			18HG07	18HG38	18IG21	18IG27
			18HG16	18IG06		
417647	112082075	TUBE, ENDOTRACH 7.5MM (PARKLAND) (10/BX) TELFLX	18FG16	19KT41	KME21C2026	KME22A2602
			18FG29	19KT48	KME21C2051	KME22A2603
			18GG13	19KT49	KME21C2066	KME22A2898
			18GG15	19LT07	KME21C2195	KME22A2937
			18GG16	19LT23	KME21C2499	KME22A3082
			18GG22	19LT31	KME21C2615	KME22A3196
			18GG24	19LT35	KME21C2680	KME22A3219
18GG26	19LT54	KME21C2766	KME22A3222			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
417647	112082075	TUBE, ENDOTRACH 7.5MM (PARKLAND) (10/BX) TELFLX	18GG27	20AT11	KME21C2908	KME22A3238
			18GT16	20AT14	KME21C2964	KME22B0266
			18GT19	20AT20	KME21C3394	KME22C2126
			18GT32	20AT22	KME21D0029	KME22C2350
			18GT34	20AT26	KME21D0058	KME22C2437
			18HG14	20AT32	KME21D0186	KME22C2439
			18HG21	20AT38	KME21D0586	KME22C2604
			18HG23	20AT47	KME21D0648	KME22C2756
			18HG24	20AT54	KME21D0775	KME22C2855
			18HG27	20AT58	KME21D1103	KME22C3145
			18HG33	20BT01	KME21D1400	KME22D0324
			18HT15	20BT09	KME21D1423	KME22D0464
			18HT18	20BT24	KME21D1705	KME22D1176
			18HT22	20BT25	KME21D1780	KME22D1481
			18HT24	20BT31	KME21D1923	KME22D1660
			18HT31	20BT37	KME21D1956	KME22D1688
			18IG08	20BT39	KME21D2014	KME22D2635
			18IG22	20BT50	KME21D2018	KME22D2853
			18IG23	20BT51	KME21D2262	KME22D3116
			18IG25	20CT42	KME21E0021	KME22D3236
			18IG29	20CT43	KME21E0288	KME22D3419
			18IT04	20CT51	KME21E0469	KME22D3420
			18IT18	20CT53	KME21E0498	KME22D3711
			18IT23	20CT54	KME21E1150	KME22E0179
			18IT25	20DT01	KME21E1178	KME22E0319
			18JG01	20DT02	KME21E1186	KME22E0321
			18JG02	20DT05	KME21E1264	KME22E0322
			18JG04	20DT08	KME21F0184	KME22E0641
			18JG07	20DT20	KME21F0351	KME22E0642
			18JG10	20DT32	KME21F0385	KME22E1861
			18JG24	20DT37	KME21F0557	KME22E2010
			18JG32	20ET12	KME21F0594	KME22E2011
			18JG38	20ET34	KME21F0725	KME22E2647
			18JT04	20ET37	KME21F1292	KME22E3002
			18JT06	20ET39	KME21F1380	KME22F0185
			18JT22	20FT27	KME21F1529	KME22F0317
			18KG05	20FT35	KME21F1583	KME22F0775
			18KG35	20FT36	KME21F1607	KME22F0776
			18KT04	20FT39	KME21F1696	KME22F1265
			18KT11	20FT43	KME21F1751	KME22F2028
18KT21	20FT49	KME21F2182	KME22F2237			
18LG19	20FT53	KME21G0114	KME22F2721			
18LG38	20FT55	KME21G0169	KME22H1514			
18LG40	20FT58	KME21G0898	KME22H1855			
18LT07	20FT61	KME21G1942	KME22H1898			
19AG21	20FT71	KME21H0829	KME22H1953			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
417647	112082075	TUBE, ENDOTRACH 7.5MM (PARKLAND) (10/BX) TELFLX	19AG36	20GG41	KME21J0200	KME22H1987
			19AT20	20GG42	KME21J0359	KME22H1993
			19BG01	20GT04	KME21J0418	KME22H2002
			19BG04	20GT10	KME21J0562	KME22H2172
			19BG06	20GT11	KME21J0647	KME22H2247
			19BG37	20GT14	KME21J0671	KME22H2300
			19BT04	20GT18	KME21J0726	KME22H2359
			19BT12	20GT25	KME21J0763	KME22H2417
			19BT18	20GT31	KME21J0854	KME22H2448
			19CG04	20GT41	KME21J1212	KME22H2607
			19CG06	20GT51	KME21J1709	KME22H2865
			19CT03	20GT57	KME21J1842	KME22J0150
			19CT04	KME20H2118	KME21J2080	KME22J0511
			19CT20	KME20H2346	KME21K0305	KME22J1107
			19CT26	KME20H2728	KME21K0986	KME22J1258
			19CT34	KME20H3115	KME21K1712	KME22J1490
			19CT42	KME20H3120	KME21K1855	KME22J3017
			19CT50	KME20J0048	KME21K2171	KME22J3143
			19CT57	KME20J1051	KME21K2475	KME22J3344
			19CT64	KME20J1455	KME21K2531	KME22K0071
			19CT69	KME20J1469	KME21K2745	KME22K0116
			19DT10	KME20J1657	KME21L0628	KME22K0138
			19DT34	KME20J1971	KME21L0880	KME22K0277
			19DT36	KME20J1997	KME21L1025	KME22K0364
			19DT42	KME20J2159	KME21L1118	KME22K0430
			19DT45	KME20J2538	KME21L1443	KME22K0580
			19DT49	KME20J2789	KME21L1581	KME22K0581
			19DT50	KME20J2836	KME21L1663	KME22L1286
			19ET02	KME20J2918	KME21L1698	KME22L1559
			19ET06	KME20J2944	KME21L1992	KME22L1716
			19ET08	KME20J2988	KME21L2222	KME22L2285
			19ET12	KME20J2992	KME21L2455	KME22M1576
			19ET22	KME20J3022	KME21L2775	KME22M1645
			19ET29	KME20J3165	KME21M0193	KME22M1646
			19ET34	KME20J3284	KME21M0310	KME22M1682
			19ET36	KME20K0052	KME21M0314	KME22M1688
			19ET49	KME20K0181	KME21M0382	KME22M1717
			19ET67	KME20K1113	KME21M0500	KME22M1730
			19ET70	KME20K3044	KME21M0690	KME22M1764
			19ET72	KME20K3577	KME21M0802	KME22M1963
			19ET77	KME20L0484	KME21M0852	KME22M2037
			19ET82	KME20L0755	KME21M0930	KME22M2165
			19FT04	KME20L1092	KME21M1006	KME22M2260
			19FT25	KME20L1156	KME21M1065	KME22M2454
			19FT33	KME20L1298	KME21M1070	KME22M2521
19GT02	KME20L1299	KME21M1085	KME22M2556			

MMS #	MFG Catalog #	Description	Affected Lot(s)			
417647	112082075	TUBE, ENDOTRACH 7.5MM (PARKLAND) (10/BX) TELFLX	19GT14	KME20L1704	KME21M1245	KME22M2565
			19GT16	KME20L1924	KME21M1268	KME22M2630
			19GT17	KME20L1983	KME21M1385	KME22M2637
			19GT25	KME20L2605	KME21M1503	KME22M2649
			19GT29	KME20L2798	KME21M1583	KME22M2722
			19GT30	KME20M0034	KME21M1613	KME22M2726
			19GT31	KME20M0321	KME21M1638	KME22M2773
			19GT32	KME20M1067	KME21M1749	KME22M2794
			19GT35	KME20M1507	KME21M2135	KME22M2875
			19GT40	KME20M1890	KME21M2478	KME22M2897
			19GT64	KME20M2212	KME21M2499	KME22M2908
			19HT03	KME21A3004	KME21M2646	KME22M2928
			19HT09	KME21A3259	KME21M2699	KME22M2943
			19HT15	KME21B0179	KME21M2708	KME22M3002
			19HT22	KME21B0823	KME21M2745	KME22M3017
			19HT71	KME21B2480	KME21M2790	KME22M3046
			19IT12	KME21B2482	KME21M2871	KME22M3088
			19IT13	KME21C0239	KME21M2967	KME22M3100
			19IT41	KME21C0264	KME22A0545	KME22M3142
			19IT42	KME21C0614	KME22A0854	KME22M3188
			19IT43	KME21C0615	KME22A0940	KME23A2932
			19IT45	KME21C0661	KME22A1136	KME23A3025
			19JT11	KME21C0702	KME22A1224	KME23B0026
			19JT14	KME21C0730	KME22A1367	KME23B0240
			19JT59	KME21C0826	KME22A1392	KME23B0703
			19JT61	KME21C0855	KME22A1754	KME23B1328
			19KT06	KME21C1093	KME22A1990	KME23B1570
			19KT09	KME21C1223	KME22A1993	KME23B1672
			19KT13	KME21C1323	KME22A2115	KME23B1673
			19KT14	KME21C1399	KME22A2129	KME23B1718
			19KT18	KME21C1601	KME22A2315	KME23B2081
			19KT19	KME21C1704	KME22A2351	KME23B2098
			19KT20	KME21C1746	KME22A2364	KME23C1322
19KT26	KME21C1944	KME22A2414	KME23C1441			
19KT35	KME21C2008	KME22A2537	KME23C1689			
19KT38						
420192	170080	TUBE, ENDOTRACH CUFF W/STYL 8.0MM (10/BX)	18FG09	19AG27	20AT53	KME21B1752
			18FG20	19AG31	20BT05	KME21B2156
			18FT32	19BG01	20BT47	KME21C0465
			18GG05	19CT41	20CT05	KME21C0533
			18GG08	19DT38	20CT12	KME21C2826
			18GG14	19ET38	20CT27	KME21D1204
			18HG07	19ET45	20CT36	KME21E1287
			18HG12	19ET50	20CT51	KME21F1329
			18HG22	19FT24	20DT09	KME21J0410
			18HG27	19GT59	20FT44	KME21L0274

MMS #	MFG Catalog #	Description	Affected Lot(s)			
420192	170080	TUBE, ENDOTRACH CUFF W/STYL 8.0MM (10/BX)	18HG35	19IT06	20GT25	KME21L0813
			18IG02	19IT25	KME20H0248	KME21L0851
			18IG09	19IT54	KME20H2019	KME22E2351
			18IG15	19JT42	KME20L0396	KME22E2354
			18JG04	19JT49	KME20L1278	KME22E2794
			18JG06	19JT57	KME20L1279	KME22E2796
			18JG14	19KT25	KME20L2830	KME22F0303
			18JG26	19LT04	KME20M0847	KME22J0007
			18LG12	20AT14	KME21A0386	KME22J3345
			18LG18	20AT25	KME21A2800	KME22L1514
			19AG13	20AT32	KME21A2801	KME23B2041
			19AG24	20AT37	KME21B1245	
			421461	170135	TUBE, ENDOTRACH UNCUFF W/STYL 3.5MM (10/BX)	18KG03
18KG07	19HT61	KME20K1629				KME22A1510
18KG18	19HT66	KME20K2547				KME22B0538
18LG22	19IT54	KME20K2583				KME22B0726
19AG18	19JT22	KME20K3590				KME22B1029
19AG21	19JT29	KME20L0632				KME22B2334
19AG24	20AT32	KME20L1486				KME22C1825
19AG37	20AT53	KME20L2401				KME22C2261
19BG04	20BT05	KME20M0811				KME22D0351
19BG12	20BT23	KME20M3276				KME22E2104
19BG40	20BT47	KME20M3277				KME22E2263
19CT49	20CT12	KME21A1898				KME22E2264
19CT58	20CT36	KME21A2813				KME22F2163
19DT06	20DT09	KME21B1686				KME22F2418
19DT20	20DT37	KME21C0018				KME22H0050
19DT38	20ET16	KME21D1850				KME22H0350
19ET22	20ET21	KME21E0464				KME22H2625
19ET28	20ET28	KME21K1250				KME22J0996
19GT45	KME20G0242	KME21L1893				KME22K2777
19GT59	KME20H1846	KME21L1894				KME23A2592
19GT65	KME20J0328	KME21L2343	KME23A3201			
19HT05						
421463	170145	TUBE, ENDOTRACH UNCUFF STR 4.5MM (10/BX)	18JG14	20CT12	KME20L1488	KME21J1149
			18LG15	20CT51	KME20L2377	KME21M1818
			19AG24	20DT29	KME21A1891	KME21M2952
			19BG12	20GG41	KME21C0019	KME22D1896
			19BG16	KME20H2929	KME21C1566	KME22E1744
			19CT49	KME20K1886	KME21C1974	KME22F0717
			19IT25	KME20K2586	KME21D0445	KME22G3114
20BT05						
421466	170055	TUBE, ENDOTRACH CUFF W/STYL STR 5.5MM (10/BX)	18FG18	20AT25	KME20K2587	KME21M0313
			18GG14	20AT32	KME20L1469	KME21M1832
			18GG31	20BT05	KME20L2378	KME22B0193
			18HG22	20BT37	KME21A0383	KME22D1787

MMS #	MFG Catalog #	Description	Affected Lot(s)			
421466	170055	TUBE, ENDOTRACH CUFF W/STYL STR 5.5MM (10/BX)	18KG23	20BT47	KME21A1901	KME22F1665
			19BG12	20ET21	KME21B2682	KME22G2296
			19ET72	KME20G0383	KME21C2708	KME22H3106
			19IT25	KME20K0936	KME21G0756	KME22K0295
421467	170060	TUBE, ENDOTRACH CUFF W/STYL STR 6.0MM (10/BX)	18GG31	19DT20	KME20H1653	KME21K0548
			18HG22	19ET22	KME20J1236	KME21K0549
			18HG32	19ET38	KME20J1237	KME21L0444
			18HG35	19ET66	KME20L1487	KME21M1392
			18IG09	19ET77	KME20L2379	KME22A1512
			18IG24	19FT18	KME20M0494	KME22A1513
			18IG27	19IT14	KME21A0384	KME22B1392
			18JG02	19IT25	KME21A1482	KME22D0595
			18JG17	19JT42	KME21C0293	KME22E0622
			18LG12	20DT29	KME21E1411	KME22E2793
			19AG37	20ET04	KME21G0757	KME22G2464
			19BG12	20FT34	KME21G2208	KME22K1055
			19CT49	20FT54	KME21G2209	KME23C2504
			19DT15	20GG41		
421468	170065	TUBE, ENDOTRACH CUFF W/STYL STR 6.5MM(10/BX)	18FG05	19CT71	20CT22	KME21A1812
			18FG25	19DT15	20DT37	KME21B1670
			18GG05	19ET60	20ET21	KME21B2499
			18GG10	19FT41	20GT14	KME21B2500
			18GG20	19GT38	KME20H0903	KME21K1127
			18GT32	19HT15	KME20H0904	KME21K1256
			18HG19	19HT61	KME20H2018	KME21K3029
			18HG32	19HT66	KME20J1239	KME21L0568
			18HG38	19IT25	KME20J1240	KME22B0064
			18IG09	19JT18	KME20K0938	KME22B0066
			18JG04	19JT22	KME20K2623	KME22D0757
			18JG06	19JT57	KME20L0487	KME22F1796
			18JG34	19JT70	KME20L1705	KME22G2468
			18KG35	19KT41	KME20M1578	KME22H1867
			19AG13	20AT25	KME20M1751	KME22J3351
			19AG37	20AT53	KME20M2709	KME22L0238
			19BG01	20BT27	KME21A0385	KME22L0239
			19CT58	20CT05		
421471	170070	TUBE, ENDOTRACH CUFF W/STYL STR 7.0MM (10/BX)	18FG09	19FT41	20ET16	KME21C1422
			18FG15	19GT05	20ET56	KME21C1592
			18FG17	19GT10	20FT10	KME21C2266
			18FG20	19GT52	20FT54	KME21C2504
			18FG21	19GT59	20GG36	KME21C2706
			18GG10	19GT65	20GG45	KME21D1043
			18GG28	19HT40	KME20H0245	KME21D1430
			18GT32	19HT66	KME20H1287	KME21E1290
			18HG07	19IT44	KME20H1288	KME21F1526
			18HG16	19IT54	KME20J0571	KME21G0541

MMS #	MFG Catalog #	Description	Affected Lot(s)			
421471	170070	TUBE, ENDOTRACH CUFF W/STYL STR 7.0MM (10/BX)	18HG35	19JT22	KME20J0572	KME21G0542
			18HG38	19JT42	KME20J0802	KME21G0543
			18IG09	19JT57	KME20J3033	KME21G1559
			18JG02	19JT64	KME20K2544	KME21H0608
			18JG17	19KT30	KME20L0053	KME21K0551
			18JG37	19KT41	KME20L0055	KME21K0552
			18KG23	19LT04	KME20L0942	KME21K1767
			18KG35	19LT17	KME20L2372	KME21K2469
			18LG05	20AT14	KME20M0381	KME21L1897
			18LG12	20AT25	KME20M0665	KME21L1898
			18LG18	20AT44	KME20M1103	KME22E2554
			19AG13	20AT53	KME20M1754	KME22E2795
			19AG31	20BT05	KME20M1755	KME22F0302
			19AG34	20BT23	KME21A0893	KME22F0755
			19AG37	20BT47	KME21A2264	KME22F0756
			19BG04	20CT05	KME21B1074	KME22F0837
			19CT17	20CT12	KME21B1075	KME22F0838
			19CT41	20CT27	KME21C0463	KME22F1638
			19CT49	20CT51	KME21C0464	KME22G1630
			19CT55	20DT09	KME21C1027	KME22H2376
19FT24	20DT37					
421477	170075	TUBE, ENDOTRACH CUFF W/STYL 7.5MM (10/BX)	18FG09	19BG07	20BT47	KME21C1895
			18FG15	19CT37	20CT22	KME21C2855
			18FG17	19CT41	20CT36	KME21C3294
			18FG22	19CT49	20CT51	KME21C3295
			18FG30	19CT55	20DT09	KME21D1785
			18FT32	19FT14	20DT37	KME21F1326
			18GG02	19FT18	20ET16	KME21F1327
			18GG05	19FT35	20ET37	KME21G0544
			18GG28	19GT10	20FT54	KME21K1410
			18GT32	19GT52	20FT61	KME21L2179
			18HG12	19GT59	KME20H0977	KME21L2346
			18HG16	19GT65	KME20H0979	KME21L2347
			18HG22	19HT10	KME20J0575	KME21M0195
			18HG27	19HT35	KME20J2715	KME21M0287
			18HG38	19HT48	KME20K1352	KME21M0869
			18IG06	19IT06	KME20L1132	KME21M0965
			18IG15	19IT14	KME20M0033	KME21M1267
			18IG17	19IT25	KME20M0382	KME21M1748
			18IG27	19IT44	KME20M1149	KME22A2896
			18JG02	19JT22	KME20M1225	KME22E2405
			18JG11	19JT49	KME20M1227	KME22E2911
			18JG26	19JT70	KME20M1537	KME22E2954
			18KG07	19KT25	KME20M2881	KME22E2987
			18KG18	19LT04	KME20M2882	KME22F0181
			18KG44	19LT32	KME21A0824	KME22F2658

MMS #	MFG Catalog #	Description	Affected Lot(s)			
421477	170075	TUBE, ENDOTRACH CUFF W/STYL 7.5MM (10/BX)	18LG12	20AT14	KME21A1161	KME22F2659
			18LG22	20AT25	KME21A1813	KME22H2852
			19AG13	20AT32	KME21B0813	KME22H3104
			19AG27	20AT44	KME21C0212	KME22J0009
			19AG34	20AT53	KME21C0213	KME22K0312
			19AG37	20BG12	KME21C1594	KME22L2286
421508	170085	TUBE, ENDOTRACH CUFF W/STYL STR 8.5MM (10/BX)	18FG09	19ET38	20BT05	KME21C2385
			18GG20	19ET60	20BT23	KME21E1296
			18GG28	19HT61	20CT51	KME21K1876
			18HG07	19HT66	20GG36	KME21L0573
			18JG11	19IT44	KME20L0398	KME22F0484
			18KG18	19KT21	KME20L2376	KME22F1069
			19AG37	19KT55	KME20M3247	KME22J0471
			19CT32	19LT32	KME21A1894	KME22J2353
421509	170090	TUBE, ENDOTRACH CUFF W/STYL STR 9.0MM (10/BX)	18FG09	19GT52	20AT25	KME22E2491
			18HG27	19GT65	KME20L0625	KME22E2912
			18IG12	19HT35	KME21E1232	KME22M0849
			18LG18			
465321	112480025	TUBE, MURPHY ENDOTRACH CUFF 2.5MM (10/BX) TELFLX	18FG31	19IT59	KME21A0348	KME22E2191
			18GT32	19LT08	KME21B0087	KME22G1878
			18IG03	19LT31	KME21B0465	KME22G2593
			18IG29	20BG15	KME21C0163	KME22H0532
			18KG08	20BT01	KME21C1805	KME22H1185
			18LG03	20CT07	KME21C2447	KME22H2178
			19AG08	20DT24	KME21C2688	KME22H2180
			19AG38	20ET53	KME21C3244	KME22J1309
			19BG12	20GT38	KME21K2877	KME22J2321
			19BG35	KME20H2416	KME22A1408	KME22K0879
			19DT13	KME20H2417	KME22A2390	KME22K1786
			19DT36	KME20J1068	KME22A2545	KME22K1787
			19ET51	KME20K2143	KME22B0845	KME23B0559
			19FT12	KME20K2598	KME22B2017	KME23B1263
			19FT28	KME20L0449	KME22B2848	KME23B1265
			19GT02	KME20L1346	KME22C2802	KME23C1148
19HT05	KME20M0127	KME22E0598	KME23C2363			
19HT79	KME20M1797	KME22E2130	KME23C2886			
465322	112480030	TUBE, MURPHY ENDOTRACH CUFF 3.0MM (10/BX) TELFLX	18GG16	19GT02	20GT38	KME22C2669
			18HG10	19IT38	KME20J3096	KME22C2687
			18HG25	19JT21	KME20K2145	KME22D3223
			18IG03	19JT55	KME20L1345	KME22E1149
			19AG23	19KT20	KME20M1508	KME22E2287
			19AG38	19KT37	KME20M1802	KME22G1233
			19BG12	19LT31	KME21B0621	KME22G1879
			19BG16	20AT40	KME21B2379	KME23B0133
	19CT42	20AT54	KME21C0168	KME23B0560		

MMS #	MFG Catalog #	Description	Affected Lot(s)			
465322	112480030	TUBE, MURPHY ENDOTRACH CUFF 3.0MM (10/BX) TELFLX	19CT64	20BT51	KME22A1495	KME23B0717
			19DT36	20DT01	KME22A2394	KME23B1351
			19ET09	20DT20	KME22A2396	KME23B2197
			19ET31	20DT24	KME22B2938	KME23C0078
			19ET53	20ET30	KME22C2668	KME23C0079
626214	1700100	TUBE, ENDOTRACH CUFFED STR 10.MM	18IG06	18KG44	KME22J2841	
678664	112080095	TUBE, MAGILL ENDOTRACH NASAL 9.5MM (10/BX)	18IT25	19HT79	20DT19	KME22C0673
			18JG13	19KT08	KME20K2310	KME22F1172
			19CT66	20BT13	KME20L1153	
1132959	170140	TUBE, ENDOTRACH UNCUFF STR 4.0MM (10/BX)	18JG17	19ET72	KME20L1491	KME22C2262
			18JG32	19ET77	KME20L2402	KME22C3452
			18JG34	19FT24	KME21A1159	KME22E1239
			18KG29	19JT29	KME21A1900	KME22F0049
			18LG02	20BG16	KME21A2810	KME22F1951
			18LG15	20CT36	KME21B0572	KME22H0913
			18LG18	20DT37	KME21B1688	KME22H0967
			19AG09	20ET28	KME21C0024	KME22H2910
			19AG21	20GG45	KME21C0493	KME22J1137
			19BG12	20GT14	KME21C2999	KME22K0018
			19CT32	KME20H2017	KME21L1921	KME22K0337
			19DT06	KME20K0772	KME21L1922	KME22L1356
			19DT20	KME20K2549	KME22B1077	KME23B1954
19ET38	KME20K2584					
1137457	150025	AIRWAY, LMA MASK REUSE SZ2.5	18GT32	18HG29	18IG17	18IG24
			18HG07	18HG38	18IG21	18IG27
			18HG16	18IG06		
1137459	150030	AIRWAY, LMA MASK REUSE SZ3	18IG15	18IG27	18JG02	18JG14
			18IG21			
1137462	150020	AIRWAY, LMA MASK REUSE SZ2	18FG01	18HG22	18IG02	18JG17
			18FG05	18HG27	18JG04	18JG21
			18HG07	18HG38	18JG06	18JG23
1201753	150040	AIRWAY, LMA MASK RUSBL SZ4	18IG12	18IG21	18IG27	18JG14
			18IG15			
1210844	170050	TUBE, ENDOTRACH PRELOADED UNCUFFED SZ5.0 (10/BX)	18FG05	19DT33	KME20J1238	KME21K1248
			18GG08	19ET28	KME20K2551	KME21M2648
			18GT32	19FT18	KME20L1089	KME22B1391
			18HG38	19IT25	KME20L2371	KME22D1768
			18JG02	19KT48	KME21A2064	KME22D1769
			18JG11	20AT14	KME21B2498	KME22E2745
			18KG15	20DT09	KME21C0488	KME22F1794
			18LG22	20ET04	KME21F1552	KME22H1168
			19AG24	20GT31	KME21G1708	KME22K1248
19DT06	KME20J1235					
1210846	170120	TUBE, ENDOTRACH PRELOADED UNCUFFED SZ2.0 (10/BX)	18JG23	19CT58	KME21B1671	KME22D1788
			18JG26	19CT71	KME21B2506	KME22E1489
			18JG34	19DT06	KME21C2475	KME22F1497

MMS #	MFG Catalog #	Description	Affected Lot(s)			
1210846	170120	TUBE, ENDOTRACH PRELOADED UNCUFFED SZ2.0 (10/BX)	18KG26	19DT20	KME21D0447	KME22F1500
			18LG05	19ET77	KME21K1723	KME22F1669
			19AG18	19FT14	KME22A0056	KME22F1670
			19AG21	KME20H0975	KME22B0179	KME22G3099
			19AG24	KME20H0976	KME22B0180	KME22H3107
			19AG37	KME20J0324	KME22B0181	KME22H3275
			19BG01	KME20K0771	KME22B0464	KME22J2354
			19BG04	KME20K2545	KME22C2255	KME22J2355
			19BG07	KME20K2568	KME22C2256	KME22J2356
			19BG12	KME20L0627	KME22D1638	KME23A2360
			19BG16	KME21A1886	KME22D1639	

McKesson Medical-Surgical Inc.
 Device Recall Reply Form: RC-2023-128
 Teleflex Medical Various Tracheal Tubes

June 19, 2023

Complete this reply form and return all pages immediately via email to MMSRecalls@McKesson.com or fax at **(866) 871-0270** should you have affected product.

To ensure timely credit to your account and support the completion of this notice, please respond within 30 days.

Date: _____ Ship to Acct Number: _____

Your Name: _____ Email Address: _____

Phone Number: _____ Fax Number: _____

Account Name: _____

Address: _____

City, State Zip: _____

I acknowledge that I DO HAVE product affected by this notification and have followed the instructions for return.

Qty	Unit of Measure	MMS #	MFG Catalog #	Description
		191450	112080060	TUBE, MAGILL ENDOTRACH NASAL 6.0MM (10/BX) TELFLX
		191305	112080070	TUBE, MAGILL ENDOTRACH NASAL 7.0MM (10/BX) TELFLX
		1137462	150020	AIRWAY, LMA MASK REUSE SZ2
		1137457	150025	AIRWAY, LMA MASK REUSE SZ2.5
		1210844	170050	TUBE, ENDOTRACH PRELOADED UNCUFFED SZ5.0 (10/BX)
		421467	170060	TUBE, ENDOTRACH CUFF W/STYL STR 6.0MM (10/BX)
		421468	170065	TUBE, ENDOTRACH CUFF W/STYL STR 6.5MM(10/BX)
		421471	170070	TUBE, ENDOTRACH CUFF W/STYL STR 7.0MM (10/BX)
		420192	170080	TUBE, ENDOTRACH CUFF W/STYL 8.0MM (10/BX)
		421508	170085	TUBE, ENDOTRACH CUFF W/STYL STR 8.5MM (10/BX)
		421509	170090	TUBE, ENDOTRACH CUFF W/STYL STR 9.0MM (10/BX)
		349776	504555	TUBE, ENDOTRACH CUFF 5.5MM (10/BX) TELFLX
		349777	504560	TUBE, ENDOTRACH CUFF W/STYLETTE 6.0MM (10/BX) TELFLX
		349781	504580	TUBE, ENDOTRACH CUFF W/STYLETTE 8.0MM (10/BX) TELFLX

Qty	Unit of Measure	MMS #	MFG Catalog #	Description
		349782	504585	TUBE, ENDOTRACH CUFF 8.5MM (10/BX) TELFLX
		349783	504590	TUBE, ENDOTRACH CUFF 9.0MM (10/BX) TELFLX
		349769	506525	TUBE, ENDOTRACH UNCUFF W/STYLETTE 2.5MM (10/BX) TELFLX
		349770	506530	TUBE, ENDOTRACH UNCUFF W/STYLETTE 3.0MM (10/BX) TELFLX
		349771	506535	TUBE, ENDOTRACH UNCUFF 3.5MM (10/BX) TELFLX
		349772	506540	TUBE, ENDOTRACH UNCUFF W/STYLETTE 4.0MM (10/BX) TELFLX
		349773	506545	TUBE, ENDOTRACH UNCUFF 4.5MM (10/BX) TELFLX
		626214	1700100	TUBE, ENDOTRACH CUFFED STR 10.MM (10/BX) TELFLX
		229807	100180030	TUBE, ENDOTRACH ORAL UNCUFF 3.0MM (10/BX) TELFLX
		229808	100180035	TUBE, ENDOTRACH ORAL UNCUFF 3.5MM (10/BX) TELFLX
		229809	100180040	TUBE, ENDOTRACH ORAL UNCUFF 4.0MM (10/BX) TELFLX
		229810	100180045	TUBE, ENDOTRACH ORAL UNCUFF 4.5MM (10/BX) TELFLX
		229811	100180050	TUBE, ENDOTRACH ORAL UNCUFF 5.0MM (10/BX) TELFLX
		229812	100180055	TUBE, ENDOTRACH ORAL UNCUFF 5.5MM (10/BX) TELFLX
		229813	100180060	TUBE, ENDOTRACH ORAL UNCUFF 6.0MM (10/BX) TELFLX
		229815	100180070	TUBE, ENDOTRACH ORAL UNCUFF 7.0MM (10/BX) TELFLX
		229787	100181030	TUBE, ENDOTRACH NASAL UNCUFF 3.0MM (10/BX) TELFLX
		229789	100181040	TUBE, ENDOTRACH NASAL UNCUFF 4.0MM (10/BX) TELFLX
		229790	100181045	TUBE, ENDOTRACH NASAL UNCUFF 4.5MM (10/BX) TELFLX
		229791	100181050	TUBE, ENDOTRACH NASAL UNCUFF 5.0MM (10/BX) TELFLX
		229792	100181055	TUBE, ENDOTRACH NASAL UNCUFF 5.5MM (10/BX) TELFLX
		229793	100181060	TUBE, ENDOTRACH NASAL UNCUFF 6.0MM (10/BX) TELFLX
		229794	100181065	TUBE, ENDOTRACH NASAL UNCUFF 6.5MM (10/BX) TELFLX
		189776	100380025	TUBE, MAGILL ENDOTRACH UNCUFF 2.5MM (10/BX) TELFLX

Qty	Unit of Measure	MMS #	MFG Catalog #	Description
		189777	100380030	TUBE, MAGILL ENDOTRACH UNCUFF 3.0MM (10/BX) TELFLX
		177517	100380035	TUBE, MAGILL ENDOTRACH UNCUFF 3.5MM (10/BX) TELFLX
		189778	100380040	TUBE, MAGILL ENDOTRACH UNCUFF 4.0MM (10/BX) TELFLX
		252280	100380045	TUBE, MAGILL ENDO UNCUFF 4.5MMTELFLX
		189780	100380050	TUBE, MAGILL ENDOTRACH UNCUFF 5.0MM (10/BX) TELFLX
		173546	100380055	TUBE, MAGILL ENDOTRACH UNCUFF 5.5MM (10/BX) TELFLX
		189781	100380060	TUBE, MAGILL ENDOTRACH UNCUFF 6.0MM (10/BX) TELFLX
		147586	100382025	TUBE, MURPHY ENDOTRACH UNCUFF 2.5MM (10/BX) TELFLX
		144479	100382030	TUBE, MURPHY ENDOTRACH UNCUFF 3.0MM (10/BX) TELFLX
		144480	100382035	TUBE, MURPHY ENDOTRACH UNCUFF 3.5MM (10/BX) TELFLX
		144481	100382040	TUBE, MURPHY ENDOTRACH UNCUFF 4.0MM (10/BX) TELFLX
		159587	100382045	TUBE, MURPHY ENDOTRACH UNCUFF 4.5MM (10/BX) TELFLX
		144482	100382050	TUBE, MURPHY ENDOTRACH UNCUFF 5.0MM (10/BX) TELFLX
		159586	100382055	TUBE, MURPHY ENDOTRACH UNCUFF 5.5MM (10/BX) TELFLX
		144483	100382060	TUBE, MURPHY ENDOTRACH UNCUFF 6.0MM (10/BX) TELFLX
		189783	100382065	TUBE, MURPHY ENDOTRACH UNCUFF 6.5MM (10/BX) TELFLX
		229795	111780035	TUBE, MURPHY ENDOTRACH CUFF 3.5MM (10/BX) TELFLX
		229797	111780045	TUBE, AGT MURPHY ENDOTRACH CUFF 4.5MM (10/BX) TELFLX
		229798	111780050	TUBE, AGT MURPHY ENDOTRACH CUFF 5.0MM (10/BX) TELFLX
		229799	111780055	TUBE, AGT MURPHY ENDOTRACH CUFF 5.5MM (10/BX) TELFLX
		229801	111780065	TUBE, AGT MURPHY ENDOTRACH CUFF 6.5MM (10/BX) TELFLX
		229802	111780070	TUBE, AGT MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX
		229804	111780080	TUBE, AGT MURPHY ENDOTRACH CUFF 8.0MM (10/BX) TELFLX
		229805	111780085	TUBE, MURPHY ENDOTRACH CUFF 8.5MM (10/BX) TELFLX
		229777	111781035	TUBE, MURPHY ENDOTRACH CUFF 3.5MM (10/BX) TELFLX

Qty	Unit of Measure	MMS #	MFG Catalog #	Description
		229778	111781040	TUBE, NASAL MURPHY ENDOTRACH CUFF 40MM (10/BX) TELFLX
		229781	111781055	TUBE, MURPHY ENDOTRACH CUFF 5.5MM (10/BX) TELFLX
		229782	111781060	TUBE, MURPHY ENDOTRACH CUFF 6.0MM (10/BX) TELFLX
		229783	111781065	TUBE, MURPHY ENDOTRACH CUFF 6.5MM (10/BX) TELFLX
		229785	111781075	TUBE, MURPHY ENDOTRACH CUFF 7.5MM (10/BX) TELFLX
		229786	111781080	TUBE, MURPHY NASAL CUFF 8.0MM (10/BX) TELFLX
		242472	112080050	TUBE, MAGILL ENDOTRACH NASAL 5.0MM (10/BX) TELFLX
		242473	112080055	TUBE, MAGILL ENDOTRACH NASAL 5.5MM (10/BX) TELFLX
		204182	112080065	TUBE, MAGILL ENDOTRACH NASAL 6.5MM (10/BX) TELFLX
		191306	112080075	TUBE, MAGILL ENDOTRACH NASAL 7.5MM (10/BX) TELFLX
		132026	112080080	TUBE, ENDOTRACH NASORAL 8.0MM (10/BX)112480 TELFLX
		204183	112080085	TUBE, MAGILL ENDOTRACH NASAL 8.5MM (10/BX) TELFLX
		149619	112082050	TUBE, MURPHY ENDOTRACH CUFF 5.0MM (10/BX) TELFLX
		146252	112082055	TUBE, MURPHY ENDOTRACH CUFF 5.5MM (10/BX) TELFLX
		144244	112082075	TUBE, MURPHY ENDOTRACH CUFF 7.5MM (10/BX) TELFLX
		417647	112082075	TUBE, ENDOTRACH 7.5MM (PARKLAND) (10/BX) TELFLX
		144477	112082090	TUBE, MURPHY ENDOTRACH CUFF 9.0MM(10/BX) TELFLX
		152534	112082095	TUBE, ENDOTRACH NASORAL 9.5MM (10/BX) TELFLX
		149620	112082100	TUBE, ENDOTRACH NASORAL 10.0MM(10/BX) TELFLX
		465321	112480025	TUBE, MURPHY ENDOTRACH CUFF 2.5MM (10/BX) TELFLX
		465322	112480030	TUBE, MURPHY ENDOTRACH CUFF 3.0MM (10/BX) TELFLX
		265413	112480040	TUBE, MURPHY ENDOTRACH CUFF 4.0MM (10/BX) TELFLX
		265414	112480045	TUBE, MURPHY ENDOTRACH CUFF 4.5MM (10/BX) TELFLX
		229796	111780040	TUBE, MURPHY ENDOTRACH CUFF 4.0MM (10/BX) TELFLX
		229800	111780060	TUBE, AGT MURPHY ENDOTRACH CUFF 6.0MM (10/BX) TELFLX

Qty	Unit of Measure	MMS #	MFG Catalog #	Description
		229803	111780075	TUBE, AGT MURPHY ENDOTRACH CUFF 7.5MM (10/BX) TELFLX
		229779	111781045	TUBE, MURPHY ENDOTRACH CUFF 4.5MM (10/BX) TELFLX
		229780	111781050	TUBE, MURPHY ENDOTRACH CUFF 5.0MM (10/BX) TELFLX
		229784	111781070	TUBE, MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX
		144476	112082060	TUBE, MURPHY ENDOTRACH CUFF 6.0MM (10/BX) TELFLX
		144243	112082065	TUBE, MURPHY ENDOTRACH CUFF 6.5MM (10/BX) TELFLX
		130596	112082070	TUBE, MURPHY ENDOTRACH CUFF 7.0MM (10/BX) TELFLX
		134105	112082080	TUBE, MURPHY ENDOTRACH CUFF 8.0MM (10/BX) TELFLX
		146251	112082085	TUBE, MURPHY ENDOTRACH CUFF 8.5MM (10/BX) TELFLX
		265412	112480035	TUBE, MURPHY ENDOTRACH CUFF 3.5MM (10/BX) TELFLX
		421477	170075	TUBE, ENDOTRACH CUFF W/STYL 7.5MM (10/BX)
		1210846	170120	TUBE, ENDOTRACH PRELOADED UNCUFFED SZ2.0 (10/BX)
		349775	504550	TUBE, ENDOTRACH CUFF 5.0MM (10/BX) TELFLX
		349778	504565	TUBE, ENDOTRACH CUFF W/STYLETTE6.5MM (10/BX) TELFLX
		349779	504570	TUBE, ENDOTRACH CUFF W/STYLETTE 7.0MM (10/BX)
		349780	504575	TUBE, ENDOTRACH CUFF W/STYLETTE 7.5MM (10/BX) TELFLX
		349774	506550	TUBE, ENDOTRACH UNCUFF W/STYLETTE 5.0MM (10/BX) TELFLX
		1137459	150030	AIRWAY, LMA MASK REUSE SZ3
		1201753	150040	AIRWAY, LMA MASK RUSBL SZ4
		421466	170055	TUBE, ENDOTRACH CUFF W/STYL STR 5.5MM (10/BX)
		1132959	170140	TUBE, ENDOTRACH UNCUFF STR 4.0MM (10/BX)
		421463	170145	TUBE, ENDOTRACH UNCUFF STR 4.5MM (10/BX)
		189782	100380065	TUBE, MAGILL ENDOTRACH UNCUFF 6.5MM (10/BX)
		678664	112080095	TUBE, MAGILL ENDOTRACH NASAL 9.5MM (10/BX)
		278648	112080100	TUBE, MAGILL ENDOTRACH CUFF 10.0MM (10/BX)

Qty	Unit of Measure	MMS #	MFG Catalog #	Description
		421461	170135	TUBE, ENDOTRACH UNCUFF W/STYL 3.5MM (10/BX)
		401083	150020	TUBE, ENDOTRACH UNCUFF W/STYL 2.0MM (10/BX) TELFLX
		401084	150025	TUBE, ENDOTRACH UNCUFF 2.5MM (10/BX) TELFLX
		191307	112080090	TUBE, MAGILL ENDOTRACH NASAL 9.0MM (10/BX) TELFLX

*Return Affected lot numbers only

*** Any product returned in addition to or in lieu of affected product will be destroyed, without issuance of a credit. The affected lot numbers are listed on the McKesson customer letter.**

If you are on a McKesson truck route, a delivery professional will pick up the affected products, otherwise you will receive UPS return label(s) via email or fax.

I am on a McKesson truck route, please schedule a delivery professional pick up.

Please send my UPS/Return label by Fax or Email. Number of UPS Parcels to be returned: _____

If you have any questions about information provided in this communication, please contact the McKesson Recall Message Center at MMSRecalls@McKesson.com or call (800) 688-8840.

May-2023

URGENT MEDICAL DEVICE RECALL

Type of Action	Recall Notice
Teleflex Reference	EIF-000537
Product Code & Batch/Lot Number	Refer to Appendix 2
Commercial Name	
TOP Endotracheal tube with Cuff	RUSCHELIT® Preformed Nasal Tracheal Tube, Two Eyes without Cuff
Slick Set ® Cuffed Endotracheal Tube and Stylet Set, oral/nasal	RUSCHELIT® Safety Clear Tracheal Tube, oral/nasal
Slick Set ® Uncuffed Endotracheal Tube and Stylet Set, oral/nasal	Safety Clear Tracheal tube (without Cuff)
Flexi-Set Cuffed Endotracheal Tube and Stylet Set, oral/nasal	RUSCHELIT® Preformed Oral Tracheal Tube, Cuffed, Oral, Murphy Eye
Flexi-Set Uncuffed Endotracheal Tube and Stylet Set, oral/nasal	AGT Orotracheal tube
Preformed AGT Oral Endotracheal Tube uncuffed/plain - Murphy	RUSCHELIT® Preformed Nasal Tracheal Tube, Cuffed, Nasal, Murphy Eye
Preformed AGT Nasal Endotracheal Tube uncuffed/plain - Murphy	AGT Nasotracheal tube
Endotracheal Tube oral/nasal uncuffed/plain - Murphy	RUSCHELIT® Safety Clear Plus Tracheal Tube, oral/nasal, Cuffed
Preformed AGT Oral Endotracheal Tube Murphy Eye, High Volume, Low Pressure Cuff	RUSCHELIT® Super Safety Clear Microlaryngeal Tube, oral/nasal
Preformed AGT Nasal Endotracheal Tube Murphy Eye, High Volume, Low Pressure Cuff	RUSCHELIT® Super Safety Clear Tracheal Tube, oral/nasal, Cuffed
Endotracheal Tube oral/nasal Magill, High Volume, Low Pressure Cuff	Flexiset Super Safety Clear Tracheal Tube, oral/nasal with Cuff and Insertion Aid
Endotracheal Tube oral/nasal Murphy Eye, High Volume, Low Pressure Cuff	RUSCHELIT ® Safety Clear Tracheal Tube, oral/nasal, Murphy
RUSCHELIT ® Safety Clear Tracheal Tube, oral/nasal, Magill	Super Safety Clear Tracheal tube
RUSCHELIT® Preformed Oral Tracheal Tube, Two Eyes without Cuff	ENDOSOFT Tracheal tube

Dear Customer,

Details of affected devices

Teleflex LLC has initiated a voluntary recall for the above listed products; refer to Appendix 2 for product code and lot number information.

Description of the problem & immediate actions required

This voluntary recall has been initiated due to reports of disconnection of the 15mm connector from the endotracheal tube (ET tube) for the affected products. There is the possibility of oxygen desaturation and in that event any immediate and long-term health consequences are dependent on the degree and duration of desaturation, which may include serious injury or death.

Where patients are undergoing mechanical ventilation in either the operating room or critical care settings, the ventilation devices to which the affected products are connected are designed to alarm immediately upon a disconnection in the breathing circuit, alerting the clinician to reattach the connector. Additional standards of care such as pulse oximetry also alert clinicians to desaturation within seconds of its occurrence, again permitting prompt reattachment of the connector.

For product in situ, Teleflex advises clinical staff to ensure the 15 mm connector is seated firmly in the ET tube to prevent disconnection during use per the product instructions for use. Should disconnection occur, reconnect the two components promptly and securely in the manner described in the product instructions for use. Clinical staff may wish to consider replacing the device, making sure to evaluate on a case by case basis the risks associated with extubation and reintubation.

As of 14-April-2023, Teleflex received 173 complaints reporting connector issues for products in scope of this recall. Of these 173 complaints, 10 reported injury, including eight reports of patient desaturation, and three reports of patient death. Two complaints reported that the patient deaths were unrelated to the disconnection of the device and one complaint reported that it was impossible to determine whether the device contributed to the patient death.

The initial investigation has identified that the disconnection results from intermittent cross-contamination of the 15mm connector with trace amounts of silicone oil. Although this does increase the lubricity of the connection, this does not present additional expected risks.

Our records indicate you have received products that are subject to this recall.

Depending on your device location please adhere to the following Action list:

Device location	Action List Number
Medical facilities (hospitals, medical staff, etc.)	1
Distributors	2

Action list number 1 – Medical facilities

1. We request that you immediately check your inventory for product within the scope of this recall. **Users should cease use and distribution of affected product and immediately quarantine the affected product.**
2. If you have affected product, mark the applicable checkbox on the Acknowledgement Form (Appendix 1) and fax it to 1-855-419-8507, including “Attn: Customer Service”, or e-mail it to recalls@teleflex.com. This will allow us to document the amount of affected product you have on hand for return. A customer service representative will contact you with a Return Goods Authorization (RGA) number and will provide instructions for the return of affected products to Teleflex.
3. If you do not have affected product, mark the applicable checkbox on the Acknowledgement Form (Appendix 1) and return the form to Teleflex at the contact details provided.
4. Teleflex (or your local dealer) will issue a credit note upon receipt of the returned affected product.

Action list number 2 – Distributors

1. Provide a copy of this recall notice to all customers who have received impacted product. Each of your customers is then required to complete the Acknowledgement Form and return it to you.

2. We request that you immediately check your inventory for affected product. **Cease use and distribution of, and immediately quarantine, the affected product.** You may then return all product in scope.
3. ~~As a distributor, you are then required to confirm to Teleflex that you have completed the field activity outlined in actions 1 and 2 of this Action List Number 2. Upon completion of your actions, please forward the completed Acknowledgement Form to the e-mail address below.~~ **Important** – Please ensure you only list batch numbers in scope of this recall notice when completing this form.
4. If you have further distributed product outside of your country, please notify Teleflex Customer Service by return e-mail to the e-mail address below.

Adverse reactions or quality problems experienced with the use of this product should be reported to Teleflex Customer Service using the contact information below or may also be reported to the FDA's MedWatch Adverse Event Reporting program either by phone at [1-888-INFO-FDA \(1-888-463-6332\)](tel:1-888-INFO-FDA) or online at <https://www.fda.gov/safety/medwatch-fda-safety-information-and-adverse-event-reporting-program>.

Transmission of this Field Safety Notice

This notice should be passed on to all persons who need to be aware within your organization or to any organization where the potentially affected devices have been transferred.

Please consider end users, clinicians, risk managers, supply chain/distribution centres, etc., in the circulation of this notice. Please maintain awareness of this notice until all required actions have been completed in your organization.

Contact reference person

Should you require further information or support concerning this issue, please contact Customer Service via email, phone, or FAX using the details provided below. Customer Service hours of operation are 8am–7pm EST.

Customer Service:

Contact: Customer Service

FAX: 1-855-419-8507

Telephone: 1-866-396-2111

Email: Recalls@teleflex.com

Teleflex is committed to providing high quality, safe and effective products. We regret any inconvenience this action may cause your operations. If you have any other questions, please contact your local Teleflex sales representative or Teleflex Customer Service.

For and on behalf of Teleflex,

Padraig Hegarty

Padraig Hegarty VP, Global QA (Manufacturing)