


Temperature Monitoring to Meet Your Needs

Basic Monitoring

- Record Min/Max Temps
- Programmable Alarm
- Calibrated (cert. included)


More Info



Data Logging

- Programmable Logging Interval
- USB Download (no software)
- Alarm History Reporting


More Info



Cloud Monitoring

- Cloud-Based Applications
- Monitor Locally or Remotely
- Secure Cloud-Based Storage



More Info



Basic Monitoring

Accurately monitor temperatures and record high/low readings.

Record Min/Max Temps
Programmable Alarm
Calibrated (cert. included)



TM-5650

[Web/Pricing](#)

- Triple display simultaneously shows minimum, maximum and current temperatures (°F/°C switchable).
- Alarm is programmable in 1° increments, sounds continuously when temperature rises above/falls below set points, and can be manually turned off.
- Glycol-filled (GRAS by the FDA) bottle insulates the probe sensor from rapid temperature changes.
- 10' micro-cable permits coolers and refrigerator doors to close on it.
- Flip-open stand and wall mount.

Temperature Range: -50° to 70°C (-58° to 158°F)

Resolution: 1°

Accuracy: ±1°

Dimensions: 5.25" x 2.75" x .75" @ 4 oz.


Probe/Bottle: 2.5" x 1" bottle w/ sealed probe; 10' cord

Certificate: multi-point calibration; individually-numbered

Data Logging

Accurately monitor & log temperatures; with USB downloadable data.

Programmable Logging Interval
USB Download (no software)
Alarm History Reporting



TM-6430

[Web/Pricing](#)

- 1 bottle (single zone)
- 525,600 data points
- memory can be cleared



TM-6431

[Web/Pricing](#)

- 2 bottles (dual zone)
- 525,600 data points
- memory can be cleared



TM-6440

[Web/Pricing](#)

- 1 bottle (single zone)
- 1,051,200 data points
- data is locked (21 CFR 11 compliant)



TM-6441

[Web/Pricing](#)

- 2 bottles (dual zone)
- 1,051,200 data points
- data is locked (21 CFR 11 compliant)

- Programmable logging interval from 1 minute to 24 hour observations.
- USB Download (use any USB flash drive/memory stick).
- Smart-Alarm™ (signals until acknowledged).
- Triple display simultaneously shows minimum, maximum and current temperatures (°F/°C switchable).

Temperature Range: -50° to 70°C (-58° to 158°F)

Resolution: .01°

Accuracy: ± .25°

Dimensions: 4.25" x 2.75" x .75" @ 4 oz


Probes/Bottles: 2.5" x 1" bottle w/ sealed probe; 10' cord

Certificate: multi-point calibration; individually-numbered

Cloud Monitoring

Accurately monitor & log temperatures; with Wi-Fi transmission and cloud storage.

Cloud-Based Applications
Monitor Locally or Remotely
Secure Cloud-Based Storage



+



TM-6500

[Web/Pricing](#)

- No local software required; web, Android, iOS data visibility.
- Connects easily to Wi-Fi and provides continuous data transmission to the cloud-based interface.
- Connect the App and view current conditions, control alarms, view history, run reports, and so much more.
- Alerts: temperature alarm, connectivity interruption, low battery.
- Secure, cloud-based data link.
- Receive mobile push/e-mail/text notifications.

Purchase



Create*



Connect



*Requires TraceableLIVE® Subscription

Temperature Range: -50° to 60°C (-58° to 140°F)

Resolution: .01°

Accuracy: ± .25°

Dimensions: 6" x 2.75" x .75" @ 6 oz

Probe/Bottle: 2.5" x 1" bottle w/ sealed probe; 9' cord

Certificate: multi-point calibration; individually-numbered



Vaccinator 500 Go-Bag

Supplies for off-site clinics...
just at staff and vaccine.

Products are available direct or thru any of our fine distributors. (GSA, Amazon, Grainger, etc.)



608-526-6901



info@vericor-med.com



Holmen, WI USA



VeriCorMed.com

Systems

- » Patient Surge
- » Medical Kits
- » Vaccination

Products

- » Organizers & Cases
- » Cool Cubes™
- » Medical Products

Services

- » Cache Conversion
- » System Training
- » System Replenishment



Keep product at the right temperature for days...
without ice or electricity!

- Vaccine (fresh or frozen)
- Blood Products
- Medicine
- Biospecimens



Cold Chain Management

Cool Cubes™ are safe, easy, container and pack-out systems for temperature sensitive product.