

Traceable® Excursion-Trac™ Dual Zone Temperature Monitoring & Data Logger Kit

Model: TM-6431

The Traceable® Excursion-Trac™ Dual Zone Temperature Monitoring & Data Logger Kit provides accurate temperature readings, ideal for monitoring two temperatures while storing vaccines or other important samples. The data logging feature monitors the high and low readings in 1 minute intervals and is capable of storing over 500,000 data points. Memory on the unit can be cleared after downloading data to a USB stick. VeriCor is a proud distributor of Traceable® products.



Features

- Traceable to NIST for accuracy
 - Certificate included
- Records high/low readings
- Signals out-of-range conditions (visual & audible)
- Stores time and dates of excursions
- USB download (CSV file)

Specifications

Temperature Range
-50° to 70°C (-58° to 158°F)

Resolution
0.01°

Accuracy
± 0.25°C

Dimensions & Weight
2¾" x 4¼" x ¾" @ 5oz

10' Probe Cables/2 Probes



USB Capability

*plug in flash drive
to download
(flash drive not
included)



Go to Web Page >>



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®