

# Traceable® Memory-Loc™ Dual Zone Temperature Monitoring & Data Logger Kit

Model: TM-6441

The Traceable® Memory-Loc™ Dual Zone Temperature Monitoring & Data Logger Kit provides accurate temperature readings, ideal for monitoring two temperatures while storing vaccines or other important samples. The data logging feature monitors the high and low readings in 1 minute intervals and is capable of storing over 1 million data points. Raw data is locked and may not be cleared or changed on base unit. VeriCor is a proud distributor of Traceable® products.



## Features

- Traceable to NIST for accuracy
  - Certificate included
- Records high/low readings
- Signals out-of-range conditions (visual & audible)
- Stores time and dates of excursions
- USB download (CSV file)

## Specifications

*Temperature Range*  
-50° to 70°C (-58° to 158°F)

*Resolution*  
0.01°

*Accuracy*  
± 0.25°C

*Dimensions & Weight*  
2¾" x 4¼" x ¾" @ 5oz

10' Probe Cables/2 Probes



USB Capability

\*plug in flash drive to download (flash drive not included)



Go to Web Page >>



info@VeriCorMed.com  
(866) 469-6019

For more information & pricing: [www.VeriCorMed.com](http://www.VeriCorMed.com)

**VERICOR** LLC  
Response-Friendly®