

Traceable® Temperature Monitoring Kit

Model: TM-4127

The Traceable® Temperature Monitoring Kit is the world's most popular refrigerator/freezer thermometer. This calibrated, digital thermometer has a detachable, GRAS, glycol-encased probe (recommended by the CDC) to insulate the readings from rapid temperature changes (i.e. when a door is opened). The triple display simultaneously shows minimum, maximum and current temperatures for effective monitoring. VeriCor is a proud distributor of Traceable® products.



Features

- Traceable to NIST for accuracy
 - Certificate included
- Records high/low readings
- Signals out-of-range conditions (visual & audible)
- 10 foot micro-thin cable
- Flip-open stand & wall mount

Specifications

Temperature Range
-50° to 70°C (-58° to 158°F)

Resolution
1°

Accuracy
± 1° C

Dimensions & Weight
2¾" x 4¼" x ¾" @ 4oz

Back Profile



Side Profile



Go to Web Page >>



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®