

TraceableLIVE® Dual Zone Temperature Monitoring & Data Logger Kit

Model: TM-6500

The TraceableLIVE® Dual Zone Temperature Monitoring & Data Logger Kit allows you to monitor sensitive products wirelessly directly from your computer, tablet or smartphone. Dual glycol filled bottle probes provide accurate temperature readings, ideal for monitoring two temperatures while storing vaccines or other important samples. Audible and visual alarms help to identify any out-of-range conditions. VeriCor is a proud distributor of Traceable® products.



*REQUIRES TRACEABLELIVE® SUBSCRIPTION

Features

- Provides continuous data transmission to the cloud-based interface via WiFi
- Traceable to NIST for accuracy
 - Certificate included
- Records high/low readings
- Signals out-of-range conditions (visual & audible)

Specifications

Temperature Range
-50° to 60°C (-58° to 140°F)

Resolution
0.01°

Accuracy
± 0.25°C

Dimensions & Weight
6" x 2¾" x ¾" @ 6oz

9' Probe Cables/2 Probes



TraceableLIVE® APP

*monitor your environment
wherever you go



Go to Web Page >>



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

VERICOR LLC
Response-Friendly®