

Upright Medical Organizer

Model: MC-UMO

The Upright Medical Organizer (UMO) epitomizes VeriCor's Response-Friendly™ medical equipment. This ultimate organizer, made of durable vinyl, facilitates the rendering of immediate care anywhere. Thanks to the dual-use EZ Hang System™, this organizer can be hung on doors, shelter frames or on the self-supporting stand (sold separately).



EZ Hang System™



Clear Pockets



Double-Zippers



Features

- Light, durable, roll-it-up design
- Clear, double-zipped pockets
- Carry handle when rolled up
- Pre-attached straps and hooks enable various hanging capabilities (EZ Hang System™)
- Highly visible ID holder
- Structural reinforcement rod

Specifications

Dimensions (L x W)
67" x 26" x 5"

Weight
7 lbs.

Storage Pockets
20 clear

Additional Accessories:
- Stand (MK-SUMO)
- Protective Bag (MK-BUMO)



1 - lay flat



2 - roll up



3 - clip and carry

* product for demo purposes only

[Go to Web Page >>](#)



info@VeriCorMed.com
(866) 469-6019

For more information & pricing: www.VeriCorMed.com

